


# INSTALLATION GUIDE

## Phone adapter





# Table of contents


Intro	4
Phone adapter overview	5
Connecting to power	6
Connecting to Phone line	7
Prepare SoundGate for phone calls	9
Using the Phone adapter and support	10
Product approval precautions and markings	11
Warranty	14
Warnings	15

## Intro



This user guide explains how to install the Phone adapter for use with the SoundGate. For instructions how to use your Phone adapter please consult the Phone adapter User Guide.

Before your SoundGate can be used for conversations over your land line phone, the Phone adapter must be connected to the phone line and a power outlet.

All needed parts for the installation are included in two numbered boxes.



## Phone adapter overview



### Status indicator

Light	Indication
Constant	SoundGate connected
Slow fading flash	Searching for SoundGate
Flashing	Pairing or Telephone ringing
No light	Power off




## Connecting to power

**Begin the installation by connecting the Phone adapter to a power outlet.** Use the Power adapter from box no. 1.

- Ensure that the wall socket is an appropriate voltage for the power adapter
- Plug the power adapter into the wall socket
- Plug the power cord into the power connection on the back of the Phone adapter
- The status indicator will flash rapidly for a few seconds to confirm that the power is on

After a few seconds the status indicator will then slowly flash and fade while it is searching for a SoundGate. When a SoundGate is found the status light will be steadily illuminated. **Note:** See page 9 for pairing instructions.





## Connecting to Phone line

**Connect the Phone adapter to the Phone line and phone.** Use the Telephone cord from box no. 2.

**Note:** The connector plugs are the same on both ends of the Telephone cord.

The existing land line connection between the wall outlet and your phone (or base station) must be disconnected at the phone end:

1. Unplug the phone line from the phone (or base station).\*
2. Plug the phone line in the Phone adapter "LINE" connector (BLACK).
3. Plug the Tele cord in the Phone adapter "PHONE" connector (YELLOW).
4. Plug the other end of the Telephone cord into the phone (or base station).\*

The installation is now complete. Before use the SoundGate and the Phone adapter must be paired (see page 9).

**Note:**

\* If your phone or base station does not have a detachable plug for the land line connection, you should either try using an alternative phone or a technician must be contacted to complete the installation.





## Prepare SoundGate for telephone usage

The very first time your SoundGate is used with the Phone adapter, the two devices must be paired together.

### Pairing

1. In order to pair the Phone adapter and the SoundGate, the devices must be placed close together (less than 30 cm/1 ft between devices).
2. Press the SoundGate Bluetooth® key for 5 seconds until the blue light starts flashing rapidly.
3. Press the red phone button on the Phone adapter for 5 seconds until the red status indicator flashes. Wait one minute.
4. After one minute the Phone adapter will be paired with the SoundGate. The red status indicator on the Phone adapter will show a solid light.



Status

Pairing should be completed only once and will normally not have to be repeated.

## Using the Phone adapter and support

For use of the land line phone with the Phone adapter and SoundGate, refer to the Phone adapter user guide.

Place the Phone adapter unit in a suitable spot near your land line phone or base station.

If you don't want to receive phone calls with your SoundGate, turn off the Phone adapter by disconnecting the power.

For normal use, the Phone adapter uses very little power so it is not necessary disconnect the power, except for long periods of time when it is not in use.

The Phone adapter has a range of approx. 10 meters (30 feet) depending on furnitures and placement. Within this range you can make hands free phone conversations with SoundGate.

For further support or tips you can check

**[www.bernafon.com](http://www.bernafon.com)**

## Product approval precautions and markings

The emission power from the Phone adapter is below international emission limits for human exposure.

This device complies with Part 15 of the FCC Rules.

This equipment complies with Part 68 of the FCC rules and the requirements adopted by the ACTA.

**CAUTION:** Any changes or modifications in construction of this device which are not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

To comply with FCC RF exposure requirements, the device and the antenna for this device must be installed to ensure a minimum separation distance of 20 cm or more from a person's body. Other operating configurations should be avoided.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced technician for help.

IC statement:

Operation is subject to the following two conditions:

- 1) this device may not cause interference.
- 2) this device must accept any interference, including interference that may cause undesired operation of the device.

IC radiation exposure statement: This equipment complies with IC radiation exposure limits set forth for an uncontrolled environment. End users must follow the specific operating instructions for satisfying RF exposure compliance. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

This Phone adapter is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC.

Declaration of conformity is available at:

Bernafon AG

Morgenstrasse 131

3018 Bern

Switzerland

[www.bernafon.com](http://www.bernafon.com)

**CE** This product is intended for connection to the public switched analogue telephone network and may be operated in EU and EFTA countries.

Only use the power adapter that was originally supplied with your Phone adapter.

 **Bluetooth®**



## International warranty

The Phone adapter has a limited warranty, covering defects in materials and workmanship. This warranty covers the Phone adapter itself, but not accessories such as charger units, cables, etc. This warranty becomes invalid if a defect is a result of misuse or maltreatment.

This warranty is void if non-authorized personnel have repaired the Phone adapter. Please review the warranty card with your hearing care professional and make sure that your warranty card is completed.

The above warranty does not affect any legal rights that you might have under applicable national legislation governing sale of consumer products. Your hearing care professional may have issued a warranty that goes beyond the clauses of this limited warranty. Please consult your hearing care professional for further information.

## Warnings

The Phone adapter is not a toy and therefore should be kept out of the reach of children and anyone else who might swallow parts or otherwise cause injury to themselves. Special attention should be paid to the smaller components in order to prevent children from swallowing them and choking.

If a part is swallowed, see a doctor immediately.

The safety of the use of the Phone adapter with an auxiliary input cable is determined by external signal source. When the input cable is connected to mains operated equipment, this equipment must comply with IEC-60065, IEC-60601 or equivalent safety standards.

The Phone adapter is designed to comply with the most stringent Standards of International Electromagnetic Compatibility. However, the Phone adapter might cause interference with other medical devices. Please contact relevant personnel to obtain a permission before using the Phone adapter in hospitals or similar venues.

The Phone adapter must never be washed or immersed in water or other liquids.











## **Australia**

Bernafon Australia Pty. Ltd.  
12/97 Castlemaine Street  
Milton QLD 4064  
Freecall 1800 809 111  
Phone +61 7 3250 0300  
Fax +61 7 3250 0372  
[www.bernafon.com.au](http://www.bernafon.com.au)

## **Canada**

Bernafon Canada Ltd.  
500 Trillium Drive, Unit 15  
Kitchener, Ontario  
Canada N2R 1A7  
Phone +1 519 748 6669  
Toll Free +1 800 265 8250  
Toll Free Fax +1 888 748 9158  
[www.bernafon.ca](http://www.bernafon.ca)

## **Distributor:**

### **Switzerland**

Bernafon AG  
Morgenstrasse 131  
3018 Bern  
Phone +41 31 998 15 15  
Fax +41 31 998 15 90  
[www.bernafon.com](http://www.bernafon.com)



## **New Zealand**

Bernafon New Zealand Ltd.  
Level 1, Building F  
27-29 William Pickering Drive  
Albany, Auckland 0632  
Toll Free 0800 44 22 57  
Telephone +64 9 415 7917  
Fax +64 9 415 7916  
[www.bernafon.co.nz](http://www.bernafon.co.nz)

## **United Kingdom**

Bernafon UK  
Cadzow Industrial Estate  
Off Low Waters Road  
Hamilton  
ML3 7QE Scotland  
Phone +44 1698 285 968  
Fax +44 1698 421 456  
[www.bernafon.co.uk](http://www.bernafon.co.uk)

## **United States of America**

Bernafon, LLC  
200 Cottontail Lane, Bldg. B  
Somerset, NJ 08873  
Phone +1 888-941-4203  
Fax +1 732-560-4877  
[www.bernafon-us.com](http://www.bernafon-us.com)

[www.hear-it.org](http://www.hear-it.org)

**For more information about  
hearing and hearing loss**

**bernafon<sup>®</sup>**

*Your hearing • Our passion*