


# INSTALLATION GUIDE

## TV adapter





# Table of contents

Purpose of this guide	4
TV adapter overview	5
Connecting to power	6
Connecting to TV	7
Optional adapter for connecting to TV	8
Prepare SoundGate for TV listening	9
Using the SoundGate and support	10
Product approval precautions and markings	11
International warranty	13
Warnings	14


## Purpose of this guide

This is an easy guide through the installation of the TV adapter. Once the installation is complete this manual is no longer required.



Before your Bluetooth® enabled instruments can receive sound from the TV, the TV adapter must be connected to the TV and a power outlet.

All needed parts for the installation are included in the packaging.

The connection parts are packed in numbered boxes from 1 to 3. The numbers are referring to each step of the installation.



## TV adapter overview




### Power indicator

Light	Meaning
Solid	Power on
Off	Power off

### Status indicator

Light	Meaning
Solid	Connected to SoundGate and streaming TV sound
Off	Stand-by
Flashing	Pairing





## Connecting to power

**Begin the installation by connecting power to the TV adapter.**

Use the power adapter from box no. 1.

- Plug the power adapter into the wall socket or another power source.
- Plug the power cord into the power connection marked **A**.
- Check that the power is turned on and there is a solid green light in the power indicator on the front of the TV adapter.





## Connecting to TV


**Connect the TV to the TV adapter.**


**Use the TV cord from box no. 2.**

- Connect the black plug on the TV cord to the connector marked **B**.
- Connect the red and white connectors to the "AUDIO OUT" terminals on the TV (typically marked "L – R" "AUDIO OUT" and found on the TV back or bottom panel).

If the TV has no "L – R" RCA Phono "OUTPUT", an adapter should be used for connecting the TV cord. Then proceed to step 3.

The installation is now complete and can be checked by using your SoundGate (refer to the user guide).






## Optional adapter for connecting to TV

**Head phone connection** using the jack adapter.

- Connect the TV cord to the jack adapter and plug it into the TV's head phone output.

The head phone output may be found on the side or front of the TV.

Note that the TV sound over the speaker may turn off while the jack adapter is plugged in.



## Prepare SoundGate for TV listening

The very first time your SoundGate is used with the TV adapter, the two devices must be paired together.

### Pairing

1. In order to pair the TV adapter and the SoundGate, the devices must be placed close together. Place the SoundGate right next to the TV adapter (less than 30 cm /1 ft between devices).
2. Press the SoundGate Bluetooth® key for 5 seconds until the blue light starts flashing rapidly.
3. The two devices are now establishing a connection with each other. Wait one minute.
4. After one minute the TV adapter connects to the SoundGate and the pairing is complete. The orange STATUS indicator on the TV adapter will show a solid light.



Status

Pairing should be completed only once and will normally not have to be repeated.

## Using the SoundGate and support

Place the TV adapter in a suitable spot close to the TV. To use the SoundGate with the TV adapter, refer to the SoundGate user guide.

### **Turn off the TV adapter**

The TV adapter uses very little power so it is not necessary to turn it off.

### **TV sound range**

The TV adapter has a range of about 10 meters / 30 feet, depending on room furnitures and placement.

For further support or tips you can check the SoundGate information on **[www.bernafon.com](http://www.bernafon.com)**

## Product approval precautions and markings

The emission power from the TV adapter is below international emission limits for human exposure.

This device complies with Part 15 of the FCC Rules.

Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

To comply with FCC RF exposure requirements, the device and the antenna for this device must be installed to ensure a minimum separation distance of 20 cm or more from a person's body. Other operating configurations should be avoided.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to

radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This TV adapter is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC.

Declaration of conformity is available at:

Bernafon AG  
Morgenstrasse 131  
3018 Bern  
Switzerland  
[www.bernafon.com](http://www.bernafon.com)

**CE 0678**



Only use the power adapter that was originally supplied with your TV adapter.

## International warranty

The TV adapter has a limited warranty, covering defects in materials and workmanship. This warranty covers the TV adapter itself, but not accessories such as charger units, cables, etc. This warranty becomes invalid if a defect is a result of misuse or maltreatment.

This warranty is void if non-authorized personnel have repaired the TV adapter. Please review the warranty card with your hearing care professional and make sure that your warranty card is completed.

The above warranty does not affect any legal rights that you might have under applicable national legislation governing sale of consumer products. Your hearing care professional may have issued a warranty that goes beyond the clauses of this limited warranty. Please consult your hearing care professional for further information.

## Warnings

The TV adapter is not a toy and therefore should be kept out of the reach of children and anyone else who might swallow parts or otherwise cause injury to themselves. Special attention should be paid to the smaller components in order to prevent children from swallowing them and choking.

If a part is swallowed, see a doctor immediately.

The safety of the use of the TV adapter with an auxiliary input cable is determined by external signal source. When the input cable is connected to mains operated equipment, this equipment must comply with IEC-60065, IEC-60601 or equivalent safety standards.

The TV adapter is designed to comply with the most stringent Standards of International Electromagnetic Compatibility. However, the TV adapter might cause interference with other medical devices. Please contact relevant personnel to obtain a permission before using the TV adapter in hospitals or similar venues.

The TV adapter must never be washed or immersed in water or other liquids.



**Distributor:**

Bernaфон AG  
Morgenstrasse 131  
3018 Bern  
Switzerland  
Phone +41 31 998 15 15  
Fax +41 31 998 15 90  
[www.bernaфон.com](http://www.bernaфон.com)

CE 0678

**Canada**

Bernaфон Canada Ltd.  
500 Trillium Drive, Unit 15  
Kitchener, Ontario  
Canada N2R 1A7  
Phone +1 519 748 6669  
Toll Free +1 800 265 8250  
Toll Free Fax +1 888 748 9158  
[www.bernaфон.ca](http://www.bernaфон.ca)

**United States of America**

Bernaфон, LLC  
200 Cottontail Lane, Bldg. B  
Somerset, NJ 08873  
Phone +1 888-941-4203  
Fax +1 732-560-4877  
[www.bernaфон-us.com](http://www.bernaфон-us.com)

[www.hear-it.org](http://www.hear-it.org)

**For more information about  
hearing and hearing loss**