

# SOUNDGATE

SoundGate for TV listening

# USER GUIDE

## Intro

This user guide explains how the SoundGate is used for TV listening.


The TV sound is transmitted by the TV adapter to the SoundGate via Bluetooth® technology.

Make sure that the TV adapter is correctly connected to the TV. Refer to the installation guide.

When the TV adapter is turned on all operations for TV listening are made from the SoundGate.

TV listening is possible within a range of approx. 10 meters / 30 feet from the TV.



## SoundGate overview



## Prepare SoundGate for TV listening

The very first time your SoundGate is used with the TV adapter, the two devices must be paired together.

### Pairing

1. In order to pair the TV adapter and the SoundGate, the devices must be placed close together. Place the SoundGate right next to the TV adapter (less than 30 cm /1 ft between devices).
2. Press the SoundGate Bluetooth® key for 5 seconds until the blue light starts flashing rapidly. 
3. The two devices are now establishing a connection with each other. Wait one minute.
4. After one minute the TV adapter connects to the SoundGate and the pairing is complete. The orange STATUS indicator on the TV adapter will show a solid light.  Status

Pairing should be completed only once and will normally not have to be repeated.

# TV listening/start and stop

## Before listening

Wear the SoundGate around your neck with the neck strap. You must be within 10 meters / 30 feet from the TV set.

Use the SoundGate audio key for TV listening.



## START sound

Press the audio key.

The SoundGate will connect to the TV adapter. After a few seconds a solid orange light will be shown on the SoundGate and TV adapter.

The hearing instruments will then receive the TV sound.



Status

## STOP sound

Press the audio key.

The SoundGate will disconnect from TV adapter. All orange lights will turn off. The hearing instruments will then return to the previous listening program.

## Notice

The Bluetooth® function on the SoundGate must be activated to allow transmission from the TV. If not active, it will automatically turn on when starting the TV listening. De-activating the Bluetooth® function must be done manually with the Bluetooth® key.

