Instructions for Use

IN-THE-EAR HEARING AIDS

Nevara, Saphira, Juna, Carista, Acriva, Chronos, Inizia







Table of Contents

Hearing Aid Description	9
Step-by-Step Instructions	
for Using Your Hearing Aid	12
Step 1: Inserting the Battery	12
Step 2: Turning the Hearing Aid ON	14
Step 3: Inserting Your Hearing Aid	15
Step 4: Changing the Volume	17
Step 4a: Changing the Volume with	
the Optional Volume Control Wheel	17
Step 4b: Changing the Volume with	
the Optional Push Button	18
Step 5: Changing the Programs	20
Step 6: Muting Your Hearing Aid	21
Step 7: Removing Your Hearing Aid	22
Step 8: Turning the Hearing Aid OFF	22
Step 9: Changing the Battery	23
Features of Your Hearing Aid	24
Caring for Your Hearing Aid	25
General Care Instructions	25
Daily Care Instructions	26
Specific Care Instructions	27
Common Problems and Solutions	28
Accessories	30

Warnings	32
International Warranty	39
Mobile Phone	41
Technical Information	43
Information and Explanation of Symbols	47

The following Bernafon hearing aid models are covered within this booklet:

Ju	na 9	7
\neg	11 10 1	TEPD/ITED, with wireless functionality TCPD/ITCD, with wireless functionality
	JU9 (JU9 (CICP/CICx, with wireless functionality CIC
	JU9 I JU9 I	ICx, with wireless functionality IC
	JU7 I JU7 I	TEPD/ITED, with wireless functionality TCPD/ITCD, with wireless functionality
	JU7 (JU7 (TC CICP/CICx, with wireless functionality CIC ICx, with wireless functionality CIC ICx, with wireless functionality IC TEPD/ITED, with wireless functionality TCPD/ITCD, with wireless functionality TC CICP/CICx, with wireless functionality CIC
Ac	riva 9	17
\neg	A RQ I	ITED, with wireless functionality ITCPD/ITCD, with wireless functionality ITC CICP/CIC IIC ITED, with wireless functionality ITCPD/ITCD, with wireless functionality
	AK/I	ITED, with wireless functionality ITEDD/ITCD, with wireless functionality ITC CICP/CIC

	ronos 9 7 5 CN9 ITED CN9 ITCPD CN9 ITCD, with wireless functionality CN9 ITCP
	CN9 CICP CN9 CIC CN7 ITED CN7 ITCPD
	CN7 ITCD, with wireless functionality CN7 ITCP CN7 CICP CN7 CIC
	CN5 ITED CN5 ITCPD CN5 ITCD, with wireless functionality CN5 ITCP CN5 CICP CN5 CICP
Sa	phira 5 3
	SA5 ITEPD/ITED, with wireless functionality SA5 ITCPD/ITCD, with wireless functionality SA5 ITC SA5 CICP/CICx, with wireless functionality
	SA5 CIC SA5 IICx, with wireless functionality SA5 IIC SA3 ITEPD/ITED, with wireless functionality
	SA3 ITCPD/ITCD, with wireless functionality SA3 ITC

SA3 CICP/CICx, with wireless functionality SA3 CIC
rista 5 3 CA5 ITED, with wireless functionality CA5 ITCPD/ITCD, with wireless functionality CA5 ITC CA5 CICP/CIC CA3 ITED, with wireless functionality CA3 ITCPD/ITCD, with wireless functionality CA3 ITC CA3 CICP/CIC
zia 3 1 IN3 ITED IN3 ITCP IN3 CICP IN3 ITCD IN3 ITC IN3 ITC IN1 ITC IN1 ITCD IN1 ITCD IN1 ITC IN1 ITC
NE1 ITED, with wireless functionality NE1 ITCD, with wireless functionality NE1 ITC NE1 CICP/CICx, with wireless functionality



This booklet guides you on how to use and maintain your new hearing aid. Please read the booklet carefully including the warning section. This will help you to achieve the full benefit of your new hearing aid.

Your hearing care professional has adjusted the hearing aid to meet your needs. If you have additional questions, please contact your hearing care professional.

Indication for Use / Intended Purpose

The hearing aids are intended to amplify and transmit sound to the ear and thereby compensate for impaired hearing. The hearing aids are intended to be used by children (> 36 months) and adults.

IMPORTANT NOTICE


The hearing aid amplification is uniquely adjusted and optimized to your personal hearing capabilities during the fitting performed by your hearing care professional.

Hearing Aid Description

For your in-the-ear instrument, there are several different styles and sizes. Please identify your style. This will make it easier for you to navigate through this booklet.

Please note that your hearing aid may look slightly different compared to the illustrations on page 10 and 11.

ITEPD, ITED, ITCPD, ITCD, ITCP, and ITC Models




5 26

3

ITEPD/ITED

ITCPD/ITCD



ITCP/ITC

- 1 Vent
- 2 Canal
- 3 Sound outlet with wax protection
- 4 Microphone opening with O-cap filter
- 5 Battery door
- 6 Push button* (optional)
- 7 Volume control (optional)

^{*}For ITEPD, ITED, ITCPD, ITCD models in Juna, Saphira, Nevara, Acriva, and Carista families the push button can be programmed for volume control use

CICP, CICx, CIC, IICx, and IIC Models





CICP/CICx/CIC

IICx/IIC

- 1 Vent
- 2 Canal
- 3 Sound outlet with wax protection
- 4 Microphone opening with T-cap filter
- 5 Battery door
- 6 Push button (optional)*
- 7 Removal (pull-out) string


^{*}For Juna, Saphira, and Nevara CICP and CICx the push button can be programmed for volume control use

Step-by-Step Instructions for Using Your Hearing Aid

Step 1: Inserting the Battery

Your hearing aid uses the following battery size:


- · ITEPD, ITED: size 13
- · ITCPD, ITCD, ITCP, ITC: size 312
- · CICP, CICx, CIC, IICx, IIC: size 10
- Gently swing the battery door fully open, but do not force it (A)
- · Remove the sticker from the new battery (B)
- Place the battery into the empty compartment.
 The + sign on the battery should face up. (C)



*

For maximum power, allow the battery to be aired for 30 seconds before placing it into the empty compartment.

The MultiTool can be used for battery change. Use the magnetic end to remove and insert batteries. The MultiTool is provided by your hearing care professional.



IMPORTANT NOTICE

Your MultiTool has a built-in magnet. Keep the MultiTool 30 cm away from credit cards and other magnetically sensitive devices.

Step 2: Turning the Hearing Aid ON

Close the battery door completely. You should notice a click. The hearing aid is now ON.



Click "ON"


Never use excessive force to open or close the battery door.

Step 3: Inserting Your Hearing Aid

Your hearing aid has been programmed individually for your right or left ear.

You will see a color marking on your hearing aid. This can help you distinguish between the left (blue) and right (red) hearing aids.

- When inserting the right hearing aid, hold it with the right hand. When inserting the left hearing aid, hold it with the left hand.
- Hold your hearing aid between your thumb and index finger with the microphone on top. If your hearing aid has a pull-out string, this must be on the bottom. (A)
- Place the canal part of your hearing aid into your ear canal (B)



 Gently pull your earlobe down with the other hand while pushing the hearing aid in until it feels secure and comfortable

*

It takes patience and practice to insert your earmold correctly. If you have difficulty, please consult your hearing care professional.

Step 4: Changing the Volume

Your hearing aid automatically adjusts the volume to the changing sound environment.

Some models allow you to adjust the volume yourself with an optional control on the hearing aid. Please ask your hearing care professional if this function has been made available on your hearing aid.

Step 4a: Changing the Volume with the Optional Volume Control Wheel

To increase the volume, turn the volume control towards the front.

To decrease the volume, turn the volume control towards the back.



The hearing aid will click to confirm that the volume has been changed and then will beep when you reach the programmed maximum or minimum. You may hear hear a beep, when you return to the preset loudness level in your hearing aid.

Step 4b: Changing the Volume with the Optional Push Button

If your hearing aid has been made with a push button, your hearing care professional may be able to program it to change volume.

Two hearing aids

	Left	Right
Short push (less than 1 sec.)	☐ volume dowr	n 🗌 volume up
One hearing aid		
Short push (less than 1 sec.)	☐ volume dowr	١
Medium push (about 1 sec.)	☐ volume up	

The hearing aid will click to confirm that the volume has been changed and then will beep when you reach the programmed maximum or minimum. You may hear a beep, when you return to the preset loudness level in your hearing aid. Please ask your hearing care professional for further information.

*

Your hearing aid automatically returns to a preset loudness level when it starts, the battery is low, or when the program is changed. If you find that the volume is not adequate, your hearing care professional may need to adjust your hearing aid settings.

Step 5: Changing the Programs

If your hearing aid has been made with a push button, your hearing care professional can program it to change programs.

The hearing aid will beep to confirm the program change. The number of beeps you hear will tell you which program you are in.

Up to 4 listening programs can be configured by your hearing care professional depending on your hearing aid. Please ask your hearing care professional about the programs that are available with your hearing aid.

Available Programs

(to be filled in by your hearing care professional)

Program	Beeps	Function
1	1 🖍	
2	2 🖍 🖍	
3	3 111	
4	4 1 1 1 1	

Step 6: Muting Your Hearing Aid

If your hearing aid has been made with a push button, your hearing care professional may be able to program it to mute your hearing aid. Please ask your hearing care professional if this function is available with your hearing aid.

A long press (about 2 seconds) on the button will cause your hearing aid to mute.

To un-mute the hearing aid, give any control on the hearing aid a push.



If you wear one hearing aid (monaural fitting), it is not possible to mute it.

IMPORTANT NOTICE


Do not use the mute function as an off switch, as the hearing aid still draws from the battery in this mode.

Step 7: Removing Your Hearing Aid

- Use your thumb to push up against the bottom (back part) of your ear to loosen the hearing aid
- Grasp the hearing aid at its edge between your thumb and forefinger. Gently remove it from your ear.
- If your hearing aid has a removal (pull-out) string, pull it gently

Step 8: Turning the Hearing Aid OFF

Place your fingernail underneath the front edge of the battery door and lift to open. The hearing aid is now OFF.



Step 9: Changing the Battery

When the battery is running out, the hearing aid will beep at regular intervals. You should be prepared to replace the battery. The time until the hearing aid stops working depends on the battery type and manufacturer. Mercury-free batteries generally cause earlier warning beeps.

Features of Your Hearing Aid

The following features may be activated in your hearing aid. Please consult your hearing care professional to find out how your hearing aid is configured.

□ Telecoil helps you hear better when using a telephone or in buildings where inductive loop systems are installed.



This symbol or a similar sign should be displayed wherever a permanent loop has been installed.

Auto Telephone can automatically activate a phone program in the hearing aid, if your telephone has a dedicated magnet.



The magnet needs to be placed at your telephone next to the sound outlet

Caring for Your Hearing Aid

Healthy ears produce a waxy substance that can clog your hearing aid. Please follow these cleaning instructions to prevent wax build-up and ensure optimal performance of your hearing aid.

For further information on the care of your hearing aid, please consult your hearing care professional.

General Care Instructions

Treat your hearing aid like any other delicate electronic device (check the safety guidelines on page 38) and make sure that the hearing aid does not come in contact with moisture or water. Always clean and dry your hands thoroughly before handling the hearing aid.

IMPORTANT NOTICE

Do not wear your hearing aid while showering or participating in water activities. Do not immerse your hearing aid in water or other liquids. Do not use the battery door as a handle to insert or remove your hearing aid. It is not designed for this purpose.

Daily Care Instructions

- Check your hearing aid for ear wax and wipe it clean with a cloth or tissue
- If necessary, use the cleaning tool to remove wax from the canal and vent openings on your hearing aid
- Open the battery door fully to allow air to circulate during the night
- A dry storage kit is recommended to remove any moisture that may have accumulated in the hearing aid. Always remove zinc-air batteries before drying your hearing aid with a drying kit. Drying out zinc-air batteries shortens their lifetime.

IMPORTANT NOTICE

Use a soft, dry cloth to clean the hearing aid. It must never be washed or immersed in water or other liquids.

Specific Care Instructions

A: Cleaning the Sound Outlet

The sound outlet is an integral part of your hearing aid and should be kept clean and free from earwax at all times.

Your hearing aid is built with a wax protection system. Please ask your hearing care professional for the instructions on how to exchange the wax protection filter on your hearing aid.

B: Microphone Protection System

The microphone is an integral part of your hearing aid and should be kept clean and free from cerumen (earwax) at all times.

Your hearing aid may have a built-in wax protection system for the microphone. Please ask your hearing care professional for the instructions on how to exchange the microphone protection filter on your hearing aid.

Common Problems and Solutions

Hearing aid whistles or squeals

Check that the earpiece is inserted correctly. If this is the case and the hearing aid still whistles or squeals, please consult your hearing care professional.

No volume, level too soft or humming noise

Check if the volume level is too soft and adjust by increasing the volume level with the control on your hearing aid. If you have no volume at all, check whether the hearing aid is in mute position. If the problem still occurs, check that the battery door is closed completely. Also check if the battery is inserted correctly. If the problem still remains, change the battery. If the problem still persists, contact your hearing care professional.

Whirring noise, fading, weak or motor-boating sounds

Open and close the battery door several times or clean the battery contacts carefully with a dry cotton swab. If the problem still occurs, change the battery. If the problem still remains, contact your hearing care professional.

Hearing aid switches from ON to OFF periodically

Your battery is running down. Please change the battery.

Hearing aid beeps without any action from you Your battery is running down. Please change the battery.

Other problems with your hearing aid

If other problems occur with your hearing aid which are not listed, contact your hearing care professional.

Accessories


that may be purchased to enhance your hearing aid. Depending on the hearing aid family and local regulations, the following accessories are available:


Remote control
SoundGate communication device (for wireless connection to mobile phones, music players, etc.)
SoundGate Mic (for wireless external microphone connection via the SoundGate)
TV Adapter (for wireless connection to your television via the SoundGate)
Phone Adapter (for wireless connection to your

Bernafon offers a wide range of optional accessories

For further information on accessories, please contact your hearing care professional.

landline phone via the SoundGate)





Warnings

You should familiarize yourself fully with the following general warnings and the entire contents of this booklet before using your hearing aid to ensure personal safety and correct use.

Consult your hearing care professional if you experience unexpected operations or events with your hearing aid.

Please note that a hearing aid will not restore normal hearing and will not prevent or improve a hearing impairment resulting from organic conditions. Furthermore, note that in most cases, infrequent use of a hearing aid does not permit a user to attain full benefit from it.

Usage of Hearing Aids

- Hearing aids should be used only as directed and adjusted by your hearing care professional. Misuse can result in sudden and permanent hearing loss.
- Never allow others to wear your hearing aid as incorrect usage could cause permanent damage to their hearing

Choking Hazards and Risk of Swallowing Batteries or Other Small Parts

- Hearing aids, their parts, and batteries should be kept out of reach of children and anyone who might swallow these items, or otherwise cause injury to themselves
- Batteries have occasionally been mistaken for pills.
 Therefore check your medicine carefully before swallowing any pills.
- Most hearing aids can be supplied with a tamperresistant battery drawer upon request. This is strongly recommended for infants, small children, and people with learning difficulties. Ask your hearing care professional if a tamper-resistant battery drawer is available for your hearing aid.

If a battery or other small parts are swallowed, see a doctor immediately.

Battery Use

- Always use batteries recommended by your hearing care professional. Batteries of low quality may leak and cause bodily harm.
- Never attempt to recharge your batteries and never dispose of batteries by burning them. There is a risk that the batteries will explode.

Dysfunction

 Be aware of the possibility that your hearing aid may stop working without notice. Keep this in mind when you depend on warning sounds (e.g., when you are in traffic). The hearing aids may stop functioning, for instance if the batteries have expired or if the tubing is blocked by moisture or ear wax.

Active Implants

- · Caution must be taken with active implants
- If your hearing aid has a wireless transmission, keep the hearing aid at least 15 cm away from the implant, e.g. do not carry it in a breast pocket. In general, please follow the guidelines recommended by the manufacturers of implantable defibrillators and pacemakers on use with mobile phones.

- Your Autophone magnet or MultiTool (which has a built-in magnet) should be kept more than 30 cm away from the implant, e.g. do not carry it in a breast pocket. In general, please follow the guidelines recommended by the manufacturer of implantable defibrillators and pacemakers on use with magnets.
- If you have an active brain implant, please contact the manufacturer of your implantable device for information about the risk of disturbance

Explosives

 The power source in your hearing aid has insufficient energy to cause fire in normal usage conditions. The hearing aid has not been tested for compliance with international standards concerning explosive environments. We recommended not to use your hearing aid in areas where there is a danger of explosions.

X-ray, CT, MR, PET Scanning and Electrotherapy

 Remove your hearing aid for example during X-ray, CT / MR / PET scanning electrotherapy or surgery as your hearing aid may be damaged when exposed to strong fields

Avoiding Heat and Chemicals

- Your hearing aid must never be exposed to extreme heat e.g., left inside a parked car in the sun
- Your hearing aid must not be dried in microwave ovens or other ovens
- The chemicals in cosmetics, hairspray, perfume, after shave lotion, suntan lotion and insect repellent can damage your hearing aid. Always remove your hearing aid before applying such products and allow time to dry before putting it on.

Power Instrument

 Special care should be exercised in selecting, fitting and using a hearing aid where maximum sound pressure capability exceeds 132 dB SPL (IEC 60318-4), as there may be risk of impairing the remaining hearing of the hearing instrument user

For information of whether your instrument is a power instrument, please ask your hearing care professional.

Possible Side Effects

- Hearing aids, moulds or domes may cause an accelerated accumulation of ear wax
- The otherwise non-allergenic materials used in hearing aids may in rare cases cause a skin irritation or other side effects

Please seek consultation with a physician if these conditions occur.

Interference

 Your hearing aid has been thoroughly tested for interference, according to the most stringent international standards. However, interference with your hearing aid and other devices may occur, (e.g., some mobile telephones, citizens band systems, and shop alarm systems). If this occurs, increase the distance between the hearing aid and the device

Connection to External Equipment

 The safety of the use of the hearing aid with an auxiliary input cable is determined by the external signal source. When the input cable is connected to equipment plugged into a wall outlet, this equipment must comply with IEC 60601, IEC 60065 or equivalent safety standards.

Safety Information

- Hearing aids should be adjusted by a trained hearing care professional
- Never insert cleaning tools into the sound outlet or microphone inlet. This could damage the hearing aid.
- · Remove your hearing aid before sleeping
- Keep your hearing aid in the case for protection when you are not wearing it
- Be aware of the possibility that the directional microphone in your hearing aid may reduce the volume of some warning sounds coming from behind you

International Warranty

Your hearing aid is covered by an international limited warranty issued by the manufacturer from the date of delivery.

This limited warranty covers manufacturing and material defects in the hearing aid itself, but not accessories such as batteries, tubing, ear wax filters etc. Problems arising from improper handling or care, excessive use, accidents, repairs made by an unauthorized party, exposure to corrosive conditions, physical changes in your ear, damage due to foreign objects entering the device, or incorrect adjustments are NOT covered by the limited warranty and may void it

The above warranty does not affect any legal rights that you might have under applicable national legislation governing sale of consumer goods. Your hearing care professional may have issued a warranty that goes beyond the clauses of this limited warranty. Please consult him/her for further information.

If you need service

Take your hearing aid to your hearing care professional, who may be able to sort out minor problems and adjustments immediately.

Date:	Model:	
Warranty period:		
Model L:	Model R:	
Serial no.:	Serial no.:	
Battery size:	Battery size:	
Hearing Centre		

Do not attempt to repair your hearing aid yourself.

Mobile Phone

Some hearing aid users have reported a buzzing sound in their hearing aid when they are using mobile phones, indicating that the mobile phone and hearing aid may not be compatible.

The ANSI C63.19 standard determines the prediction of compatibility between a specific hearing aid and a mobile phone by adding the numerical value of the rating for the hearing aid immunity to the numerical value of the rating for the mobile phone emissions. A sum of 4 would indicate that the combination of wireless device and hearing aid is usable; a combined rating that equals 5 would provide normal use; a combined rating of 6 or greater would indicate excellent performance.

Whereas all hearing aids have acoustic coupling, only the larger instruments have the physical space for telecoil (inductive) coupling. These two types of coupling have different rating scales (M1-M4 for acoustic coupling and T1-T4 for telecoil coupling, respectively) and both ratings are therefore relevant when predicting the compatibility of a particular instrument.

For a hearing aid with both acoustic coupling and telecoil coupling with a rating of M4/T2 and with a telephone rating of M3/T3, the combined rating is 7 (M4 + M3) for the acoustic coupling and 5 (T2 + T3) for the telecoil coupling. According to the guideline given above, both types of coupling will thereby be acceptable, with the acoustic coupling indicating excellent performance and the telecoil coupling indicating normal use. The above equipment performance measurements, categories and system classifications are based upon the best information available, but it cannot be guaranteed that all users will be satisfied.

The immunity of all Bernafon hearing instrument models covered by these instructions for use is at least M2/T2.

IMPORTANT NOTICE

The performance of individual hearing instruments may vary with individual cell phones. Therefore, please try this hearing instrument with your mobile phone or, if you are purchasing a new phone, be sure to try it with your hearing instrument prior to purchase. For additional guidance, please ask your cell phone provider for the booklet entitled "Hearing Aid Compatibility with Digital Wireless Cell Phones".

Technical Information

The hearing aid contains a radio transmitter using short range magnetic induction technology working at 3.84 MHz. The magnetic field strength of the transmitter is < -42 dBµA/m @ 10m.

The emission power from the radio system is well below international emission limits for human exposure. For comparison, the radiation of the hearing aid is lower than unintended electromagnetic radiation from for example halogen lamps, computer monitors, dishwashers, etc.

The hearing aid complies with international standards concerning Electromagnetic Compatibility.

Due to the limited space available on the hearing aid all relevant approval markings are found in this document.

Juna hearing aid models: JU9 ITEPD, JU9 ITED, JU9 ITCPD, JU9 ITCD, JU7 ITEPD, JU7 ITED, JU7 ITCPD, JU7 ITCD;

Acriva hearing aid models: AR9 ITED, AR9 ITCPD, AR9 ITCD, AR7 ITED, AR7 ITCPD, AR7 ITCD; Saphira hearing aid models: SA ITC/ITE; Carista hearing aid models: CA5 ITED, CA5 ITCPD, CA5 ITCD. CA3 ITCD. CA3 ITCD:

Nevara hearing aid models: NE ITC/ITE

FCC ID: U6XF2ITE01 IC: 7031A-F2ITE01

Juna hearing aid models: JU9 CICP, JU9 CICx, JU9 IICx, JU7 CICP, JU7 CICx; Saphira hearing aid models: SA5 CICP, SA5 CICx, SA5 IICx, SA3 CICP, SA3 CICx; Nevara hearing aid models: NE1 CICP, NE1 CICx

Contains FCC ID: U6XF2CIC01 Contains IC: 7031A-F2CIC01

Chronos hearing aid models: CN9 ITCD, CN7 ITCD, CN5 ITCD

FCC ID: U6XFUITE01 IC: 7031A-FUITE01

The device complies with Part 15 of the FCC rules and with Industry Canada's licence-exempt RSSs.

Operation is subject to the following two conditions:

- 1. This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This Class B digital apparatus complies with Canadian ICES-003.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try

to correct the interference by one or more of the following measures:

- · Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the manufacturer or an experienced radio/ TV technician for help.

Hereby, Bernafon AG declares that the hearing aids covered in this booklet are in compliance with Directives 93/42/EEC, 1999/5/EC and 2011/65/EU.

The full text of the declaration of conformity is available at:

Bernafon AG Morgenstrasse 131 3018 Bern Switzerland

C€ 0543



equipment must be handled according to loca regulations.





Information and Explanation of Symbols



The CE marking indicates compliance to all applicable European Directives. The 4-digit numbers after the CE marking correspond to the identification number of notified bodies.



This symbol indicates that the products described in these instructions for use adhere to the requirements for an applied part of Type B of EN 60601-1. The surface of the hearing aid is specified as applied part of Type B.*

^{*} not for Inizia



The crossed-out wheeled bin indicates the European Directive 2011/65/EU on waste of electronic equipment applies. Please recycle your hearing aid and batteries according to your local regulations or return them to your hearing care professional for disposal.



This symbol indicates to follow the instructions for use in this booklet



EMC and Radio communications compliance label Australia and New Zealand



Manufacturer



Keep dry

REF

Catalogue number



Serial number

	Temperature	Humidity	
		<u></u>	
Operating Conditions	+1°C - +40°C	5% – 93%	
Storage and transportation conditions	–25°C – +60°C	Non condensing	



Australia

Bernafon Australia 629 Nudgee Road Nundah QLD 4012 Freecall 1800 809 111 Phone +61 7 3250 0300 info@bernafon.com.au

Canada

Bernafon Canada 500 Trillium Drive, Unit 15 Kitchener, Ontario Canada N2R 1A7 Canada +1 800 265 8250 contact@bernafon.ca



Bernafon AG Morgenstrasse 131 3018 Bern Switzerland Phone +41 31 998 15 15 info@bernafon.com www.bernafon.com

New Zealand

Bernafon New Zealand Level 1, Building F 27-29 William Pickering Drive North Harbour, Albany Auckland 0632 Toll Free 0800 442 257 Phone +64 9 415 7917 Fax +64 9 415 7916

South Africa

Bernafon South Africa (Pty) Ltd 39 Van Vuuren Street Constantia Kloof 1709 Phone +27 11 675 6104

United Kingdom

Bernafon UK Cadzow Industrial Estate Off Low Waters Road Hamilton ML3 7QE Scotland Phone +44 1698 285 968 info@bernafon.co.uk



