



- Audiometry
- REM/Live  
Speech Mapping
- Hearing Instrument Testing  
& Evaluation
- Tinnitus Assessment
- Video Otoscopy
- Hearing Instrument Vacuum



[www.medrx-usa.com](http://www.medrx-usa.com)

**MedRx**<sup>®</sup>  
*Precision Audiometric Instruments*

**MedRx Gives You Choices.**  
Choose the fit that's right for  
your practice.

**The AVANT™ line**

**of products consists of:**

- 3 Dual Channel Audiometers  
ARC, Stealth and A2D+
- 3 REM/LSM Systems  
ARC, REM Speech+ and REMsp
- 1 Hearing Instrument Test Chamber

**MedRx also offers:**

- 2 Video Oscopes
- 1 Tinnitus Assessment System
- 3 Tympanometers
- 1 UltraVac

**MedRx**

[www.medrx-usa.com](http://www.medrx-usa.com)



**MedRx, Inc**

1200 Starkey Rd. #105

Largo, FL 33771 U.S.A.

888-392-1234

727-584-9600 • FAX 727-584-9602

Email: [sales@medrx-usa.com](mailto:sales@medrx-usa.com)

EC REP

MedRx's Authorized Representative in Europe

(Regulatory affairs only)

Emergo Europe, Prinsessegracht 20

2514 AP The Hague, The Netherlands



#0086

HID (Human Interface Device)

technology - choose any USB port

and your computer will recognize the device

consistently after initial installation



# Contents

## **New MedRx Tinnometer**

Revolutionary Tinnitus Assessment..... 1

## **AVANT ARC**

Audiometry/REM/LSM Combination ..... 2

## **AVANT Stealth Audiometer**

Clinical Audiometer..... 3

## **AVANT A2D+ Audiometer**

Diagnostic Audiometer ..... 4

## **AVANT REM Speech+**

Binaural REM/LSM System ..... 5

## **Ultimate Office+**

REM+ and A2D+ in Portable Case ..... 6

## **AVANT REMsp**

Monaural REM/LSM System ..... 7

## **Otowave Tymps**

Wireless, Hand-Held Tympanometers ..... 8

## **AVANT HIT+**

Hearing Instrument Test Chamber..... 9

## **USB Video Oscopes**

External or LED Light Source ..... 10


## **UltraVac+**

Hearing Instrument Vacuum ..... 11

## **Technical Data**


Product Specs..... 12





Confidently track your patient's tinnitus  
with tools designed specifically for tinnitus

# MedRx Tinnometer



## Tinnitus Features

- Control the Level, Shape & Frequency
- Customized Stimulus and Reports
- Test Time Under 10 Minutes
- Save and Recall Sessions
- Quickly Pitch Match
- Testing Up To 16 kHz
- NOAH™, TIMS® and Sycle.Net™ Compatible

## Audiometer Features

- Air Threshold Testing
- Automatic Air Audiometry
- 125 Hz - 8 kHz
- USB Connection to Computer
- Small Footprint -  
Approx. 5" x 5" x 1.25"

## Revolutionary Tinnitus Assessment

The MedRx Tinnometer provides a whole new approach to tinnitus assessment. Confidently track your patient's tinnitus with tools designed specifically for tinnitus. Add recurring revenue with yearly tinnitus assessments. Track changes in tinnitus easily with NOAH™ sessions. Generate customized reports specific to tinnitus assessments meeting Medicare requirements.

### Air Audiometry

The MedRx Tinnometer has air threshold capabilities allowing additional versatility. Use this product for community screenings where air thresholds are needed.

Call Today! **888-392-1234**



Fully Test, Fit & Counsel Patients  
with One Sleek, Lightweight Device

# AVANT ARC



Combining the Power  
of Audiometry,  
Real Ear  
Measurement  
& Live Speech  
Mapping in One  
Compact System

The AVANT ARC offers pure tone audiometry via earphones and bone conduction, masking and speech audiometry with SRT (Speech Recognition Threshold), WR (Word Recognition), SISI (Short Increment Sensitivity Index), ABLB (Alternate Binaural Loudness Balance) and Tone Decay Tests. Additional features are HLS (Hearing Loss Simulator), MHA (Master Hearing Aid), QuickSIN™ testing and automated audiometry.

The AVANT ARC features Live Speech Mapping plus all traditional Real Ear Measurements and includes targets for MSS (Modified Speech Spectrum), DSL v5.0, NAL-NL1 and NAL-NL2 as well as HLS (Hearing Loss Simulator) and MHA (Master Hearing Aid) Modules.

## Features

- Complete Air, Bone, Speech and Masking Audiometry
- Binaural Real Ear Measurement and Live Speech Mapping
- Powerful 3rd Party Counseling Tools
- Built-In Special Tests, Word Lists and Auto-Scoring
- Integrated QuickSIN
- PC-based and Portable
- 2 x 20 W built-in Amplifiers
- Small Footprint - Approx. 8" x 5" x 1.25"
- USB Connection to Computer
- NOAH™, TIMS® and Sycle.Net™ Compatible



Add the Revolutionary Tinnitus Assessment Module. (See page 1)




# AVANT Stealth

## Compact and Powerful PC-Based Dual-Channel Clinical Audiometer

The AVANT Stealth clinical audiometer is a compact PC-based 2-channel audiometer, allowing multiple signal routing options utilizing cutting edge sound processing and sound generating technology.

This audiometer has an incredibly small footprint (approx. 8" x 5" x 1.25" - L x W x H) and contains 2 x 20 watt built-in amplifiers and can be upgraded to include high frequency testing up to 20,000 Hz.

The AVANT Stealth is a powerful PC-based audiometer that allows fast and accurate air, bone and speech testing, has dedicated transducer outputs and offers an intuitive user interface for data collection, patient monitoring and counseling.



### Features

- Dual-Channel Clinical Audiometer with User-Selectable Signal Routing
- Complete Air, Bone, Speech and Masking Audiometry
- Built-in Special Tests, Word Lists and Auto-Scoring
- Integrated QuickSIN
- Automated Audiometry
- HLS (Hearing Loss Simulator) & MHA (Master Hearing Aid) for 3rd Party Demonstration
- PC- based and Portable
- USB Connection to Computer
- NOAH™, TIMS® and Sycle.Net™ Compatible
- High frequency option allows testing up to 20 kHz




*Add the Revolutionary Tinnitus Assessment Module. (See page 1)*

**Call Today! 888-392-1234**



The AVANT Audiometer software can run stand-alone or within NOAH



# AVANT A2D+



## Compact PC-Based Dual Channel Diagnostic Audiometer

The AVANT A2D+ is a Dual Channel Diagnostic audiometer including Air, Bone, Speech and Masking functions.

It offers pure tone audiometry via earphones and bone conduction, masking and speech audiometry with SRT (Speech Recognition Threshold), WR (Word Recognition), SISI (Short Increment Sensitivity Index), ABLB (Alternate Binaural Loudness Balance), Tone Decay Tests, QuickSIN and Automated Audiometry.

The product is compact, Approx 6.5" x 5" x 1.25" (L x W x H), and when combined with a laptop is portable and easily configured for any office layout.

### Features

- Dual Channel Diagnostic Audiometer
- Complete Air, Bone, Speech and Masking Audiometry
- Built-in Special Tests, Word Lists and Auto-Scoring
- Integrated QuickSIN
- Automated Audiometry
- HLS (Hearing Loss Simulator) & MHA (Master Hearing Aid) for 3rd Party Demonstration
- PC- based and Portable
- USB Connection to Computer
- NOAH™, TIMS® and Sycle.Net™ Compatible



Add the Revolutionary Tinnitus Assessment Module. (See page1)



# AVANT REM Speech+



## Complete PC-Based Real Ear Measurement and Live Speech Mapping System



The AVANT REM Speech+ represents a new era of precision in-situ verification of hearing aids on both ears simultaneously.

In addition to Binaural Live Speech Mapping, the REM software supports all traditional Real Ear Measurements

The REM Speech+ system is cost effective, easy to connect and utilizes today's proven technology for maximum performance and accuracy.

### Features

- Binaural Real Ear Measurement and Live Speech Mapping System
- The AVANT REM Software includes targets for MSS (Modified Speech Spectrum), DSL v5.0, NAL-NL1 and NAL-NL2
- HLS (Hearing Loss Simulator) & MHA (Master Hearing Aid) for 3rd Party Demonstration
- 2 x 20 W built-in amplifiers
- PC-based and Portable
- USB Connection to Computer
- NOAH™, TIMS® and Sycle.Net™ Compatible

**Call Today! 888-392-1234**

# Ultimate Office+



*The REM+ and the A2D+ Audiometer in one convenient custom carrying case*



## Ultimate Office+ Package

The AVANT Ultimate Office combines the portability of the AVANT A2D+ air, bone and speech diagnostic audiometer with the AVANT REM Speech+ REAL EAR Measurement and Live Speech Mapping system in one custom carrying case.

Storage for all devices and accessories makes this “portable office” the preferred choice of professionals who conduct off-site testing.

Complete testing and verification in one compact case makes set-up quick, efficient and organized.



## The Complete Audiological Suite






Plug In To  
Live Speech Mapping

# AVANT REMsp

## The smallest full Real Ear Measurement and Live Speech Mapping System

The REMsp performs all of the traditional  
Real Ear Measurements and features  
monaural Live Speech Mapping for an  
affordable price.



### Features

- Monaural Real Ear Measurement and Live Speech Mapping System
- Extremely compact design  
Approx. 3.25" x 1" x .75"
- The AVANT REM Software includes targets for MSS (Modified Speech Spectrum), DSL v5.0, NAL-NL1 and NAL-NL2
- HLS (Hearing Loss Simulator) & MHA (Master Hearing Aid) for 3rd Party Demonstration
- PC-based and Portable
- USB Connection to Computer
- NOAH™, TIMS® and Sycle.Net™ Compatible



Call Today! **888-392-1234**



# Otowave Tymps

The MedRx Otowave 102, 202 & 202H offer a comprehensive range of screening and diagnostic Tympanometers

## Features

- Fast, accurate middle ear measurements
- Programmable 4-Frequency reflex test
- PC Interface; NOAH™ and MedRx AVANT Impedance module
- Ergonomic design lightweight & portable
- Data Transfer to a PC via the MedRx software module and IR dongle

## Otowave 102 Screening Tympanometry

This wireless hand-held Tympanometer with large screen provides 226 Hz probe tone impedance measurement including 4 acoustic reflex frequencies. Results are stored internally (up to 30 records). All data can be transferred to a PC via the MedRx software module and IR dongle.

The standard **Otowave 202** provides 226 Hz probe tone impedance measurements together with a user programmable range of both ipsi and contra lateral reflex test measurements at 500 Hz, 1 kHz, 2 kHz and 4 kHz.

The **Otowave 202-H** option features a comprehensive range of test functions including user defined 226 Hz and 1 kHz probe tone impedance measurements and a user programmable range of both ipsi and contra lateral reflex test measurements at 500 Hz, 1 kHz, 2 kHz and 4 kHz. A Scalar mode, Vector mode and Component mode display is available when using the 1 kHz Probe Tone.





# AVANT HIT+

## Hearing Instrument Test Chamber

The AVANT HIT+ represents a new era of precision Hearing Instrument Testing for your office. Compact yet rugged, this PC-based system is USB powered and performs 10 automated, selectable ANSI & IEC Tests. The device is small enough to be discreetly situated in any office setting. Its modern design complements the full line of AVANT instrumentation. The test results are automatically stored in NOAH™ for convenient retrieval and can be printed on any printer.

Automated, selectable ANSI & IEC tests: OSPL-90; Full-On Gain; Reference Test Gain; Frequency Response; Equivalent Input Noise; Harmonic Distortion; Battery Current; Input/Output (AGC); Attack/Release; Induction Coil (SPLITS).

### Features

- Small footprint (Approx 9.5" x 9.8" x 6")
- 10 standard ANSI and IEC Hearing Instrument Tests
- Coupler & Reference Microphones, 2cc Coupler with Adapters, Set of 4 Battery Pills
- Customizable User/Test Protocols
- HID device - True Plug and Play
- USB connection to computer
- Modern design
- NOAH™, TIMS® and Sycle.Net™ Compatible



Call Today! **888-392-1234**



Store pictures in NOAH, retrieve results easily for followup visits showing the patients their progress.

# Video Otoscopes

MedRx Video Otoscopes Consist of: Otoscope Probe, High Resolution Color Camera, LED Light Handle or an External Light Source and Video Otoscope Software. User Friendly, NOAH™ & TIMS® Compatible. An Integrated Video Module is Standard in all MedRx Software.

Choose from two MedRx Video Otoscopes with a patented lens/probe design that present ultra-clear, crisp images and require no focusing.

Choose from two MedRx Light Sources. The powerful external Standard Light Source or the built-in LED Light Handle.

**Analog Camera with Probe** - Plug the Video Otoscope straight into a monitor or connect to your computer via a USB Capture Dongle.



The image shows the analog camera with probe and a USB capture dongle. The camera has a long, thin probe and a black handle. The USB capture dongle is a small black box with a USB port and a video output port. Labels 'BACK' and 'USB Capture Dongle' are present.

**The external (SLS) Standard Light Source**

The SLS utilizes a 21 volt, 150 watt halogen bulb, fiber optic cable and a variable light intensity dial.



The image shows the external standard light source, which is a white rectangular box with a black fiber optic cable and a variable light intensity dial. A label 'MedRx' is visible on the box.

**USB Camera with Probe** - One USB Cable Connection To Your Computer. Capture, View and Store Digital Images.



The image shows the USB camera with probe and a USB cable. The camera has a long, thin probe and a black handle. The USB cable is a standard black cable. Labels 'BACK' and 'USB Cable' are present.

**The built-in LED Light Handle**

The Video Otoscope Handle produces virtually no heat and comes with two sets of rechargeable batteries for portable use.



The image shows the built-in LED light handle, which is a black handle with a long, thin probe. A label 'MedRx' is visible on the handle.

The MedRx Video Otoscope sets the clinical standard for precise imagery and archiving.

**DIGITAL  
DISPLAY**

*New separate  
pressure and  
vacuum wands  
minimize debris  
in system*

# UltraVac+

## Small, Powerful & Effective Hearing Instrument Vacuum

The UltraVac+ was designed to give dispensers and technicians a functional tool for repairing, restoring and maintaining hearing aids.

The device is small and effective in removing cerumen and debris from hearing aid tubing and ports. It offers vacuum and positive air flow simultaneously to get even the most stubborn wax out of hearing aids.

### Features

- Digital display of chamber pressure and time remaining
- Automatic detection of operating mode
- Separate pressure and vacuum wands minimize debris in system
- In-line vacuum filter for added reliability
- 4 minute drying chamber
- Auto shut off

*The UltraVac+ System  
comes standard with  
6 Assorted Vacuum Wand Tips  
plus a 4 Piece Toolkit*



**Call Today! 888-392-1234**



# MedRx Tinnometer



## Technical Specifications

**Standards:** Screening Audiometer as per ANSI S3.6-2010, IEC 60645-1:2012, Type 4; Tone Audiometry, Tinnitus Assessment

**Outputs:** Insert Earphones, DD45 Or DD450

**Frequency Range:** Air: 125 Hz – 8 kHz

**Tinnometer Frequency Range:** Air: 125 Hz – 16 kHz

**Max Output:** Air Conduction: 100 dB HL For Mid-Range Frequencies,

**Attenuation:** 1 dB Step Or 5 dB Step, User Selectable

**Optional Accessories:** DD45 Headphones & Insert Earphones

**Compatible with:** NOAH™, TIMS® and SycleNet™

**Power Requirements:** USB-powered

**Dimensions:** Approx. 5" x 5" x 1" (L x W x H)  
Approx 12 cm x 12 cm x 3 cm (W x D x H)

**Net Weight:** < 1 lbs • < 500 g

**Standard Accessories:** DD450, Patient Response Switch, Talkback Microphone, Operator Mic, Monitor Headset, Software & Manuals, Carrying Case



# AVANT ARC

## Technical Specifications

### REAL EAR MEASUREMENT

The Device Meets or exceeds all tests required in the ANSI S3.46-1997 Methods of Measurement of Real-Ear Performance Characteristics of Hearing Aids, along with the requirements of IEC/EN 61669:2001.

<b>Probe Microphones (L/R):</b>	Dual Electret Microphone Elements
<b>Probe Microphone Tube:</b>	Silicone 1.0 mm Nominal Diameter
<b>Measurement Range:</b>	40-120 ± 3 dB SPL
<b>Measured Frequency Range:</b>	125-8000Hz
<b>Test Stimuli:</b>	Broadband Noise and Synthesized Random Noise - Pink, White, Byrne LTASS and ANSI weighted; ICRA; ISTS Microphone, File, CD-ROM for Live Speech Mapping, Chirp
<b>Test Stimulus Levels at 1m:</b>	40-90 dB SPL in 1 dB Steps – 200Hz through 8K Hz (depending on speaker wattage and efficiency)
<b>Test Stimulus Accuracy:</b>	± 3dB SPL
<b>Equalization:</b>	Pressure Method
<b>Analysis Mode:</b>	User Selectable 1/3, 1/6, 1/9, 1/12, 1/24, 1/48 Octave Bands
<b>ANSI S3.46-1997 Test Available IEC/EN 61669:2001:</b>	Real Ear Unaided Response, Real Ear Unaided Gain, Real Ear Insertion Gain, Real Ear Occluded Response, Real Ear Occluded Gain, Real Ear Aided Response, Real Ear Aided Gain
<b>Other Test Available:</b>	Live Speech Mapping™ with Peaks and LTAS analysis; Real Ear to Coupler Difference, Occlusion Effect, Percentile Analysis
<b>Prescription Methods:</b>	NAL-RP, 1/3 Gain, 1/2 Gain, Berger, Pogo 1, Pogo 2, FIG6, DSL m[I/O], NAL-NL1, NAL-NL2
<b>Probe Monitoring:</b>	Available with Operator Headset

### REM EXTERNAL CONNECTIONS

<b>Power connection:</b>	USB 2.0 input 5.0 Volt Bus
<b>USB 2.0 input:</b>	Standard USB "B" socket
<b>Line-Output jack (REM or Audiometry Speakers):</b>	3.5mm Stereo Jack
<b>Speaker Output (Internal Amplifier) (2):</b>	3.81mm Pluggable Spring Clamp
<b>Probe Microphones inputs (2):</b>	8 pin Mini-DIN
<b>Operator Headset Jack (REM or Audiometry):</b>	3.5mm Stereo Jack
<b>Patient Headset Jack (Client):</b>	3.5mm Stereo Jack
<b>Power Jack:</b>	2.1mm X 5.5mm

### HEARING LOSS SIMULATOR AND HEARING AID SIMULATOR

Software based sound equalization with available Live Speech Mapping functionality. Frequency Range 125Hz – 8000 Hz, 13 Band Equalizer

### AUDIOMETRY

**Standards:** Clinical Audiometer as per ANSI S3.6-2010 Type 2 AE (IEC 60645-1 & 2), Tone Audiometry, Speech Audiometry, Stenger Test, QuickSIN™, ABLB, SISI, Tone Decay, Automated Audiometry

<b>Channels:</b>	Two channels
<b>Outputs:</b>	IP30 Insert Earphones, EAR 3A® Insert Earphones or TDH 39 Headphones (DD45), B71 Bone Conductor, Free Field-Line Level Output or Internal Amplifier
<b>Tone Stimuli:</b>	Pure Tone, Warble Tone, Continuous or Pulsed, Warble modulation frequency and Pulse period are user adjustable.
<b>Masking Signals:</b>	Tone Audiometry: Narrow Band Noise (default), Speech Weighted Noise, White Noise. Speech Audiometry: Speech Weighted Noise (default), White Noise, External Recorded (opposite channel).
<b>Frequency Range</b>	
<b>USB Power only:</b>	Air: 125Hz – 8000Hz (limited 8000Hz to 12500Hz available) Bone: 250Hz – 8000Hz
<b>Sound Field:</b>	125Hz – 8000Hz (Line Level)
<b>Acoustic Distortion:</b>	< 1.0% at 500 Hz, 100dB SPL
<b>Noise Floor:</b>	< -10dB HL from 125 Hz – 8000 Hz (12500 Hz)
<b>Attenuation:</b>	1dB or 5dB steps, user selectable
<b>Minimum / Maximum Output:</b>	-10 dB to 120 dB HL at 1 KHz – Air (¼ inch mono jacks), -10 dB to 75 dB HL at 1 KHz – Bone (¼ inch mono jack)
<b>Free Field Output:</b>	Frequency Range 125-8,000 Hz, Dynamic Range 60-90+ dB SPL at 1 meter distance, (Using 50 watt stereo amplifier with 89 dB sensitivity speakers)
<b>Speech Input:</b>	Microphone (3.5 mm stereo jacks)
<b>I/O Jacks - 3.5mm:</b>	Operator Headphones (output shared with REM), Operator Talk Forward Microphone, Patient Talk Back Microphone, Free Field (Line Out shared with REM)
<b>I/O Jacks – 1/4":</b>	Left Air Conduction, Right Air Conduction, Bone Conduction, Patient Response Switch

### POWER (FOR BOTH REM AND AUDIOMETRY)

<b>USB 2.0 Input:</b>	5.0 Volt Bus
<b>Max Power Consumption:</b>	Less than 500 mA at 5.0 volts
<b>Power Supply Internal Speaker Amp:</b>	15V DC, 2A
<b>Optional Powered Speakers:</b>	120V, 60 Hz or 100V – 240V, 50/60 Hz available
<b>Dimensions:</b>	Approx 8" x 5" x 1.25" (L x W x H) Approx 20cm x 12cm x 3cm (L x W x H)
<b>Weight:</b>	< 2 lbs • < 1kg

Call Today! **888-392-1234**



## Technical Specifications

<b>Standards:</b>	2-Channel Clinical Audiometer as per ANSI S3.6-2010, IEC 60645-1:2012, IEC 60645-2:1993, IEC 60645-4:1994, Type 1 HFAE; Tone Audiometry, Speech Audiometry, Stenger Test, QuickSIN™, ABLB, SISI, Tone Decay, Automated Audiometry
<b>Options:</b>	Automated Audiometry, High Frequency Audiometry
<b>Outputs:</b>	Insert Earphones, TDH 39, DD45 Or HDA 300 Headphones, Bone Conductor, Free Field via High Power Internal Amplifiers, 2x20 Watts Into 4 Ohms
<b>Frequency Range:</b>	Air: 125 Hz - 8000 Hz, Bone: 250 Hz - 8000 Hz, Optional: High Frequency Range With Sennheiser HDA 300 Headphones: 8000 Hz - 20,000 Hz
<b>Maximum Output:</b>	Air Conduction: 120 dB HL For Mid-Range Frequencies, Bone Conduction: 70 dB HL, Sound Field: 95 dB HL (depends on speakers)
<b>Attenuation:</b>	1 dB Step Or 5 dB Step, User Selectable
<b>Speech Input:</b>	Live Microphone, MP3/Wave Files, CD
<b>Communication Port:</b>	USB 2.0 (Backward Compatible With 1.1)
<b>Masking Signals</b>	
<b>Tone Audiometry:</b>	Narrow Band Noise (default), Speech Weighted Noise, White Noise
<b>Speech Audiometry:</b>	Speech Weighted Noise (default), White Noise, CD/File, Opposite Channel
<b>Hearing Loss Simulator and Hearing Instrument Simulator:</b>	Frequency Range: 125 Hz - 8000 Hz, 13 Band Equalizer.
<b>Optional Accessories:</b>	TDH 39 or DD45 Headphones & DD450 (High Freq. Headphones)
<b>Compatible:</b>	NOAH™, TIMS® and Sycle.Net™
<b>Power Requirements:</b>	USB-powered or External Power DC 15 V/2A
<b>Power Supply:</b>	100V - 240V, 50/60 Hz
<b>Dimensions:</b>	Approx 8" x 5" x 1.25" (L x W x H) • Approx 20cm x 12cm x 3cm (L x W x H)
<b>Weight:</b>	< 2 lbs • < 1kg
<b>Standard Accessories:</b>	Insert Earphones or DD45 Headphones, Bone Conductor, Patient Response Switch, Talkback Microphone, Operator Mic / Monitor Headset, Speaker Outputs, Auditec Sound File License USB Cable, Software & Manuals, Carrying Case



# AVANT A2D+

## Technical Specifications

<b>Standards:</b>	ANSI S3.6-2010, Type 2 AE (IEC 60645-1&2) Tone Audiometry, Speech Audiometry, Stenger Test, QuickSIN™, ABLB, SISI, Tone Decay, Automated Audiometry
<b>Channels:</b>	Two Channels
<b>Outputs:</b>	Insert Earphones or TDH39 Headphones (DD45). B71 Bone, Conductor, Free Field - Line Level Output
<b>Tone Stimuli:</b>	Pure Tone, Warble Tone, Continuous Or Pulsed. Warble Modulation Frequency And Pulse Period Are User Adjustable
<b>Masking Signals:</b>	Tone Audiometry: Narrow Band Noise (Default), Speech Weighted Noise, White Noise. Speech Audiometry: Speech Weighted Noise (Default), White Noise, External Recorded (Opposite Channel)
<b>Frequency Range USB Power Only:</b>	Air: 125Hz – 8000 Hz (limited 8000 Hz to 12500 Hz available) Bone: 250Hz – 8000Hz Sound Field: 125 Hz - 8000 Hz (Line Level)
<b>Acoustic Distortion:</b>	< 1.0% At 500 Hz, 100dB SPL
<b>Noise Floor:</b>	< -10dB HL From 125 Hz-8000 Hz
<b>Attenuation:</b>	1dB Or 5dB Steps, User Selectable
<b>Minimum / Maximum Output:</b>	-10 dB To 120 dB HL At 1 KHz – Air (¼ Inch Mono Jacks), -10 dB To 75 dB HL At 1 KHz – Bone (¼ Inch Mono Jack)
<b>Free Field Output:</b>	Frequency Range 125 Hz - 8000 Hz, Dynamic Range 60-90+ dB SPL At 1 Meter Distance, (Using 50 Watt Stereo Amplifier With 89 dB Sensitivity Speakers)
<b>Speech Input:</b>	Microphone (3.5 mm Stereo Jacks)
<b>I/O Jacks - 3.5mm:</b>	Operator Headphones (Output), Operator Talk Forward Microphone, Patient Talk Back Microphone, Free Field (Line Out)
<b>I/O Jacks – 1/4”:</b>	Left Air Conduction (2), Right Air Conduction (2), Bone Conduction, Patient Response Switch
<b>Communication Port:</b>	USB (Provides All Device Power)
<b>Power Requirements:</b>	USB Power +5 Volts DC, Less Than 500mA
<b>Dimensions:</b>	Approx 6.5” x 5” x 1.25” (L x W x H) • Approx 16cm x 12cm x 3cm (L x W x H)
<b>Weight:</b>	< 1 lb • < 500 g
<b>Standard Accessories:</b>	Insert Earphones or DD45 Headphones, Bone Conductor, Operator Mic / Monitor Headset, Patient Response Switch, Talkback Microphone, Auditec Sound File License, USB Cable, Software & Manuals, Carrying Case

Call Today! **888-392-1234**

# AVANT REM Speech+



## Technical Specifications

<b>Probe Microphones (L/R):</b>	Dual Electret Microphone Elements (2 Probe Microphones)
<b>Probe Microphone Tube:</b>	Silicone 1.0 mm Nominal Diameter
<b>Measurement Range:</b>	40 - 120 $\pm$ 3 dB SPL
<b>Measurement Frequency Range:</b>	125-8000Hz
<b>Test Stimuli:</b>	Broadband Noise and Synthesized Random Noise - Pink, White, Byrne LTASS and ANSI weighted; ICRA; ISTS; Microphone, File, CD-ROM for Live Speech Mapping, Chirp
<b>Test Stimulus Levels at 1m:</b>	45 - 90 dB SPL in 1 dB Steps (depending on speaker wattage & efficiency)
<b>Test Stimulus Accuracy:</b>	$\pm$ 3 dB SPL
<b>Equalization:</b>	Pressure Method
<b>Analysis Mode:</b>	User Selectable 1/3, 1/6, 1/12, 1/24, 1/48 Octave Bands
<b>ANSI S3.46-1997 Test Available IEC/EN 61669:2001:</b>	Real Ear Unaided Response, Real Ear Unaided Gain; Real Ear Insertion Gain; Real Ear Occluded Response; Real Ear Occluded Gain; Real Ear Aided Response; Real Ear Aided Gain
<b>Other Test Available:</b>	Live Speech Mapping with Peaks and LTASS analysis; Real Ear to Coupler Difference, Occlusion Effect, Percentile Analysis, 3D Speech Mapping and MPO Testing
<b>Prescription Methods:</b>	NAL-RP; 1/3 Gain; 1/2 Gain; Berger; Pogo 1; Pogo 2; FIG6; DSL m[l/0] NAL-NL1; NAL-NL2
<b>External Connections:</b>	Power Connection USB 2.0 Input 5.0 Volt Bus; Line Output Jack (Speakers) 3.5 mm Stereo Jack; Speaker Output (Internal Amplifier) (2) 3.81mm Pluggable Spring Clamp; Probe Microphones Inputs (2) 8 Pin Mini-DIN; Operator Headset Jack 3.5 mm Stereo Jack; Patient Headset Jack 3.5 mm Stereo Jack; Power Jack 2.1 mm X 5.5 mm.
<b>Dimensions:</b>	Approx 6.5" x 5" x 1.25" (L x W x H) • Approx 16cm x 12cm x 3cm (L x W x H)
<b>Weight:</b>	< 1 lb • <500 g
<b>Standard Accessories:</b>	Sure -Probe™ Microphone System with Lighted Visual Cues (2 Probe Microphones), 2 Headphones, 1 Powered Speaker, Carrying Case, USB Cable, Probe Tubes, Protective Probe Microphone Storage Box, External Power Supply for Internal Stereo Speaker Amplifier. Optional RECD Coupler





# AVANT REMsp

## Technical Specifications

<b>Probe Microphones (L/R) :</b>	Dual Electret Microphone Elements (One Probe Microphones)
<b>Probe Microphone Tube :</b>	Silicone 1.0 mm Nominal Diameter
<b>Measurement Range :</b>	45 - 110 $\pm$ 3 dB SPL
<b>Measurement Frequency Range:</b>	125 - 8000Hz
<b>Test Stimuli :</b>	Broadband Noise and Synthesized Random Noise - Pink, White, Byrne LTASS and ANSI weighted; ICRA; ISTS; Microphone, File, CD-ROM for Live Speech Mapping, Chirp
<b>Test Stimulus Levels at 1m:</b>	45 - 90 dB SPL in 1 dB Steps (depending on speaker wattage & efficiency)
<b>Test Stimulus Accuracy:</b>	$\pm$ 3 dB SPL
<b>Equalization :</b>	Pressure Method
<b>Analysis Mode :</b>	User Selectable 1/3, 1/6, 1/12, 1/24, 1/48 Octave Bands
<b>ANSI S3.46-1997 Test Available IEC/EN 61669:2001:</b>	Real Ear Unaided Response, Real Ear Unaided Gain; Real Ear Insertion Gain; Real Ear Occluded Response; Real Ear Occluded Gain; Real Ear Aided Response; Real Ear Aided Gain
<b>Other Test Available :</b>	Live Speech Mapping with Peaks and LTASS analysis; Real Ear to Coupler Difference, Occlusion Effect, Percentile Analysis, 3D Speech Mapping and MPO Testing
<b>Prescription Methods :</b>	NAL-RP; 1/3 Gain; 1/2 Gain; Berger; Pogo 1; Pogo 2; FIG6; DSL m[I/O] NAL-NL1; NAL-NL2
<b>External Connections :</b>	Power Connection USB 2.0 Input 5.0 Volt Bus $\pm$ 0.2 Volts, "A" Plug; Line-Output Jack (Speakers) 3.5 mm Stereo Jack (0.15 AC Volts RMS, Min. at 1 kHz); Probe Microphone Permanently Connected.
<b>Dimensions :</b>	Approx 3.25" x 1" x .75" (L x W x H) • Approx 8.3 cm x 2.5 cm x 2.0 cm (L x W x H)
<b>Weight:</b>	< 4 oz • < 120 g
<b>Standard Accessories:</b>	Sure-Probe™ Microphone with Lighted Visual Cue and Adjustable Loop, 1 Powered Speaker, Carrying Case, USB Extender Cable, Probe Tubes, Protective Probe Microphone Storage Box, Optional RECD Coupler

Call Today! **888-392-1234**



## Technical Specifications

<b>Otowane Tymp 102-4 Tympanometry Measurements:</b>	Probe Tone Levels: 226 Hz $\pm 2\%$ , 85 dB SPL $\pm 2$ dB, over range 0.2 ml to 5 ml
<b>Pressure Range:</b>	+200 daPa to -400 daPa $\pm 10$ daPa
<b>Ear Volume Measurement Range:</b>	0.2 ml to 5 ml $\pm 0.01$ ml or 10 % (whichever is larger) over entire range
<b>Reflex Measurement Range Reflex DB Range:</b>	500 Hz, 1 kHz, 2 kHz, 4 kHz Frequency $\pm 0.1$ % 85 to 100 dBHL (programmable in 5 or 10 dB steps)
<b>Reflex Measurement Range:</b>	0.01 ml to 0.5 ml $\pm 0.01$ ml configurable in 0.01 ml steps
<b>Standards:</b>	Safety IEC 60601-1; EMC IEC 60601-1-2 Impedance IEC 60645-5 Type 2 Tympanometer ANSI 3.39 CE Mark Complies to EU Medical Device Directive
<b>Power:</b>	Battery 4 x Alkaline AA cells or 4 x NiMH ( $\leq 2.3$ Ah) rechargeable
<b>Dimensions:</b>	Approx. 8" x 3" x 1" (L x W x H) • Approx. 21 cm x 8 cm x 4 cm (L x W x H)
<b>Weight:</b>	Approx < 1 lb • < 380 Gram
<b>Standard Accessories 102-4:</b>	Rugged Transportable Carrying Case, Test Cavities, Set of Disposable Eartips, 4 x 1.5V Rechargeable NiMH Batteries, Battery Charging System, Operator's Manual. Optional Accessories: Wireless IR Printer, Infra-red USB Adapter
<b>Otowane Tymp 202-202H Tympanometry Measurements:</b>	Probe Tone Levels: 226 Hz $\pm 2\%$ , 85 dB SPL $\pm 2$ dB 1000 Hz $\pm 2\%$ ; 79 dB SPL $\pm 2$ dB (202-H only) over ear canal volume range
<b>Pressure Range:</b>	+200 daPa to -400 daPa $\pm 10$ daPa or $\pm 10\%$ (whichever is larger) over range 0.1ml to 5ml
<b>Ear Volume Measurement Range:</b>	226 Hz: 0.2 ml to 5 ml ; 1000 Hz: 0.1 ml to 5ml $\pm 0.1$ ml or $\pm 5\%$ (whichever is larger)
<b>Reflex Measurement Range Reflex DB Range:</b>	Ipsilateral and contralateral or both User-selectable from: 500 Hz, 1 kHz, 2 kHz and 4 kHz ( $\pm 2$ %) Reflex levels: Ipsilateral: 70 dBHL to 110 dBHL ( $\pm 3$ dB) Contralateral: 70 dBHL to 110 dBHL ( $\pm 3$ dB)
<b>Reflex Measurement Range:</b>	0.01 ml to 0.5 ml $\pm 0.01$ ml configurable in 0.01ml steps
<b>Standards:</b>	Safety IEC 60601-1; EMC IEC 60601-1-2 Impedance IEC 60645-5 Type 2 Tympanometer ANSI S3.39 CE Mark Complies to EU Medical Device Directive
<b>Power:</b>	Mains: 100-240 VAC; 50/60 Hz via mains adapter (approved in medical safety standards); Batteries: 4 x AA (either Alkaline or NiMH, the latter recharged external to the instrument )
<b>Dimensions:</b>	Base unit: Approx. 8" x 3" x 1" (L x W x H) • Approx 19 x 8.5 x 4 cm (excluding connections) Probe: Approx. 5" x 1" • Approx 13 x 2.5 cm (max.) diameter
<b>Weight:</b>	Base unit: 11.6 oz • 330 g (without batteries, using mains power), 15 oz • 430 g (with batteries), Probe: 3.8 oz • 110 g (incl. connecting cable)
<b>Standard Accessories 202-202H:</b>	Test cavities, 4 x 1.5 V Alkaline AA batteries, Contralateral transducer, Set of disposable ear tips, Carrying Case, MedRx NOAH impedance module, USB cable, Mains adapter, Operating Manual



# AVANT HIT+

## Technical Specifications

**Standards:** Hearing Aid Analyzer As Per ANSI S3.22-2009, IEC 60118-7:2005

**Speaker Output:** Max 95 dB SPL

**Frequency Range:** 125 Hz - 8000 Hz  $\pm$  0.5%

**Coupler Mic:** Max Input Level: 140 dB SPL

**Reference Mic:** Omni Directional Microphone

**Battery Simulator Output:** Adjustable Output Voltage: 0.1V - 1.5V In 100mV Steps  
Accuracy  $\pm$  5% Current Measurement: 20uA - 20mA Accuracy  $\pm$  5%  
Battery Pills Provided: 10A, 13, 312 And 675

**Magnetic Loop:** 31.6mA/m Magnetic Strength, Per ANSI Standard

**Communication Port:** USB

**Power Requirements:** USB Power

**Dimensions:** Approx 9.5" x 9.8" x 6" (W x D x H)  
Approx 24 cm x 25 cm x 15c m (W x D x H)

**Net Weight:** < 7 lbs • < 4 kg

**Enclosure:** Acoustically Dampened Enclosure, 15dB Minimum Attenuation

**Environmental Requirements:** Working Temperature Range From 50°F To 95°F • 10° To 35° C

**Standard Accessories:** Battery Pill Set, 2cc Coupler with Adapters, USB Cable, Hex Wrench, O-Ring, 25mm Tube, Foam Pads, Blue Tack

**Automated, Selectable ANSI & IEC tests:** OSPL-90; Full-On Gain; Reference Test Gain; Frequency Response; Equivalent Input Noise; Harmonic Distortion; Battery Current; Input/Output (AGC); Attack/Release; Induction Coil (SPLITS)

# Video Otoscope



## Technical Specifications

<b>USB Camera</b>	
<b>Image Sensor:</b>	1/3" Interline XGA Color Progressive CCD: ICX204AK (Sony)
<b>Effective Picture Element:</b>	1024 (H) x 779 (V)
<b>Chip Size:</b>	5.80 (H) x 4.92 (V) mm
<b>Cell Size:</b>	4.65 (H) x 4.65 (V) $\mu$ m
<b>Scanning System:</b>	Progressive
<b>Resolution:</b>	1024 (H) x 768 (V) (Full Scanning) 1024 (H) x 344 (V) (1/2)
<b>Maximum Frame Rate</b>	
<b>Full Scanning:</b>	29.18 fps (Normal) / 14.59 fps (1/2 clock) / 7.295 fps (1/4 clock)
<b>Pixel Frequency:</b>	29.5 MHz (Normal) / 14.75 MHz (1/2 clock) / 7.375 MHz (1/4 clock)
<b>Video Output:</b>	8bit / 10bit / 12bit
<b>Minimum Scene Illumination</b>	
<b>21 Lux @ F1.2:</b>	
<b>Sync System:</b>	Internal
<b>Electronic Shutter:</b>	Auto / Manual (Software Selectable)
<b>Normal:</b>	1/29,500,000 ~ 1/29.18 seconds
<b>Gain:</b>	Auto / Manual (Software Selectable)
<b>Gamma:</b>	Manual (Software Selectable)
<b>White Balance:</b>	Auto / Manual / One shot (Software Selectable)
<b>Input / Output:</b>	USB 2.0 High Speed
<b>Power</b>	
<b>Input Voltage:</b>	+5 Vdc through USB connector (+4.4 ~ +5.25V)
<b>Consumption:</b>	< 300 mA
<b>Dimensions:</b>	28 mm x 28 mm x 37 mm (W x H x D)
<b>Lens Mount:</b>	CS Mount
<b>Environmental</b>	
<b>Weight:</b>	Approx. 45g
<b>Interface Connector :</b>	
	USB: mini-B USB connector IO signal: 6pin connector (HR10A-7R-6PB or equivalent)
<b>Temperature Operational:</b>	0 ~ 40°C
<b>Storage:</b>	-30 ~ 65°C
<b>RoHS:</b>	RoHS Compliant

<b>Analog Camera</b>	
<b>Image Sensor:</b>	1/3" Interline CCD ICX638BKA
<b>Active Picture Element:</b>	768 (H) x 494 (V)
<b>Signal Format:</b>	NTSC
<b>Scanning System:</b>	2:1 Interlace
<b>Scanning Frequency:</b>	Horizontal Frequency 15.734kHz, Vertical Frequency 59.94Hz
<b>Sync System:</b>	Internal / External
<b>Horizontal Resolution:</b>	480 TV Lines
<b>S/N Ratio:</b>	More than 48 dB (AGC = off)
<b>Video Output Format:</b>	VBS 1.0 Vp-p, 75 $\Omega$ , Y/C
<b>Minimum Scene Illumination:</b>	0.11 lx, F1.2 (AGC=ON)
<b>Electronic Shutter:</b>	
<b>Dip Switch:</b>	1/60 (1/50PAL), 1/125, 1/250, 1/500, 1/1,000, 1/2,000, 1/4,000, 1/10,000 usec
<b>Auto/Control Software:</b>	High Speed Shutter: 1/60(1/50:PAL), 1/125, 1/250, 1/500, 1/1000, 1/2000, 1/4000, 1/10000, 1/100000 sec
<b>Auto/Control Software:</b>	Low Speed Shutter: 1 to 256FLD (Configurable through the control software)
<b>White Balance:</b>	Auto / Push to Set/ Manual (Configurable through the control software)
<b>Auto IRIS:</b>	Not supported
<b>AGC (Automatic Gain Control):</b>	On / Off
<b>Gamma:</b>	0.45/1.0 (Switchable, Con- figurable through the control software), Default:0.45
<b>Image Rotation:</b>	Normal (Default), Horizontal Flip, Vertical Flip, Horizontal Vertical Flip
<b>Still Image:</b>	Supported
<b>Lens Mount:</b>	CS Mount
<b>Optical LPF:</b>	IR cut filter
<b>Input Voltage:</b>	DC9V ~ 15V
<b>Consumption:</b>	80mA $\pm$ 20mA
<b>Operating Temperature:</b>	-10 ~ 50°C
<b>Storage Temperature:</b>	-30 ~ 65°C
<b>Dimensions:</b>	51 (W) x 51 (H) x 60.5 (D) mm
<b>Weight:</b>	Approx. 190g
<b>RoHS:</b>	RoHS Compliant



# MedRx UltraVac+

## Technical Specifications

<b>Power Requirements:</b>	100V - 240V, 50/60 Hz
<b>Power Consumption:</b>	<100 Watts
<b>Fuse:</b>	1.25 amp SB 5 x 20 mm
<b>Operating Temperature:</b>	32° TO 120° F • 0° TO 50° C
<b>Weight:</b>	Approx 11 lbs • Approx 5 kg
<b>Dimensions:</b>	Approx 11" x 10.5" x 6.5" (L x W x H) • Approx 28 cm x 26 cm x 16 cm (L x W x H)
<b>Debris Filter:</b>	I-35S
<b>Vacuum Tips:</b>	#14, #15, 2-#18 and 2-#20 gauge with Safety Lock attachment
<b>Drying Chamber Cycle Time:</b>	4 minutes with auto shut-off
<b>Vacuum Wand:</b>	Constant
<b>Standard Accessories:</b>	Drying Chamber Cup, Particulate Filter, 6 Vacuum Wand Tips, Four Piece Tool Kit





# Innovation, Continuity, Reliability.

MedRx, Inc. is a U.S. based global manufacturer and developer of advanced computerized diagnostic and hearing instrument fitting technologies, specifically designed for the hearing care professional. MedRx has created a remarkable new generation of discreet, yet powerful PC-based instrumentation for Audiometry, Real Ear Measurement, Live Speech Mapping, Hearing Instrument Testing & Evaluation and Video Otoscopy.

## **MedRx** Precision Audiometric Instruments

1200 Starkey Rd. Suite105

Largo, FL 33771

727-584-9600

888-392-1234

Fax: 727-584-9602

Email: [sales@medrx-usa.com](mailto:sales@medrx-usa.com)

Website: [www.medrx-usa.com](http://www.medrx-usa.com)

### **MedRx Minimum Computer Specs:**

Windows® - PC compatible computer, Intel™ Dual Core, 1.8 GHz or better. 2 GB RAM. 5 GB free hard drive space. Available 2.0 USB port. Windows 7, 8 or 10 Professional (32 or 64-bit).