# natus.


# Bio-logic® ABBaer® hearing screening system

🛔 🗱 🛝 🛄 😯				LOX	
1		Baer Test Comp	plete	1	
Patent Internation		6			
2		{::	) ≈		NetCeller P
					Creat
Inder	inger 1		-	-	
	Right	halle to an I	Left	1	

The Power of Choice in Newborn Hearing Screening

The ABaer<sup>®</sup> System is an all-in-one ABR, DPOAE & TEOAE screener — providing unsurpassed flexibility in hearing screening & data management

# THE ABAER HEARING SCREENER PROVIDES FLEXIBLE SOLUTIONS FOR MEETING PROGRAM REQUIREMENTS & ACHIEVING CONTINUOUS QUALITY IMPROVEMENT GOALS



## ALL-IN-ONE ABR, DPOAE & TEOAE SCREENING SYSTEM

### ABR technology

Utilizes proprietary Point Optimized Variance Ratio (POVR) Algorithm

- Developed by the House Ear Institute, a leader in auditory research since 1946
- Efficient screening
- Test automatically stops if the probability of achieving a PASS result is very low
- Statistically proven
- 99.96 % theoretical statistical bilateral sensitivity
- 95 % specificity

# FAST & SIMPLE OPERATION

- Keyboard/mouse functions provide easy navigation
- Colored, graphical displays guide you easily through
- the screening process
- Patient ABR waveforms are available for viewing
- Automated pass/refer results
- Pre-set screening parameters
- no user adjustments necessary
- Tutorial videos provide helpful tips on how to conduct a screen

		1		CONE   P OPONE	# 18 m   3	
	Latitume		Bittan	Hatelynam	Ra C Per C La	Barthand
Test Date	Rock (Pan)		Ranto PCMR Secon K (2.31	Text Exte TU-0/205 325-52-64 0.6		100 Anifacta POM 100 130 4.25
04 03 04 03 04 03 04 04 04 04 04 04 04 04 04 04 04 04 04				0.4 0.5 0.1 0.1 0.1 0.0 0.0 0.0 0.0 0.0 0.0 0.0	n	

Waveforms available for viewing

# FLEXIBLE DATA TRANSFER & DATA MANAGEMENT

- Stores an expanded set of patient demographic data
- User-definable data entry fields are available to meet program requirements
- Provides customizable reports for parents & physicians
- Access to screening database & statistical reports can be restricted to authorized individuals via password
- Includes built-in data management software
- Compatible with 3rd party data management systems such as OZ SIMS & Hi\*Track

# Distortion Product Otoacoustic Emissions (DPOAE)


- High sensitivity for cochlear hearing loss
- Frequency specific screening: 2, 3, 4, and 5 kHz
- Fast test times
- Approximately 10 seconds per ear
- In-the-ear stimulus calibration and ear probe stability check prior to testing

## Transient Evoked Otoacoustic Emissions (TEOAE)

- High sensitivity for cochlear hearing loss
- Testing frequency range: 1.2 to 3.5 kHz
- Fast test times
- Approximately 10 seconds per ear
- Ongoing ear probe stability check throughout the test

# Multiple Technology Configurations

- ABR, DPOAE & TEOAE technologies are available in different combinations
- Ideal for one-stage, two-stage or combined screening protocols



Customizable Patient Report

- Export results to personal computer via CD or USB
- Barcode reader available for input of medical record numbers



Automated PASS Result

# BUILT-IN CONTINUOUS QUALITY IMPROVEMENT

- Generates statistical reports on program outcome measures
  - Hospitals can monitor the effectiveness of their screening program on an ongoing basis

# **CONVENIENT SCREENING WORKSTATION**

- Cart allows easy transport of ABaer System to different locations
- Spacious drawer & basket are available for storing screening supplies

a ot



#### The ABaer Hearing Screening System includes:

- ABaer screener with panel PC Insert earphones or laptop computer
- Screening cart
- Data management software (built-in)
- Electrode (patient) cable
- OAE probe (optional for ABR, included with OAE)

**Ordering Information:** 

 Seiko label printer (optional)

Deskjet printer (optional)

Supply starter kit

with adapters

Instructional materials

#### ITEM PART NUMBER ABaer, ABR, Label Printer AB1001 ABaer, ABR, Regular "PC" Printer AB1002 ABaer, ABR, Label Printer, OAE Probe AB1003 ABaer, ABR, Regular "PC" Printer, OAE Probe AB1004 ABaer, ABR, DPOAE, Label Printer AB1005 ABaer, ABR, DPOAE, Regular "PC" Printer AB1006 ABaer, ABR, TEOAE, Label Printer AB1007 ABaer, ABR, TEOAE, Regular "PC" Printer AB1008 ABaer, ABR, DP & TEOAE, Label Printer AB1009 ABaer, ABR, DP & TEOAE, Regular "PC" Printer AB1010

Natus Medical Incorporated 1501 Industrial Road San Carlos, CA 94070 USA Ph: ((847) 949-5200 Fax: (847) 949-8615

**Global Sales & Support** 1-800-303-0306

www.natus.com ©2010 Natus Medical Incorporated P/N 006570A

# **Technical Specifications:**

#### PHYSICAL

ABaer interface device: **Dimensions:** 

> 7.6 x 4.0 x 1.7 inches 19.4 x 10.1 x 4.4 cm (H x W X D)

> > 14 oz (400 grams)

Weight:

**POWER SUPPLY** 6 V DC

#### **ABR DATA ACQUISITION**

Analysis Time A/D Resolution Artifact Rejection Points Per Trace **Electrode Montage**  (window length) 21.33 msec 16 Bit >13 uV 256 Forehead to Test Ear or Forehead to Nape of Neck (dependent on transducer choice)

1 Optically Isolated

>110 dB. at 50/60 Hz

Internal, 1000 Hz sine wave

x30,000 100 Hz

1500 Hz

#### ABR AMPLIFIERS

Channels Gain **High Pass Filter** Low Pass Filter CMR Ratio Impedance Test

#### **ABR STIMULUS**

Type Polarity Rate Intensity 100 µs click alternating 37.1/sec 35 dB nHL default (option for 30 or 40 dB nHL as an alternative)

#### **ABR TRANSDUCER CHOICES**

ABaer Probe (OAE-type) Insert Earphones with in-the-ear tips Insert Earphones with Halo Ear Muffin® TDH-39 headphones

#### COMPUTER

System requires PC with dedicated USB and Windows® XP Professional operating system. Bio-logic-supplied PC specifications are available upon request.

#### CART

Five-caster cart accommodates ABaer module, laptop or panel PC, deskjet printer, isolation transformer and supplies.

#### SAFETY

Designed to meet the following standards: IEC 60601-1 Class II UL 2601-1 CSA-C22.2 No. 601.1 AAMI-ES1 **CE** Certified

Contact your local customer service representative for more information. Note: Specifications are subject to change without notice.

> Trusted Reliability. Flexible Innovation. Proven Excellence.

natus hearing diagnostics