ADAPTABLE MID-LEVEL AUDIOMETER

PELLO

INSTRUCTIONS FOR USE

Table of Contents

EN – English

CS – Czech

EL - Greek

FR - French

HU - Hungarian

IT – Italian

PL – Polish

SK - Slovak

SL - Slovenian

SV – Swedish

ADAPTABLE MID-LEVEL AUDIOMETER

PELLO

INSTRUCTIONS FOR USE

Title: GSI Pello Instructions for Use

Manufacturer


Grason-Stadler, Inc.
10395 West 70th Street
Eden Prairie, MN 55344
USA

Copyright © 2021 Grason-Stadler. All rights reserved. No part of this publication may be reproduced or transmitted in any form or by any means without the prior written permission of Grason-Stadler. The information in this publication is proprietary to Grason-Stadler.

This product may incorporate software developed by others and redistributed according to license agreements. Copyright notices and licenses for these external libraries are provided in a supplementary document that is included with the other software applications distributed with this product.

Compliance

The CE 0123 mark identifies compliance with the Medical Device Regulation (EU) 2017/745. Grason-Stadler is an ISO 13485 certified corporation.

European Authorized Representative

Grason-Stadler
c/o DGS Diagnostics A/S
Audiometer Alle 1
5500 Middelfart
Denmark

0123

Caution: US Federal law restricts this device to sale by or on the order of a physician or licensed hearing care professional.

TABLE OF CONTENTS

Table of Contents	1
Preface.....	6
Manual Conventions.....	6
Regulatory Symbols.....	7
Device Symbols	8
Warnings and Cautions.....	10
Important Safety Instructions.....	10
Safety Precautions	10
Fluids	11
Damage.....	11
Connecting Additional Equipment.....	11
EXPLOSION	12
ENVIRONMENTAL FACTORS.....	12
ELECTROMAGNETIC COMPATIBILITY (EMC) CONSIDERATIONS.....	12
Introduction	13
Indications for Use.....	13
Intended Use	13
Contraindications.....	13
Essential Performance.....	13
Unpacking the system.....	14
Standard Components	15
Getting Started.....	16
System Assembly	16
Base Unit Rear Panel	16
Connections.....	16
Right Side Panel	18
Power Switch	18
Left Side Panel.....	18

USB Port	19
A/B Cable.....	19
Front Panel	20
Monitor	20
Elimination of Ambient Noise	21
Maximum Ambient Noise (dB re: 20 μ Pa to the nearest 0.5 dB)	21
Operating the Pello	22
Preliminary Checks	22
Prepare to Test.....	22
Placement of the Earphones.....	23
Placement of the Insert Phone	23
Placement of the High Frequency Transducer (optional license & DD450 headphone required)	23
Placement of the Bone Vibrator.....	23
Loud Auditory Output.....	23
Test Type Screens	24
Title Bar	24
Test Type Information.....	24
Navigation Menu	25
Time and Date	25
Common Icons	25
Tone Test Type	26
High Frequency Testing (Optional High Frequency License).....	27
TEN Test Type (Optional Special Test License)	27
Presentation Level.....	27
Test Instructions.....	28
Test Procedure	28
Interpretation	28
Speech Test Type Button.....	28
Integrated Word Files.....	29
More Test Type button (Optional Licenses).....	30

BKB-SIN (Optional Speech Plus License).....	31
Presentation Level.....	31
Test Instructions.....	31
Test Procedure	32
QuickSIN (Optional Speech Plus License)	33
Presentation Level.....	33
Test Instructions.....	33
Test Procedure	34
Configuration and Update	35
Update.....	35
Bone	35
Print.....	36
Export Logs.....	36
License.....	36
Settings.....	36
Cleaning and Maintenance	37
Routine Maintenance	37
Biological Calibration Check	37
Periodic Checks	37
Earphone and Bone Vibrator Cords	37
Hum and Noise.....	37
Distortion and Frequency Shift	37
Speech Level Check.....	38
Internal Controls Check	38
Bone Vibrator Check	38
Masking Level Check.....	38
Talk Forward Check.....	38
Calibration	38
Cleaning the System	39
Cleaning and Disinfecting Agents	39

Single-use Accessories	39
Status and Error Messages	40
Error Messages	40
Please try another selection.....	40
No test data stored	40
Printer communication error	40
Record test result in comments.....	40
The session comments have been updated with the results of the SDT test.....	40
Not supported in speech	40
Speech data limit exceeded, speech tables limited to 6 test results per ear. Latest test result will not be saved	40
Error	40
Error	40
Supplies and Accessories.....	41
Appendix A – Technical Data	42
Appendix B – Specification of Input/output Connections	48
Mains Cable	48
Headphones/Inserts (2 connectors, R, L).....	48
Bone Oscillator.....	48
Talk Forward, Talk Back Mic.....	49
Monitor Headset.....	49
External A+B.....	49
Speakers (FF1, FF2).....	49
USB A (Host) x4	50
USB B (device)	50
LAN Ethernet	50
Appendix C - Calibration Reference & Maximum Levels.....	51
Earphones - Pure Tone RETSPL.....	52
Earphones - ANSI Speech RETSPL.....	53
Earphones - IEC Speech RETSPL	54
Earphones - Pure Tone max HL.....	55

Earphones - NB noise effective masking level.....	56
Earphones - NB noise max HL	57
Earphones - ANSI Speech max HL	58
Earphones - IEC Speech max HL.....	58
Insert Earphones - Pure Tone RETSPL.....	59
Insert Earphones - ANSI Speech RETSPL.....	60
Insert Earphones - IEC Speech RETSPL	60
Insert Earphones - Pure Tone max HL	61
Insert Earphones - NB noise effective masking level.....	62
Insert Earphones - NB noise max HL	63
Insert Earphones - ANSI Speech max HL	64
Insert Earphones - IEC Speech max HL.....	64
Bone Vibrator - Pure Tone RETFL	65
Bone Vibrator - ANSI Speech RETSPL	66
Bone Vibrator - IEC Speech RETSPL.....	66
Bone Vibrator - Pure Tone max HL.....	67
Bone Vibrator - NB noise effective masking level	68
Bone Vibrator - NB noise max HL.....	69
Bone Vibrator - ANSI Speech max HL.....	70
Bone Vibrator - IEC Speech max HL	70
Free Field Speakers – ANSI RETSPL and Max HL	71
Appendix D - Warranty and Repair	72
Appendix E – Recycling and Disposal.....	73
Appendix F – Electromagnetic Compatibility (EMC).....	74
Cautions regarding EMC.....	74
Guidance and manufacturer’s declaration Electromagnetic emissions.....	75
Recommended separation distances between portable and mobile RF communications equipment	76
Guidance and Manufacturer’s Declaration Electromagnetic Immunity.....	77

PREFACE

This user manual provides information about the GSI Pello. This manual is intended for technically qualified personnel. This manual contains information pertinent to the use of the Pello system including safety information as well as maintenance and cleaning recommendations. Proper use of this device depends on careful reading of all instructions and labeling. Also be sure to follow all safety standards set by your facility.

Please note: This User Manual is not intended as a training manual. The reader should consult standard audiology texts for the theory and application of the tests provided by this instrument.

READ THIS ENTIRE MANUAL BEFORE ATTEMPTING TO USE THIS SYSTEM!

MANUAL CONVENTIONS

Throughout this manual, the following meaning of warnings, cautions and notices are used.

WARNING

The WARNING symbol identifies conditions or practices that may present danger to the patient and/or user.


CAUTION

The CAUTION Symbol identifies conditions or practices that could result in damage to the equipment.

NOTE: Notes help you identify areas of possible confusion and avoid potential problems during system operation.

REGULATORY SYMBOLS

Symbol	Description
	Conforms to Medical Device Regulation (EU) 2017/745.
	Symbol for "SERIAL NUMBER."
	GSI Part Number.
	Return to Authorized Representative, Special disposal required.
	Medical Equipment Classified by Intertek Testing Services NA Inc. with respect to electric shock, fire, and mechanical hazards only, in accordance with UL 60601-1. Classified under the Medical Device Regulation (EU) 2017/745 as a Class IIa device.
	China RoHS symbol for product compliance
	Symbol for "European Authorized Representative."
	Symbol for "Manufacturer."
	Symbol for "Date of Manufacture."
	On/Off
	Keep Dry

Symbol	Description
	<p>This side up.</p>
	<p>Consult the operating instructions/directions for use. A copy of the operating manual is available on this website: www.grason-stadler.com A printed copy of the operating instructions can be ordered from Grason-Stadler for shipment within 7 days; or you can contact your local representative.</p>
	<p>Consult the operating instructions/directions for use. A copy of the operating manual is available on this website: www.grason-stadler.com A printed copy of the operating instructions can be ordered from Grason-Stadler for shipment within 7 days; or you can contact your local representative.</p>

DEVICE SYMBOLS

The following symbols appear on the instrument.

Definition: Type B applied part – IEC 60601-1 uses the term applied part to refer to the part of the medical device which comes into physical contact with the patient for the device to carry out its intended function. This icon indicates that the GSI Pello is Type B according to the international standard IEC 60601-1. The applied parts are the headphones, insert earphones, and the bone oscillator.

Symbol	Description
	Consult Operating Instructions.
	Follow Instructions for Use.
	Connection for headphones. The red markings around the connector indicate the right ear and the blue marking indicate the left connector.
	Connection for the insert earphones. The red markings around the connector indicate the right ear and the blue marking indicate the left connector.
	Connection for the bone conductor transducer.
	Connection for the sound field speakers. The red markings around the connector indicate the right ear and the blue marking indicate the left connector.
	Patient response switch
	On/Off - Next to power mains.

WARNINGS AND CAUTIONS

IMPORTANT SAFETY INSTRUCTIONS

The following safety precautions must always be observed. General safety precautions must be followed when operating electrical equipment. Failure to observe these precautions could result in damage to the equipment and injury to the operator or patient.

The instrument must only be used by hearing health care professionals qualified to perform hearing threshold and auditory tests such as an audiologist, otolaryngologist, researcher, or a technician under direct supervision by the specialist. Users should use their professional skills when interpreting the results and this should be done in conjunction with other testing as deemed appropriate given their professional skills. Incorrect use could lead to wrong results.

The employer should instruct each employee in the recognition and avoidance of unsafe conditions and the regulations applicable to his or her work environment to control or eliminate any hazards or other exposure to illness or injury. It is understood that safety rules within individual organizations vary. If a conflict exists between the material contained in this manual and the rules of the organization using this instrument, the more stringent rules should take precedence.

In case of death or serious incident in relation to the use of the device, the incident must immediately be reported to Grason-Stadler and the local national competent authority.

SAFETY PRECAUTIONS

This product and its components will perform reliably only when operated and maintained in accordance with the instructions contained in this manual, accompanying labels, and/or inserts. A defective product should not be used.

Make sure all connections to external accessories are snug and secured properly. Parts which may be broken or missing or are visibly worn, distorted, or contaminated should be replaced immediately with clean, genuine replacement parts manufactured by or available from GSI.

This instrument is suitable in hospital environments except for near active high frequency (HF) surgical equipment and radio frequency (RF) shielded rooms of systems for magnetic resonance imaging, where the intensity of electromagnetic disturbance is high.

In case of emergency, disconnect the instrument from the supply mains by pulling the plug out of the power socket.

Before the first use of the instrument each day, or if suspect or inconsistent results are apparent, the checks specified in the Routine Maintenance section should be carried out. If the system is not functioning properly, do not operate it until all necessary repairs are made and the unit is tested and calibrated for proper functioning in accordance with Grason-Stadler published specifications.

No parts of the equipment can be serviced or maintained while in use with the patient.

Latex is not used anywhere in the manufacturing process. The base material for the earphone cushions is made from natural and synthetic rubber.

Fluids

The device is not intended to be used in environments exposed to fluid spills. No means is specified for fluid protection (not IP classed). Do not use the device in the presence of fluid that can contact any of the electronic components or wiring. Should the user suspect fluids have contacted the system components or accessories, the unit should not be used until deemed safe by an authorized service technician. Do not immerse the unit in any fluids. See the Routine Maintenance Section of this manual for the proper cleaning procedure for the instrument and its accessories and the function of single-use parts.

Damage

Do not drop or otherwise impact this instrument. If the instrument is dropped or damaged, return it to the manufacturer for repair and/or calibration. Do not use the instrument if any damage is suspected.

Do not attempt to open, modify or service the instrument. Return the instrument to the manufacturer or distributor for all repair and servicing requirements. Opening the instrument will void the warranty. See Warranty and Repair in the Appendix

Connecting Additional Equipment

This equipment is intended to be connected to other equipment thus forming a Medical Electrical System. External equipment intended for connection to signal input, signal output or other connectors shall comply with the relevant product standard e.g., IEC 60950-1 for IT equipment and the IEC 60601-series for medical electrical equipment. In addition, all such combinations – Medical Electrical Systems – shall comply with the safety requirements stated in the general standard IEC 60601-1, edition 3, clause 16. Any equipment not complying with the leakage current requirements in IEC 60601-1 shall be kept outside the patient environment i.e., at least 1.5 m from the patient support or shall be supplied via a separation transformer to reduce the leakage currents.

Any person who connects external equipment to signal input, signal output or other connectors has formed a Medical Electrical System and is therefore responsible for the system to comply with the requirements. If in doubt, contact a qualified medical technician or your local representative. If the instrument is connected to a PC (IT equipment forming a system) ensure not to touch the patient while operating the PC. If the instrument is connected to a PC (IT equipment forming a system) assembly and modifications shall be evaluated by qualified medical technician according to safety regulations in IEC 60601.

EXPLOSION

This system is not explosion proof.

Do NOT use in the presence of flammable gaseous mixtures. Users should consider the possibility of explosions or fire when using this device near flammable anesthetic gases.

Do NOT use the system in a highly oxygen-enriched environment, such as a hyperbaric chamber, oxygen tent, etc.

ENVIRONMENTAL FACTORS

Use and store the instrument indoors only. It is recommended that the instrument be operated within an ambient temperature range of 15 °C / 59 °F to 35 °C / 95 °F and in relative humidity between 30 % and 90 % (non-condensing).

Transport the instrument in temperatures between 20 °C / -4 °F to +50 °C / +122 °F and store the instrument in temperatures between 0 °C / 32 °F to 50 °C / 122 °F.

ELECTROMAGNETIC COMPATIBILITY (EMC) CONSIDERATIONS

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information in the Appendix. Although the instrument fulfills the relevant EMC requirements, precautions should be taken to avoid unnecessary exposure to electromagnetic fields, e.g., from mobile phones etc. This appendix provides guidance on the electromagnetic environment in which to operate the instrument.

Portable and mobile radio frequency (RF) communications equipment can affect medical electrical equipment. Use of this instrument adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this instrument and the other equipment should be observed to verify that they are operating normally.

INTRODUCTION

Thank you for purchasing the GSI Pello™ system. The GSI Pello is a stand-alone, 1.5 channel, clinical audiometer. It supports data transfer to a computer but does not rely on the computer to perform tests and will not lose data if the computer connection fails. User login and password controls provide security for patient data. Test data can be transferred to software such as GSI Suite and Noah and can be integrated into your EMR/EHR system.

INDICATIONS FOR USE

The Pello is intended to be used for the identification and etiology of hearing loss in patients of any age. It is intended to be used by an audiologist, ENT, hearing healthcare professional, or trained technician in a hospital, clinic, healthcare facility or other suitable quiet environment as defined in ANSI S3.1 or equivalent.

INTENDED USE

The GSI Pello is intended to be used for the purposes of determining patient hearing sensitivity. It is intended to quantify the patient's level of hearing by presenting pure tones stimuli through specific transducers at different frequencies and at different sound pressure levels.

CONTRAINDICATIONS

There are no known contraindications for use.

ESSENTIAL PERFORMANCE

Essential Performance for this instrument as per IEC 60601-1 is defined by the manufacturer as:

This instrument does not have an Essential Performance.

Absence or loss of Essential Performance cannot lead to any unacceptable immediate risk. Final diagnosis shall always be based on clinical knowledge.

UNPACKING THE SYSTEM

Although this GSI Pello Clinical Audiometer was carefully tested, inspected, and packed for shipping, it is good practice after receiving the instrument to immediately examine the outside of the container for any signs of damage. Notify the carrier if any damage is observed.

- ✓ It is recommended that you unpack your Pello carefully to make sure that all components are removed from the packing materials.
- ✓ Verify that all components are included as shown on the packing slip included with your shipment.
- ✓ If any component is missing, contact your distributor immediately to report the shortage.
- ✓ If any component appears to be damaged in shipment, contact your distributor immediately to report it. Do not attempt to use any component or device that appears to be damaged.
- ✓ Please check that all accessories listed below have been received in good condition.
- ✓ All standard components are part of the medical device and suitable for operation in a patient environment.

Notify the carrier immediately if any mechanical damage is noted. This will ensure that a proper claim is made. Save all packing material so the claim adjuster can inspect it as well. Notify your dealer or GSI when the adjuster has completed the inspection.

Save all the original packing material and the shipping container so the instrument can be properly packed if it needs to be returned for service or calibration.

STANDARD COMPONENTS

- Pello Base Unit and power cord
- Subject Response Hand Switch
- Monitor Headset with boom mic
- IP30 Insert earphones, 10 ohm
- Bone Vibrator Assembly, B81, 10 ohm
- Wireless keyboard and mouse
- GSI Pello Software and Manuals Bundle, USB Flash Drive
- GSI Suite Software Bundle, USB Flash Drive
- Cable, USB A/B, 2 meters
- Earpads for monitor headset with boom mic
- Microphone windscreen (Sennheiser PS 01)
- GSI Cleaning Cloth
- Patch Cord, 6 feet, red, 1 each
- Patch Cord, 6 feet, blue, 1 each
- Patch Cord, 6 feet, black, 1 each
- Talk Back Microphone (EM400)

GETTING STARTED


SYSTEM ASSEMBLY

1. Connect transducers to audiometer base unit
2. Connect USB cable from base unit to PC if using GSI Suite
3. Connect power cable from base unit to outlet or isolation transformer


BASE UNIT REAR PANEL

Connections

The connectors on the rear panel of the GSI Pello are shown below. The label and jacks are visible by turning the instrument around on a flat, stable surface.

	Connection	Description	Graphic
R1	Mains Power Input	IEC 14	
R2	USB Computer Connection	USB B style connector	
R3	LAN Connections	Ethernet Connection RJ45 Currently not supported	
R4	USB Connectors	USB A style plug	

	Connection	Description	Graphic
R5	FF Speaker RCA Connections Output	2 RCA jacks	
R6	Left and Right Headphone Outputs	6.35 mm mono jacks Left (blue) and Right (red)	
R7	Left and Right Insert Phone Outputs	6.35 mm mono jacks Left (blue) and Right (red)	
R8	Bone vibrator	6.35 mm mono jack	
R9	Patient Response Inputs	6.35 mm stereo jack 1 hand switch can be used	
R10	Talkback Microphone Input	6.35 mm stereo jack NOTE: Microphone inputs are between .25 mV and 5 mV for a 0 dB reading on a VU indicator; the input impedance is 3,200 ohms	
R11	Ext. A and B	3.5 mm stereo jack Input jack for optional digital music player or CD player input NOTE: External A and B inputs are between 15 mV and 500 mV for a 0 dB reading on a VU indicator; the input impedance is 50,000 ohms	

RIGHT SIDE PANEL

Power Switch



The power switch is located on the right-side panel. To turn on/off the GSI Pello, use the power switch.

NOTE: Do not block access to the power switch.

LEFT SIDE PANEL

The following connectors are found on the left side panel of the GSI Pello.

Connection	Description	Graphic
Monitor Microphone	3.5 mm stereo jack Monitor microphone	
Monitor Headphones	3.5 mm stereo jack Monitor headphone	

USB Port

The Pello is equipped with four (4) USB A and one (1) USB B ports. It is possible to connect external devices such as mouse, keyboard, or external printer to be used with the audiometer. Additionally, a memory stick can be inserted into a USB port for updating software, adding additional sound files, license updates, printing to PDF or exporting diagnostic log files.

NOTE: Scan files on a USB drive for viruses prior to installing the drive into the instrument.


A/B Cable

Remote connection to an external computer is achieved using a standard A/B USB cable.

NOTE: It is recommended to always have the USB ports enabled on the PC. Disable the “suspend USB” option on the PC.

FRONT PANEL

The front panel controls are used to operate the device. The blue shaded areas have controls that operate the 2 channels. Channel 1 is used for the stimulus and channel 2 for masking. The areas outside the blue shaded area control the functions of the device that are channel independent. The controls are labeled to indicate their function. The four navigation buttons, and the middle select button, can be used to make selections from the on-screen menus as well as navigate through the internal files for speech testing.

MONITOR

The Pello comes standard with an LCD display. The LCD is hinged to the GSI Pello and is used to display all the testing information from the instrument. The information displayed on the Pello LCD varies depending on the Test Type.

ELIMINATION OF AMBIENT NOISE

Excessive noise in the test environment reduces test validity because it masks the test signals. Maximum permissible ambient noise levels (MPANLs) are defined in the standard ANSI S3.1-1991 (R2008) American National Standard Maximum Permissible Ambient Noise Levels for Audiometric Test Rooms.

The following table shows the maximum background levels that can be present inside the room while a valid hearing test from 125 Hz to 8000 Hz is being conducted. These values apply for hearing threshold measurements to 0 dB HL. The standard specifies that MPANLs should be adjusted appropriately when hearing thresholds for pure tones are measured above and below 0 dB HL.

Maximum Ambient Noise (dB re: 20 μ Pa to the nearest 0.5 dB)

	Supra-aural Earphone	Insert Earphone	Ears not covered i.e., Sound Field
125	35	59	29
250	25	53	21
500	21	50	16
1000	26	47	13
2000	34	49	14
4000	37	50	11
8000	37	50	14

NOTE: Live voice testing requires a separate sound attenuated room for the patient to avoid feedback and direct transmission of the test stimuli.

Sound Attenuation (dB) for Earphones per ISO 4869-1

	DD45 Earphone	Insert Earphone	DD450
125	3	33.5	14.5
250	5	34.5	16.0
500	7	34.5	22.5
1000	15	35.0	28.5
2000	26	33.0	32.0
4000	32	39.5	45.5
8000	24	43.5	44.0

OPERATING THE PELLO

PRELIMINARY CHECKS

Before starting any procedures using the GSI Pello Clinical Audiometer, ensure that the power cord is plugged into a properly grounded receptacle.

Check also that all cords from the transducers, patient response hand switch (if used), and printer fit securely in their connectors on the rear and side panels. Inspect all cords for fraying and damage. If there is any damage to any cord, do not use the Pello. If speech testing with recorded voice from an external source is to be performed, check that the CD or digital music player device is connected and operating properly.

It is recommended that all parts that come into direct contact with the patient (e.g., earphone cushions) are subjected to standard disinfecting procedures between patients. This includes physically cleaning and using a recognized disinfectant. Individual manufacturer's instructions should be followed for use of any disinfecting agent to provide an appropriate level of sterilization.

Handle earphones, bone vibrator, and insert earphones with care. Do not drop them nor allow them to be banged together. Severe mechanical shock can alter their operating characteristics or change the output levels, which may require that the transducers be replaced.

PREPARE TO TEST

Once you have performed the preliminary checks and any routine maintenance checks prepare to test patients.

- √ Turn on the instrument and allow it to come to operating temperature (approximately 10 minutes).
- √ Check that the transducers and other system components are operating properly.
- √ Seat the patient comfortably in the test area.
- √ Place the selected transducers on the patient.

Do not connect or disconnect Earphones, Insert Phones, Bone Vibrator, High Frequency transducers, or any other accessories while in contact with the patient.

Placement of the Earphones

Prior to positioning the earphones on the patient's head, inspect the ear canals for any blockage due to cerumen or foreign objects. Recognize that soft-walled ear canals may collapse under the earphones, and this can lead to incorrect threshold levels. Insert phones might be used in these cases. Eliminate all obstructions, such as glasses, hair, or hearing aid, between the earphone and the patient.

Center the earphones over both ears and adjust the headband so that it rests solidly on the crown of the head and exerts pressure on both ears. Place the earphone with the red connector over the patient's right ear and the earphone with the blue connector over the left ear.

Placement of the Insert Phone

Be sure there is an eartip attached to the insert phone before inserting into the patient's ear. Inserting the insert phone without an eartip could cause harm to the patient. Push the correctly sized eartip onto the earphone and then place the insert phone securely into the patient's ear. When using the insert phones, follow the manufacturer's recommended procedure for eartip preparation, placement, and insertion.

Insert eartips are single use only. Using disposable eartips ensures sanitary conditions for each patient. Failure to comply with the single-use requirement can lead to cross-contamination and infection.

Placement of the High Frequency Transducer (optional license & DD450 headphone required)

Remove eyeglasses and earrings if possible and position the transducer directly on the head of the patient. Place the rubber cushions so that the earphone diaphragm is aimed directly at the opening into the ear canal. Adjust the headband for a tight fit. If the cushions are not tight to the ears, the test result will be false, especially at lower frequencies.

Placement of the Bone Vibrator

The bone vibrator can be placed on the promontory of the mastoid process or on the forehead, whichever has been selected in the configuration application or modified in the Configure screen.


Loud Auditory Output

The Pello audiometer is capable of outputting loud sounds and is designed to provide a visual indicator to the examiner when presenting sounds at or above 100 dB. The level display is shown with a yellow background when loud output is selected for presentation.

TEST TYPE SCREENS

The information displayed on the Pello LCD varies depending on the Test Type. There are common elements found on all screens such as the Channel 1 and 2 level settings, the Navigation menu, and the Title Bar.

Title Bar

The title bar is located at the top of the display. The title bar displays the test type in the middle. The patient's name will appear on the left side of the title bar if a patient name has been entered (or selected from an imported patient list). The right side of the title bar displays the examiner name if examiners have been entered. The examiners can be entered from the Configuration application.

Test Type Information

Under the title bar test specific information will be displayed. On the left and right side, the current output in dB HL for Channel 1 and Channel 2 will be displayed. The other

information displayed will depend on the test type and is described as part of the individual test type displays.

Navigation Menu

This menu is located at the bottom of the display. It utilizes the on-board navigation buttons or an external mouse to access the menu options. The menu is specific to the test type selected.

Time and Date

The date and time are displayed in the bottom right corner of the screen. Using the Configuration Application, the Time can be configured in a 12- or 24-hour format and the Date can be configured in any order (dd/mm/yyyy, etc.). It is also possible to update the date and time in the configuration screen of the instrument. Use the navigation buttons and the attenuator dials to set the appropriate time and date. The date and time will be updated when a new session is started.

NOTE: The time does not change automatically for daylight savings time. The operator must manually change the time using the configure button on the front panel of the instrument or the configuration application.

Common Icons

These icons are found in the test information area and common to the different test types.

Talk Forward – When pressed, a head with a headset icon will appear. This icon will remain active if the talk forward button is depressed.

Store – When either of the store buttons is pressed, a floppy disc icon flashes and the result is then displayed.

Interlock – When interlock is active, a padlock icon will appear.

Tracking – When tracking is selected, a railroad track icon will appear.


Data Transfer – When there is an active connection between the Pello and an external computer, communication will be indicated by the blue arrows.

Pencil Icon - This icon opens a comments window (must use external keyboard to utilize comment section). Comments can be entered from any test screen, and it is possible to review and edit comments from any test screen.

TONE TEST TYPE

Pressing the Tone Test Type button prepares the GSI Pello™ for pure tone air and bone conduction testing from 125 – 8,000 Hz or 125-20,000 Hz with the High Frequency option.

Each selection on the blue Navigation menu is specific to Pure Tone Testing. Pressing this button will set the defaults from the configuration application to start the test type. Auto Hz can be defined in the Config App and automatically moves the next designated test frequency when “store” is pressed. If the High Frequency license option is on the device, use the Range menu to select the desired frequency range to be tested: Standard (125-8,000 Hz), High (8,000- 20,000 Hz), or Full (125 – 20,000 Hz). If the High Frequency license option is not on the device, the range button does not appear on the Navigation menu.

NOTE: If you select a different transducer or stimulus or ear and leave the tone test type, when you return to tone, the Pello will default to the last settings that were selected.

- Press the Tone Test Type Button.
- Verify that the transducers and signals are correct.
- Verify the frequency range is correct (if licensed for High Frequency).
- Perform air conduction threshold testing.
- Press “store” after each threshold is obtained.

When the pure tone evaluation is complete, move to the next test type in your typical testing sequence.

HIGH FREQUENCY TESTING (OPTIONAL HIGH FREQUENCY LICENSE)

Use the navigation buttons or the mouse to access the Range menu. Available frequency ranges are Standard (125-8,000 Hz), High (8,000-20,000 Hz) or Full (125 – 20,000 Hz). Auto Hz can be defined in the Config App and automatically moves the next designated test frequency when “store” is pressed.

- Press the Tone Test Type Button.
- Select the desired frequency range (High or Full).
- Verify that the transducers and signals are correct.
- Perform High Frequency Testing.
- Press “store” after each threshold is obtained.

When the high frequency evaluation is complete, move to the next test type in your typical testing sequence

TEN TEST TYPE (OPTIONAL SPECIAL TEST LICENSE)

The TEN test option is available from the More test button if the Special Test license has been installed on the audiometer.

Presentation Level

- If the audiometric threshold is 60 dB or less, set the TEN level to 70 dB.
- If the audiometric threshold is 65 dB or greater, set the TEN level 10 dB higher than the threshold.
- If the TEN noise is reported to be too loud, set the TEN level at 90 dB (maximum) or the same level as audiometric threshold.

Test Instructions

When the starting level has been determined, instruct the patient in the same manner as when measuring pure tone thresholds with masking.

Test Procedure

Using the blue navigation buttons, access the TEN test type to prepare the Pello for TEN testing. The procedure for determining thresholds in the TEN is identical to the manual pure tone audiometry except that a 2 dB final step size should be used for maximum accuracy. The TEN will take approximately 4 minutes per ear (to complete all test frequencies).

Interpretation

The accepted rule is that a dead region is present when the TEN-masked threshold is at least 10 dB above the audiometric threshold and the TEN Threshold is 10 dB above the TEN noise.

SPEECH TEST TYPE BUTTON

Pressing the Speech Test type button prepares the Pello for Speech testing. The internal .Wav files can be presented by either using the present button or by a single click of a wireless mouse. The correct/incorrect/clear buttons can be used to score. It is also possible to automatically play wordlists by pressing the interrupt button. These preferences are determined in the Config App. It is critical that the test type be carefully selected as the reporting/storing is dependent upon test type. If you wish to do a PIPB rollover, you can select the speech audiogram page.

Integrated Word Files

When Speech Test Type is selected, the Pello defaults to internal speech files. These can be presented for consistent recorded speech testing.

- Define word list favorites in the Config App for fast selection of common lists.
- Utilize the navigation menu or external mouse to select the test type and the wordlist.
- Select Word Nav and use the navigation buttons to highlight word stimulus. Press the present bar to present the word.

–OR–

- Utilize the navigation buttons external mouse to present the words (single click to present).
- When the speech stimulus is being presented, the word will be highlighted yellow.

- When the patient responds (and the yellow highlight disappears), the stimulus word/sentence can be scored correct or incorrect.
- The stimulus word/sentence will turn green for correct or orange/red for incorrect. The center area of the display will indicate the % correct/#words presented.
- After the completion of each speech test type, press store to save the results in the speech results table.

When your speech evaluation is complete, move to the next test type in your typical test sequence.

MORE TEST TYPE BUTTON (OPTIONAL LICENSES)

Pressing the “More” test type button calls up a menu of the following special tests if the Special Test and/or Speech Plus license has been installed: ABLB, BKB-SIN, QuickSIN, SISI, TEN, and Tone Decay. Use the on-board navigation buttons or an external mouse to select the special test.

BKB-SIN (OPTIONAL SPEECH PLUS LICENSE)

The screenshot displays the BKB-SIN software interface. At the top, it shows 'Patient ...', 'BKB-SIN', and 'Examiner 1'. The interface is divided into two main channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Below the channels are sections for SNR Loss Averages (Group 1 and Group 2) with input fields for Right (R), Back (B), and Left (L) ears. There are also fields for PTA AC: BC and SII for both Right and Left ears. The central part of the screen features three 'Test Results Group 1' tables, each with columns for Ear, Word List, dB HL, SNR 50, and SNR Loss. Below these is a 'BKB-SIN : List Pair 1' table with 10 items (A1-A10) and their corresponding sentences. The interface also includes a 'Score' table on the right and a bottom navigation bar with buttons for Word Lists, Word Nav, Age, Aided, dB Step (set to 5), Group (set to 1), Research, and a timestamp of 10:03 AM on 02-26-2019.

Presentation Level

For standard SNR Loss testing the BKB-SIN Test should be presented at a relatively high level (loud, but below discomfort). Normative data on normal-hearing adults and normal-hearing children were collected using binaural presentation via insert earphones, at a presentation level of 70 dB HL (83 dB SPL). Normative data on adult cochlear implant users were collected using a 65 dB SPL presentation level in sound field (equivalent to 50 dB HL at 0 degrees azimuth).

Test Instructions

“You will hear a man talking to you through the earphones (or loudspeaker). He is going to say “Ready” and then he’ll say a sentence. Repeat the sentence the man says. You will hear other talkers in the background. Don’t pay any attention to them; just repeat what the man says. The background talkers will get louder, and then it will be hard for you to hear the

man's voice. When that happens, it is OK to guess; repeat anything you think you heard the man say."

Test Procedure

- Select BKB-SIN from the More Tests Menu.
- Select the proper transducer and intensity levels for each channel.
- Select the appropriate age from the Navigation Menu
- Select the appropriate word list pair. Lists 1-8 are appropriate for all patients. List pairs 9-16 are appropriate for CI patients or those with significant hearing loss.
- Using the Word Nav and front panel navigation buttons or an external mouse, select the first sentence.
- Press the present bar or click the first sentence.
- Score the four/three key words highlighted in each sentence by pressing the INCORRECT button for each word repeated incorrectly by the patient.
- Complete the entire list pair (20 sentences) and press "store".

To interpret the SNR loss score for adults, see table below.

Refer to the BKB-SIN user manual for how to interpret test results for children.

SNR Loss	Degree of SNR Loss	Expected Improvement with Directional Mic
0-3 dB	Normal/near normal	May hear better than normal hearing individuals hear in noise
3-7 dB	Mild SNR loss	May hear almost as well as normal hearing individuals hear in noise
7-15 dB	Moderate SNR loss	Directional microphones help; consider array mic
>15 dB	Severe SNR loss	Maximum SNR improvement is needed; consider FM system

QUICKSIN (OPTIONAL SPEECH PLUS LICENSE)

The screenshot displays the QuickSIN software interface. At the top, it shows 'Patient ...' and 'QuickSIN'. Below this, there are two main channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Each channel has a frequency scale from -20 to 3 dB HL. The interface also includes sections for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages', each with a grid for 'Basic', 'HFE', and 'HFE-LP' across 'R', 'B', and 'L' ears. Below these are 'Test Results Group 1' tables for each ear, with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. A 'QuickSin : List 1 (Track 3)' section contains several sentences for the test. To the right of the sentences is a 'Score' table with columns for 'S/N' and 'Score'. The bottom of the interface has a control bar with buttons for 'Word Lists', 'Word Nav', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a clock showing '10:53 AM 03-01-2019'.

Presentation Level

For pure-tone average (PTA) less than or equal to 45 dB HL, set the attenuator (s) in Channel 1 and Channel 2 to 70 dB HL. For PTA of 50 dB HL or greater, set the attenuators to a level that is judged to be “loud, but okay.” The sound should be perceived as loud, but not uncomfortably loud. Tests can be performed in the right, left or binaural conditions.

Test Instructions

“Imagine that you are at a party. There will be a woman talking and several other talkers in the background. The woman’s voice is easy to hear at first, because her voice is louder than the others. Repeat each sentence the woman says. The background talkers will gradually become louder, making it difficult to understand the woman’s voice, but please guess and repeat as much of each sentence as possible.”

Test Procedure

- Select QuickSIN from the More Tests Menu.
- Ensure the proper transducer and intensity levels are selected.
- Select the appropriate word list. Lists 1-12 are standard lists used to determine SNR loss.
- Using the Word Nav and front panel navigation buttons or an external mouse, select the first sentence.
- Press the present bar or click the selected sentence.
- Score the five key words highlighted in each sentence by pressing the INCORRECT button for each word repeated incorrectly by the patient.
- Press Store
- The SNR Loss score will appear in the SCORE/WORD window.

NOTE: It is recommended that at least 2 lists be presented in each condition. The average score will appear in the Group 1 SNR averages window.

To interpret the SNR loss score, see the table above in the BKB-SIN section.

Refer to the QuickSIN user manual for advanced testing options.

CONFIGURATION AND UPDATE

The configure button presents a Configuration screen where it is possible to view the instrument information such as serial number, software version and last calibration date. This screen displays setup options to update the Pello software, configure bone conduction symbol settings, export log files, set the date and time, and adjust the screen brightness.

The Configuration screen displays the following information:

Instrument			
App Version	1.5.3 (Build 6)	Login	No Login Required
DSPAud Version	2.16 (12/7/2015 2:50 PM)	Language	English (United States)
OS Version	1047, 5, 1.6, 15	Symbol Set	USA
Serial Number	PROTO0798	Bone	Mastoid
License Number	51815792	Display Brightness	50
Free Disk Space	2.22 MB	Mode	AWL
Free SD Card Space	10.10 GB		
Last Calibration	11/16/2016	Next Calibration	09/09/2022

Transducer		Right	Left
Headset Phone	DD450	WT000195	WT000195
Insert Phone	IP30	ID040383	ID040383
Bone Conductor	B81	003737	

At the bottom of the screen, there is a navigation bar with buttons for Update, Bone, Print, and Export Logs. On the right side of the bar, there are icons for a key and a gear, and a timestamp showing 11:46 PM 9/9/2021.

Update

Place a USB drive with the appropriate update loaded into one of the four USB ports. Select Update and then select from device, sound files, or fonts to update the instrument. Software and sound files updates must be obtained from GSI or an authorized GSI representative.

Bone

Select the symbol scheme for bone conduction testing. Choose between MAS (Mastoid) and FOR (Forehead.) This selection will be active throughout the current session. When a new session is started, the symbol scheme will revert to the configured preference.

Print

Select to change the printing format.

Export Logs

Select to obtain and send the instrument logs to an attached USB drive. The instrument logs may be requested to help troubleshoot any problems.

License

Select the License icon to see the licensed options on the instrument. In the dialog that is presented, the current options are indicated and there are buttons to allow the update or export of the license key. An update to the license key can be made manually by typing the key code into the text box or the key can also be imported from an attached USB drive.

NOTE: One license key contains all the instrument options. When a successful update occurs, the option to add Another Key refers to licensed word lists, if applicable.

Settings

Select the settings icon opens the settings dialog. Select to configure audiogram orientation, masked and unmasked symbols, printer protocol, display brightness, and to set date and time

CLEANING AND MAINTENANCE

ROUTINE MAINTENANCE

Biological Calibration Check

The design of the GSI Pello audiometer should provide trouble-free service for a long time. It is recommended to routinely record and file the audiogram of one person for the purpose of biologic calibration. This person (or group of persons) should have a known stable audiometric curve that does not exceed 25 dB HL at any frequency. This procedure should start when the GSI Pello is first installed and then be continued. Remember that individual thresholds can shift by as much as 5 dB from day to day; however, variations that exceed this range can point to difficulties which require attention.

Periodic Checks

The routine maintenance checks described below can detect possible instrument problems. If any issues are found, the instrument should receive technical service before further use. The checks should be made at periodic intervals, even if biologic checks reveal no problems.

Earphone and Bone Vibrator Cords

With extended use, all transducer cords tend to fray internally at the connectors. To evaluate the cord status, turn on the GSI Pello. Set the HL to a comfortably audible level. Place the transducer on your head. Activate both Interrupt buttons. Bend the cord next to the plug at both ends of each earphone. Listen for an intermittent signal, abrupt changes in the signal level, or a scratchy sound that coincides with the flexing of the cord. The presence of any of these conditions signifies that the cord should be replaced. Repeat this check for all transducers.

Hum and Noise

Set the GSI Pello to Tone test type with the standard earphones selected and the Channel 1 Interrupt button in the ON mode. Turn the Channel 1 Hearing Level control from 0 to 60 dB HL. Listen for low frequency hum (60 or 120 Hz) and any other noise (hiss or low rushing sound) at all attenuator levels through the earphone. Some audible noise at levels above 70 dB is permissible. If these noises are detected below 70 dB, the audiometer should be scheduled for maintenance. Repeat for Channel 2.

Distortion and Frequency Shift

Check for distortion and frequency shift by listening to the GSI Pello's output through the earphones at each frequency (in the 125 Hz to 12,000 Hz range) at a loud, but not uncomfortable level (70 to 80 dB HL for normal ears). Listen also to ensure that the signal frequencies change appropriately when the Frequency up arrow (>) and down arrow (<)

pushbuttons are operated. If distortion is heard in one earphone but not the other, the chances are high that the earphones are at fault and should be replaced. In any case, the audiometer should be scheduled for immediate maintenance.

Speech Level Check

To check the speech level with recorded speech, select the Speech test type button. Place the earphones on a person with normal hearing and present a word list at 40 dB. If intelligible speech is not heard, with the Channel 1 Hearing Level control set at 40 dB, the audiometer should be scheduled for technical service.

Internal Controls Check

Should the front panel controls lock into one state, and it is not possible to change any of the parameters, turn off the power. Wait one minute and then power on.

Bone Vibrator Check

This check must be performed in a quiet environment or in a sound room. With the frequency set to 2,000 Hz, the Channel 1 level set at 40 dB HL and the bone vibrator positioned properly, the tone should be clearly audible to a person with normal hearing – less than 25 dB. When a bone vibrator fails this test, the calibration should be verified.

Masking Level Check

Select the Tone test type. Ensure the stimulus is narrow band noise on Channel 2. Activate the Channel 2 Interrupt button and listen for a smooth, even hiss.

Talk Forward Check

Speech should be clearly audible (in the earphones) when spoken in a normal tone with the Talk Forward dB HL control set at 45 dB HL.

CALIBRATION

To guarantee that the Pello works properly, the instrument should be checked and calibrated at least once per year. Have an authorized service technician perform electrical safety checks on the unit to maintain continued compliance to IEC and UL 60601-1. The use of non-calibrated devices can lead to incorrect test results.

The service and calibration must be performed by an authorized service center. If these checks are not done, Medical Device Regulation (EU) 2017/745 and other regulations may be violated, and warranties may be void.

CLEANING THE SYSTEM

Non-disposable parts of the system, including the transducers and insert phone tubes which are in direct contact with the patient may need to be disinfected between patients.

Turn **OFF** the system and disconnect power before cleaning the instrument. Use a soft cloth lightly dampened with cleaning solution to clean all exposed surfaces. DO NOT allow liquid to contact the metal parts inside the transducers (e.g., earphones / headphones). Do not permit solutions or disinfecting agents to seep into the electronic portions of the system. Take special care around controls, connectors, and panel edges. Remove any dust from the exterior of the system with a soft brush or cloth. Use a brush to dislodge any dirt on or around the connectors and panel edges. Remove stubborn dirt with a soft cloth slightly dampened with mild detergent and water. Wipe surfaces dry afterward. Do not use instrument or transducers until they are completely dry.

Recommendations for cleaning and disinfection of the Pello presented in this document are not intended to replace or contradict policies in effect or procedures required for infection control at the facility.

Cleaning and Disinfecting Agents

According to the recommendations from the CDC, audiometric equipment is non-critical medical equipment and typically requires cleansing followed by low to intermediate level disinfecting, depending on the nature of the contamination. Cleaning should be done with a mild soapy detergent (such as dishwashing liquid) and a damp cloth, or an Endozime Sponge followed by an application of EPA-registered hospital disinfectant. Do not use any abrasive cleaners.

Use of a non-alcohol-based disinfectant is recommended for larger areas and headphones. Non-alcohol-based products that contain the active ingredient referred to as quaternary ammonia compound or hydrogen peroxide-based cleaner such as Oxivir Disinfectant Wipes are recommended to clean the ear cushions, headset, and to wipe down the machine. The quaternary ammonia compound and hydrogen peroxide are specifically designed to disinfect rubber, plastic, silicone, and acrylic products which are commonly used in hearing evaluation instruments.

Single-use Accessories

Disposable or single use supplies such as eartips for the insert earphones are intended to be discarded immediately after use. Use of disposable items on multiple patients can lead to cross-contamination and infection.

STATUS AND ERROR MESSAGES

ERROR MESSAGES

Please try another selection

Indicates an incorrect selection. This could include actions such as incompatible transducers, incompatible routing, or no calibration data stored for the selected transducers.

No test data stored

Indicates that there is no test data available to be erased, printed or transferred.

Printer communication error

Occurs during printing if communication problems with the printer are encountered.

Record test result in comments

Test results of the ABLB and Tone Decay are not recorded directly on a report. This message indicates that the results should be documented in the comments.

The session comments have been updated with the results of the SDT test

This message indicates that the stored speech detection threshold results will appear in the comments section and will be printed directly or transferred electronically.

Not supported in speech

The selected action is not supported in the speech test type.

Speech data limit exceeded, speech tables limited to 6 test results per ear. Latest test result will not be saved

Up to six speech tests can be stored for each ear. This message indicates that the maximum number of tests have been stored and the latest test cannot be added.

Error

If there are general system errors, a dialog box with "Error" in the title will be shown with the given error.

Error

The startup configuration for this test type is not fully calibrated; a search for a different configuration that is calibrated has found the currently displayed configuration

This message indicates that the selected transducers have not been calibrated.

SUPPLIES AND ACCESSORIES

To order consumables, additional accessories and to replace detachable parts that have been damaged, please contact GSI or your GSI distributor for current prices and delivery charges. Before use check the packaging for any shelf-life expirations.

Part Number	Description
8500090	Foam Eartips for Insert Phones, 13 mm (Adult), 50 Each
8106226	Foam Eartips for Insert Phones, 13 mm (Adult), 100 Each
8500110	Foam Eartips for Insert Phones, 10 mm (Pediatric), 50 Each
8106227	Foam Eartips for Insert Phones, 10 mm (Pediatric), 100 Each
8500100	Foam Eartips for Insert Phones, Jumbo, 24 Each
8500140	Connectors (nubbins) for Insert Phones to Silicone Tubes, 10 Each
8501072	Sound tubes for Insert Phones, 1 red, 1 blue
8109155	AMBAND - Calibrated Bone Conduction Headbands (10)
8522201	Earpads for monitor headset with boom mic (2)
8504476	Microphone windscreen (Sennheiser PS 01 for PC 131 headset)
8107419	Ear cushion covers (500) 3 inch (for DD45)
8513404	Ear cushion covers (100) 4 inch (for DD450 and DD65 V2)

APPENDIX A – TECHNICAL DATA

The Pello is an active, diagnostic medical product. The device is classified as a class IIa device according to the Medical Device Regulation (EU) 2017/745 and a class II device according to the US FDA. Manufactured, designed, developed, and marketed under ISO 13485 certified quality systems.

Standards Compliance	
Safety and Electromagnetic compatibility (EMC)	IEC 60601-1, IEC 60601-1-2
Calibration and Test Signal	ANSI S3.6 (2018) Type 2 IEC 60645-1 (2017) Type 2 IEC 60645-2 (2017) ISO 389-2 ISO 389-3 ISO 389-4 ISO 389-5 ISO 389-7 ISO 389-8
General Specifications	
ENVIRONMENTAL	
Transport and Handling 	Transport package shall be kept away from rain and in dry conditions
Temperature 	Operation: + 15° C to + 40° C (+ 59° F to + 104° F) Transport: - 20° C to + 50° C (- 4° F to + 122° F) Storage: 0° C to + 50° C (+32° F to + 122° F)
Humidity 	Operation: Maximum relative humidity 95 %, non-condensing Transport and Storage: Maximum relative humidity 93 %, non-condensing
Ambient air pressure	98 kPa – 104 kPa
Location	Indoor Use
Expected Service Life	7 years

POWER	
Power Supply	100 V AC – 240 V AC, 0.5 A max 50/60 Hz, 90 W
PHYSICAL	
Dimensions	37.5 x 26.7 x 35.1 cm (L x W x H) 14.8 x 10.5 x 13.8 inches (LCD raised)
Weight	3.7 kg (8.18 lbs)
System	
Channels	1.5
Pure Tone Stimuli	
Frequency Range	125-20,000 Hz* (transducer dependent)
Frequency Accuracy	± 1 %
Total Harmonic Distortion	< 2% (earphones and insert phones) < 5.5% (B81 bone vibrator)
Level Range	-20 – 120 dB HL (transducer dependent)
Signal Format	Steady: Tone continuously present Pulsed: Tone pulsed 200 msec ON, 200 msec OFF FM: Modulation Rate: 5 Hz, Modulation depth +/- 5% Pediatric Noise*: Continuously presented or pulsed
*Note: Optional license required. High Frequency license required for extended frequency range. Special Test license required for Pediatric Noise.	

TRANSDUCERS	
RadioEar IP30 Insert Earphones	Frequency Range: 125 Hz – 8000 Hz Level: -10 dB HL to 120 dB HL
RadioEar DD45 Headphones	Frequency Range: 125 Hz – 8000 Hz Level: -10 dB HL to 120 dB HL
RadioEar DD450 Headphones	Frequency Range: 125 Hz – 20,000 Hz Level: -20 dB HL to 100 dB HL
RadioEar B81 Bone Vibrator	Frequency Range: 250 Hz – 8000 Hz Level: -10 dB HL to 90 dB HL (mastoid) Level: -10 dB HL to 80 dB HL (forehead)
Sound Field	Frequency Range: 125 Hz – 8000 Hz Level: -10 dB HL to 90 dB HL (amplified speakers) Level: -10 dB HL to 102 dB HL (high performance speakers and external booster amplifier)
Note: The maximum HL values are applicable to the middle frequencies only. Optional licenses or configuration required for High Frequency, some transducers, and Sound Field speakers.	
Speech	
Microphone	For live voice testing and communications
INT/EXT A & INT/EXT B	Can be utilized for internal speech files or recorded speech material from an external digital device
Level Range	Air Conduction: -10 dB HL to 100 dB HL (DD45 Linear), Inserts (IP30): -10 dB HL to 95 dB HL Bone Conduction (B81): -10 dB HL to 60 dB HL (mastoid), -10 dB HL to 50 dB HL (forehead) Sound Field: -10 dB HL to 85 dB HL (amplified speakers)
Masking	

Narrow Band Noise (effective masking)	<p>Maximum dB HL is 15 dB below tone</p> <p>Frequency Range: 3 dB down bandwidth, 1/3 octave minimum; ½ octave maximum of test signal</p>
Speech Noise	<p>Air Conduction: -10 dB HL to 95 dB HL (DD45) Bone Conduction: -10 dB HL to 50 dB HL (mastoid B81), -10 dB HL to 40 dB HL (forehead B81) Sound Field: -10 dB HL to 85 dB HL</p> <p>Frequency Range: Equal energy per frequency 100 to 1,000 Hz with a 12 dB/octave roll-off from 1,000 to 6,000 Hz</p>
White Noise	<p>Air Conduction: -10 dB HL to 95 dB HL (DD45) Bone Conduction: -10 dB HL to 60 dB HL (mastoid B81), -10 dB HL to 50 dB HL (forehead B81) Sound Field -10 dB HL to 80 dB HL</p> <p>Frequency Range: 125 Hz to 12,000 Hz with constant bandwidth</p>
Input Specification	
External CD	7 mV RMS at max. gain for 0 dB reading; input impedance 47 kOhm
Talk Forward	100 uV RMS at max. gain for 0dB reading; input impedance 3.2 kOhm
Output Specification	
Headphones	7 V RMS at 10 Ohm load; 60-20,000 Hz -3 dB
Insert Phones	7 V RMS at 10 Ohm load; 60-20,000 Hz -3 dB
Bone	7 V RMS at 10 Ohm load; 60-10,000 Hz -3 dB

Sound field	7 V RMS at min. 2 kOhm load; 60-20,000 kHz -3 dB
Monitor	2x 3 V RMS at 32 Ohm/1.5 V RMS at 8 Ohm load; 60-20,000Hz -3dB
Free Field Specification	
Power Amplifier and Speakers	With an input of 7 V rms-Amplifier and loudspeakers must be able to create a sound pressure level of 100 dB at 1 meter and meet the following requirements
Frequency Response	125-250Hz, (+0/-10 dB) 250-4000Hz (± 3 dB) 4000-6300Hz (± 5 dB)
Total Harmonic Distortion	80 dB SPL <3% 100 dB SPL <10%
Communication and Monitoring	
Talk Forward	Permits the tester to speak through the test microphone into the selected transducer at approximately the level set by the front panel controls
Talk Back	Allows the tester to listen to comments from the patient in the testing booth
Monitor	The monitor headset can be used by the tester to listen to Channel 1, Channel 2, and/or Talk Back signals
Special Tests	
QuickSIN*	Six (6) sentences with five (5) key words per sentence are presented in four-talker babble noise. The sentences are presented at pre-recorded signal-to-noise ratios. The SNR's used are 25, 20, 15, 10, 5, and 0.


BKB-SIN*	18 List Pairs. The sentences are presented at prerecorded signal-to-noise ratios that decrease in 3-dB steps. Each list in the pair is individually scored, and the results of the two lists are averaged to obtain the List Pair score. Results are compared to normative data to obtain the SNR Loss
TEN*	TEN masking noise will be presented to the test ear. Pure tone stimuli between 500 and 4000 Hz can be used at 1, 2, or 5 dB increments to obtain TEN thresholds
ALT (ABLB)*	Tone alternating between Channel 1 and Channel 2: Channel 1 is 400 msec ON, 400 msec OFF followed by Channel 2, 400 msec ON, 400 msec OFF
SISI*	A level increment is added to a tone in the selected channel for 200 msec, every 5 seconds. The HL increments are in 1, 2 or 5 dB
Others (User defined)	MLB Lombard test Pure Tone Stenger Speech Stenger SAL
*Note: Optional license required. Special Test license required for ABLB, SISI, TEN, and Tone Decay. Speech Plus license required for BKB-SIN and QuickSIN.	

APPENDIX B – SPECIFICATION OF INPUT/OUTPUT CONNECTIONS


MAINS CABLE

Pin	Description	
1	Live	 <p>IEC C6</p>
2	Neutral	
3	Earth	


HEADPHONES/INSERTS (2 CONNECTORS, R, L)

Pin	Description	
1	Ground	 <p>6.3mm Mono</p>
2	Signal	


BONE OSCILLATOR

Pin	Description	
1	Ground	 <p>6.3mm Mono</p>
2	Signal	


Patient Response

Pin	Description	
1		 <p>6.3mm Stereo</p>
2	Contact Closure 	
3		


TALK FORWARD, TALK BACK MIC

Pin	Description	
1	Ground	 <p>3.5mm Stereo</p>
2	DC Bias	
3	Signal	


MONITOR HEADSET

Pin	Description	
1	Ground	 <p>3.5mm Stereo</p>
2	Signal 2	
3	Signal 1	

EXTERNAL A+B

Pin	Description	
1	Ground	 <p>3.5mm Stereo</p>
2	External B (CD 2)	
3	External A (CD 1)	


SPEAKERS (FF1, FF2)

Pin	Description	
1	Ground	 <p>RCA</p>
2	Signal	

USB A (HOST) x4

Pin	Description	
1	+ 5 V DC	
2	Data	
3	Data	
4	GND	

USB B (DEVICE)

Pin	Description	
1	+ 5 V DC	
2	Data	
3	Data	
4	GND	

LAN ETHERNET

Not Supported

APPENDIX C - CALIBRATION REFERENCE & MAXIMUM LEVELS

The Pello is supplied from the factory calibrated for the transducers that were purchased with it. The exception is the speakers, as those must be calibrated in the environment where they will be used. The calibration data supplied from the factory is only valid for GSI supplied transducers and cannot be applied to non-GSI supplied transducers.

It is recommended that calibration of the instrument and transducers be performed annually by authorized GSI Representatives using appropriate calibration instrumentation. If periodic checks are also desired, the tables in this section provide the SPL values per frequency for each transducer. If the measured values are not within ± 5 dB at 125, 6,000, and 8,000 Hz in the earphones, the GSI Pello should be scheduled for immediate maintenance.

It is not possible to select a dB HL value outside the limits for a particular transducer/frequency combination. An attempt to change or select a hearing level control that is outside of the limit will cause the dB HL display to flash momentarily and then the test channel value will be replaced with NR (No Response). If an audiogram is displayed and the limits for a frequency/transducer are reached, the symbol for no response is displayed in the audiogram.

It is not possible to select a test frequency that is invalid for a particular transducer.

The hearing levels listed in the Max HL tables are maximum levels. These levels are achievable only if ANSI, ISO or GSI reference threshold levels, and not customized calibration values, are used. At no time will the hearing level limit exceed 120 dB HL.

EARPHONES - PURE TONE RETSPL

Transducer	DD45	DD45	TDH50	TDH50	DD450
Impedance	10 Ω	10 Ω	60 Ω	60 Ω	23 Ω
Coupler	318-3	318-1	318-3	318-1	318-1
	RETSPL	RETSPL	RETSPL	RETSPL	RETSPL
125 Hz	47.5	45	47.5	45	30.5
160 Hz	40.5	38.5	40	38.5	26
200 Hz	33.5	32.5	33.5	32.5	22
250 Hz	27	27	26.5	27	18
315 Hz	22.5	22	22	22	15.5
400 Hz	17.5	17	17.5	17	13.5
500 Hz	13	13.5	13.5	13.5	11
630 Hz	9	10.5	10.5	10.5	8
750 Hz	6.5	9	8.5	9	6
800 Hz	6.5	8.5	8.5	8.5	6
1000 Hz	6	7.5	7.5	7.5	5.5
1250 Hz	7	7.5	7.5	7.5	6
1500 Hz	8	7.5	7.5	7.5	5.5
1600 Hz	8	8	8.5	8	5.5
2000 Hz	8	9	11	9	4.5
2500 Hz	8	10.5	10	10.5	3
3000 Hz	8	11.5	9.5	11.5	2.5
3150 Hz	8	11.5	9.5	11.5	4
4000 Hz	9	12	10.5	12	9.5
5000 Hz	13	11	12	11	14
6000 Hz	20.5	16	13.5	16	17
6300 Hz	19	21	13.5	21	17.5
8000 Hz	12	15.5	13	15.5	17.5
9000 Hz					19
10000 Hz					22
11200 Hz					23
12000 Hz			17.5	11	0
12500 Hz					27.5
14000 Hz					35
16000 Hz					56
18000 Hz					83
20000 Hz					105

- DD45 6ccm uses IEC60318-3 or NBS 9A coupler and RETSPL comes from PTB – DTU report 2009-2010. Force 4.5N \pm 0.5N
- DD45 Artificial ear uses IEC60318-1 coupler and RETSPL comes from ANSI S3.6 2018 and ISO 389-1 1998. Force 4.5N \pm 0.5N
- TDH50 6ccm uses IEC60318-3 or NBS 9A coupler and RETSPL comes from ANSI S3.6 2018. Force 4.5N \pm 0.5N
- TDH50 Artificial ear uses IEC60318-1 coupler and RETSPL comes from ANSI S3.6 2018 and ISO 389-1 1998 Force 4.5N \pm 0.5N
- HDA200 Artificial ear uses IEC60318-1 coupler with type 1 adaptor and RETSPL comes from ANSI S3.6 2018 and ISO 389-8 2004. Force 9N \pm 0.5N

EARPHONES - ANSI SPEECH RETSPL

Transducer	DD45	DD45	TDH50	TDH50	DD450
Impedance	10 Ω	10 Ω	60 Ω	60 Ω	23 Ω
Coupler	318-3	318-1	318-3	318-1	318-1
	RETSPL	RETSPL	RETSPL	RETSPL	RETSPL
Speech	18.5	20	20	20	19
Speech Equ.FF.	18.5	19.5	17	18	18.5
Speech Non-linear	6	7.5	7.5	7.5	5.5
Speech noise	18.5	20	20	20	19
Speech noise Equ.FF.	18.5	19.5	17	18	18.5
Speech noise Non-linear	6	7.5	7.5	7.5	5.5
White noise in speech	21	22.5	22.5	22.5	21.5

- DD45 (G_F-G_C) PTB-DTU report 2009-2010.
- TDH50 (G_F-G_C) ANSI S3.6 2018.
- HDA200 (G_F-G_C) ANSI S3.6 2018 and ISO 389-8 2004.
- ANSI Speech level 12.5 dB + 1 kHz RETSPL ANSI S3.6 2018 (acoustical linear weighting)
- ANSI Speech Equivalent free field level 12.5 dB + 1 kHz RETSPL - (G_F-G_C) from ANSI S3.6 2018(acoustical equivalent sensitivity weighting)
- ANSI Speech Not linear level 1 kHz RETSPL ANSI S3.6 2018 (DD45-TDH50-HDA200) and IP30- B81 12.5 dB + 1 kHz RETSPL ANSI S3.6 2018 (no weighting)

EARPHONES - IEC SPEECH RETSPL

Transducer	DD45	DD45	TDH50	TDH50	DD450
Impedance	10 Ω	10 Ω	60 Ω	60 Ω	23 Ω
Coupler	318-3	318-1	318-3	318-1	318-1
	RETSPL	RETSPL	RETSPL	RETSPL	RETSPL
Speech	20	20	20	20	20
Speech Equ.FF.	3.5	4.5	2	3	3.5
Speech Non-linear	6	7.5	7.5	7.5	5.5
Speech noise	20	20	20	20	20
Speech noise Equ.FF.	3.5	4.5	2	3	3.5
Speech noise Non-linear	6	7.5	7.5	7.5	5.5
White noise in speech	22.5	22.5	22.5	22.5	22.5

- DD45 (G_F-G_C) PTB-DTU report 2009-2010.
- TDH50 (G_F-G_C) ANSI S3.6 2018.
- HDA200 (G_F-G_C) ANSI S3.6 2018 and ISO 389-8 2004.
- IEC Speech level IEC60645-2 2017 (acoustical linear weighting)
- IEC Speech Equivalent free field level (G_F-G_C) from IEC60645-2 2017 (acoustical equivalent sensitivity weighting)
- IEC Speech Not linear level 1 kHz RETSPL (DD45-TDH50-HDA200) and IP30 - B81 IEC60645-2 2017 (no weighting)

EARPHONES - PURE TONE MAX HL

Transducer	DD45	DD45	TDH50	TDH50	DD450
Impedance	10 Ω	10 Ω	60 Ω	60 Ω	23 Ω
Coupler	318-3	318-1	318-3	318-1	318-1
	Max HL	Max HL	Max HL	Max HL	Max HL
Tone 125 Hz	90	90	85	85	100
Tone 160 Hz	95	95	90	90	105
Tone 200 Hz	100	100	95	95	105
Tone 250 Hz	110	110	105	105	110
Tone 315 Hz	115	115	110	110	115
Tone 400 Hz	120	120	115	115	115
Tone 500 Hz	120	120	120	120	115
Tone 630 Hz	120	120	120	120	120
Tone 750 Hz	120	120	120	120	120
Tone 800 Hz	120	120	120	120	120
Tone 1000 Hz	120	120	120	120	120
Tone 1250 Hz	120	120	120	120	110
Tone 1500 Hz	120	120	120	120	115
Tone 1600 Hz	120	120	120	120	115
Tone 2000 Hz	120	120	120	120	115
Tone 2500 Hz	120	120	120	120	115
Tone 3000 Hz	120	120	120	120	115
Tone 3150 Hz	120	120	120	120	115
Tone 4000 Hz	120	120	120	120	115
Tone 5000 Hz	120	120	115	115	105
Tone 6000 Hz	115	115	115	110	105
Tone 6300 Hz	115	110	110	105	105
Tone 8000 Hz	110	110	100	100	105
Tone 9000 Hz					100
Tone 10000 Hz					100
Tone 11200 Hz					95
Tone 12000 Hz			90	90	
Tone 12500 Hz					90
Tone 14000 Hz					80
Tone 16000 Hz					60
Tone 18000 Hz					30
Tone 20000 Hz					15

EARPHONES - NB NOISE EFFECTIVE MASKING LEVEL

Transducer	DD45	DD45	TDH50	TDH50	DD450
Impedance	10 Ω	10 Ω	60 Ω	60 Ω	23 Ω
Coupler	318-3	318-1	318-3	318-1	318-1
	EM	EM	EM	EM	EM
NB 125 Hz	51.5	49	51.5	49	34.5
NB 160 Hz	44.5	42.5	44	42.5	30
NB 200 Hz	37.5	36.5	37.5	36.5	26
NB 250 Hz	31	31	30.5	31	22
NB 315 Hz	26.5	26	26	26	19.5
NB 400 Hz	21.5	21	21.5	21	17.5
NB 500 Hz	17	17.5	17.5	17.5	15
NB 630 Hz	14	15.5	15.5	15.5	13
NB 750 Hz	11.5	14	13.5	14	11
NB 800 Hz	11.5	13.5	13.5	13.5	11
NB 1000 Hz	12	13.5	13.5	13.5	11.5
NB 1250 Hz	13	13.5	13.5	13.5	12
NB 1500 Hz	14	13.5	13.5	13.5	11.5
NB 1600 Hz	14	14	14.5	14	11.5
NB 2000 Hz	14	15	17	15	10.5
NB 2500 Hz	14	16.5	16	16.5	9
NB 3000 Hz	14	17.5	15.5	17.5	8.5
NB 3150 Hz	14	17.5	15.5	17.5	10
NB 4000 Hz	14	17	15.5	17	14.5
NB 5000 Hz	18	16	17	16	19
NB 6000 Hz	25.5	21	18.5	21	22
NB 6300 Hz	24	26	18.5	26	22.5
NB 8000 Hz	17	20.5	18	20.5	22.5
NB 9000 Hz					24
NB 10000 Hz					27
NB 11200 Hz					28
NB 12000 Hz			22.5	16	
NB 12500 Hz					32.5
NB 14000 Hz					40
NB 16000 Hz					61
NB 18000 Hz					88
NB 20000 Hz					110
White noise	0	0	0	0	0
TEN noise	25	25	24.5	24.5	31.5

EARPHONES - NB NOISE MAX HL

Transducer	DD45	DD45	TDH50	TDH50	DD450
Impedance	10 Ω	10 Ω	60 Ω	60 Ω	23 Ω
Coupler	318-3	318-1	318-3	318-1	318-1
	Max HL	Max HL	Max HL	Max HL	Max HL
NB 125 Hz	75	75	65	65	75
NB 160 Hz	80	80	70	70	80
NB 200 Hz	90	90	80	80	80
NB 250 Hz	95	95	85	85	85
NB 315 Hz	100	100	90	90	90
NB 400 Hz	105	105	95	95	95
NB 500 Hz	110	110	100	100	95
NB 630 Hz	110	110	100	100	95
NB 750 Hz	110	110	105	105	100
NB 800 Hz	110	110	105	105	100
NB 1000 Hz	110	110	105	105	100
NB 1250 Hz	110	110	105	105	95
NB 1500 Hz	110	110	105	105	100
NB 1600 Hz	110	110	105	105	100
NB 2000 Hz	110	110	100	100	100
NB 2500 Hz	110	110	100	100	100
NB 3000 Hz	110	110	100	100	100
NB 3150 Hz	110	110	100	100	100
NB 4000 Hz	110	110	100	100	100
NB 5000 Hz	110	110	100	100	95
NB 6000 Hz	105	105	95	95	90
NB 6300 Hz	105	100	95	90	90
NB 8000 Hz	100	100	90	85	90
NB 9000 Hz					85
NB 10000 Hz					85
NB 11200 Hz					80
NB 12000 Hz			75	75	
NB 12500 Hz					75
NB 14000 Hz					70
NB 16000 Hz					50
NB 18000 Hz					20
NB 20000 Hz					0
White noise	120	120	120	120	115
TEN noise	110	110	100	100	121

EARPHONES - ANSI SPEECH MAX HL

Transducer	DD45	DD45	TDH50	TDH50	DD450
Impedance	10 Ω	10 Ω	60 Ω	60 Ω	23 Ω
Coupler	318-3	318-1	318-3	318-1	318-1
	Max HL	Max HL	Max HL	Max HL	Max HL
Speech	110	105	100	100	90
Speech Equ.FF.	100	100	90	90	85
Speech Non-linear	120	120	115	115	110
Speech noise	100	100	95	90	85
Speech noise Equ.FF.	100	95	85	85	80
Speech noise Non-linear	115	115	110	110	105
White noise in speech	95	95	95	95	90

EARPHONES - IEC SPEECH MAX HL

Transducer	DD45	DD45	TDH50	TDH50	DD450
Impedance	10 Ω	10 Ω	60 Ω	60 Ω	23 Ω
Coupler	318-3	318-1	318-3	318-1	318-1
	Max HL	Max HL	Max HL	Max HL	Max HL
Speech	110	105	100	100	90
Speech Equ.FF.	115	115	105	105	100
Speech Non-linear	120	120	115	115	110
Speech noise	100	100	95	90	85
Speech noise Equ.FF.	115	110	100	100	95
Speech noise Non-linear	115	115	110	110	105
White noise in speech	95	95	95	95	90

INSERT EARPHONES - PURE TONE RETSPL

Transducer	IP30	IP30
Impedance	10 Ω	10 Ω
Coupler	2ccm	711
	RETSPL	RETSPL
125 Hz	26	28
160 Hz	22	24.5
200 Hz	18	21.5
250 Hz	14	17.5
315 Hz	12	15.5
400 Hz	9	13
500 Hz	5.5	9.5
630 Hz	4	7.5
750 Hz	2	6
800 Hz	1.5	5.5
1000 Hz	0	5.5
1250 Hz	2	8.5
1500 Hz	2	9.5
1600 Hz	2	9.5
2000 Hz	3	11.5
2500 Hz	5	13.5
3000 Hz	3.5	13
3150 Hz	4	13
4000 Hz	5.5	15
5000 Hz	5	18.5
6000 Hz	2	16
6300 Hz	2	16
8000 Hz	0	15.5
9000 Hz		
10000 Hz		
11200 Hz		
12000 Hz		
12500 Hz		
14000 Hz		
16000 Hz		
18000 Hz		
20000 Hz		

- IP30 2ccm uses ANSI S3.7-1995 IEC60318-5 coupler (HA-2 with 5mm rigid Tube) and RETSPL comes from ANSI S3.6 2018 and ISO 389-2 1994
- IP30 Ear simulator uses ANSI S3.25- IEC60318-4 coupler and RETSPL comes from ANSI S3.6 2018 and ISO 389-2 1994

INSERT EARPHONES - ANSI SPEECH RETSPL

Transducer	IP30	IP30
Impedance	10 Ω / 50 Ω	10 Ω / 50 Ω
Coupler	2ccm	711
	RETSPL	RETSPL
Speech		
Speech Equ.FF.		
Speech Non-linear	12.5	18
Speech noise		
Speech noise Equ.FF.		
Speech noise Non-linear	12.5	18
White noise in speech	15	20.5

- ANSI Speech level 12.5 dB + 1 kHz RETSPL ANSI S3.6 2018 (acoustical linear weighting)
- ANSI Speech Equivalent free field level 12.5 dB + 1 kHz RETSPL – (GF-GC) from ANSI S3.6 2018(acoustical equivalent sensitivity weighting)
- ANSI Speech Not linear level 1 kHz RETSPL ANSI S3.6 2018 (DD45-TDH50-HDA200) and IP30- B81 12.5 dB + 1 kHz RETSPL ANSI S3.6 2018 (no weighting)

INSERT EARPHONES - IEC SPEECH RETSPL

Transducer	IP30	IP30
Impedance	10 Ω / 50 Ω	10 Ω / 50 Ω
Coupler	2ccm	711
	RETSPL	RETSPL
Speech		
Speech Equ.FF.		
Speech Non-linear	20	20
Speech noise		
Speech noise Equ.FF.		
Speech noise Non-linear	20	20
White noise in speech	22.5	22.5

- IEC Speech level IEC60645-2 2017 (acoustical linear weighting)
- IEC Speech Equivalent free field level (G_F-G_C) from IEC60645-2 2017 (acoustical equivalent sensitivity weighting)
- IEC Speech Not linear level 1 kHz RETSPL (DD45-TDH50-HDA200) and IP30 -B81 IEC60645-2 2017 (no weighting)

INSERT EARPHONES - PURE TONE MAX HL

Transducer	IP30	IP30
Impedance	10 Ω	10 Ω
Coupler	2ccm	711
	Max HL	Max HL
Tone 125 Hz	90	90
Tone 160 Hz	95	95
Tone 200 Hz	100	100
Tone 250 Hz	105	105
Tone 315 Hz	105	105
Tone 400 Hz	110	110
Tone 500 Hz	110	110
Tone 630 Hz	115	115
Tone 750 Hz	115	115
Tone 800 Hz	115	115
Tone 1000 Hz	120	120
Tone 1250 Hz	120	120
Tone 1500 Hz	120	120
Tone 1600 Hz	120	120
Tone 2000 Hz	120	120
Tone 2500 Hz	120	120
Tone 3000 Hz	120	120
Tone 3150 Hz	120	120
Tone 4000 Hz	115	115
Tone 5000 Hz	105	105
Tone 6000 Hz	100	100
Tone 6300 Hz	100	100
Tone 8000 Hz	95	95
Tone 9000 Hz		
Tone 10000 Hz		
Tone 11200 Hz		
Tone 12000 Hz		
Tone 12500 Hz		
Tone 14000 Hz		
Tone 16000 Hz		
Tone 18000 Hz		
Tone 20000 Hz		

INSERT EARPHONES - NB NOISE EFFECTIVE MASKING LEVEL

Transducer	IP30	IP30
Impedance	10 Ω	10 Ω
Coupler	2ccm	711
	EM	EM
NB 125 Hz	30	32
NB 160 Hz	26	28.5
NB 200 Hz	22	25.5
NB 250 Hz	18	21.5
NB 315 Hz	16	19.5
NB 400 Hz	13	17
NB 500 Hz	9.5	13.5
NB 630 Hz	9	12.5
NB 750 Hz	7	11
NB 800 Hz	6.5	10.5
NB 1000 Hz	6	11.5
NB 1250 Hz	8	14.5
NB 1500 Hz	8	15.5
NB 1600 Hz	8	15.5
NB 2000 Hz	9	17.5
NB 2500 Hz	11	19.5
NB 3000 Hz	9.5	19
NB 3150 Hz	10	19
NB 4000 Hz	10.5	20
NB 5000 Hz	10	23.5
NB 6000 Hz	7	21
NB 6300 Hz	7	21
NB 8000 Hz	5	20.5
NB 9000 Hz		
NB 10000 Hz		
NB 11200 Hz		
NB 12000 Hz		
NB 12500 Hz		
NB 14000 Hz		
NB 16000 Hz		
NB 18000 Hz		
NB 20000 Hz		
White noise	0	0
TEN noise	16	25

INSERT EARPHONES - NB NOISE MAX HL

Transducer	IP30	IP30
Impedance	10 Ω	10 Ω
Coupler	2ccm	711
	Max HL	Max HL
NB 125 Hz	90	90
NB 160 Hz	95	95
NB 200 Hz	100	100
NB 250 Hz	105	105
NB 315 Hz	105	105
NB 400 Hz	105	105
NB 500 Hz	110	110
NB 630 Hz	110	110
NB 750 Hz	110	110
NB 800 Hz	110	110
NB 1000 Hz	110	110
NB 1250 Hz	110	110
NB 1500 Hz	110	110
NB 1600 Hz	110	110
NB 2000 Hz	110	110
NB 2500 Hz	110	110
NB 3000 Hz	110	110
NB 3150 Hz	110	110
NB 4000 Hz	110	110
NB 5000 Hz	105	105
NB 6000 Hz	100	100
NB 6300 Hz	100	100
NB 8000 Hz	95	90
NB 9000 Hz		
NB 10000 Hz		
NB 11200 Hz		
NB 12000 Hz		
NB 12500 Hz		
NB 14000 Hz		
NB 16000 Hz		
NB 18000 Hz		
NB 20000 Hz		
White noise	110	110
TEN noise	100	100

INSERT EARPHONES - ANSI SPEECH MAX HL

Transducer	IP30	IP30
Impedance	10 Ω	10 Ω
Coupler	2ccm	711
	Max HL	Max HL
Speech (Non-linear)	110	105
Speech noise Non-linear	100	100
White noise in speech	95	85

INSERT EARPHONES - IEC SPEECH MAX HL

Transducer	IP30	IP30
Impedance	10 Ω	10 Ω
Coupler	2ccm	711
	Max HL	Max HL
Speech (Non-linear)	100	105
Speech noise Non-linear	90	100
White noise in speech	85	85

BONE VIBRATOR - PURE TONE RETFL

Transducer	B81	B81
Impedance	10 Ω	10 Ω
Coupler	Mastoid	Forehead
	RETFL	RETFL
125 Hz		
160 Hz		
200 Hz		
250 Hz	67	79
315 Hz	64	76.5
400 Hz	61	74.5
500 Hz	58	72
630 Hz	52.5	66
750 Hz	48.5	61.5
800 Hz	47	59
1000 Hz	42.5	51
1250 Hz	39	49
1500 Hz	36.5	47.5
1600 Hz	35.5	46.5
2000 Hz	31	42.5
2500 Hz	29.5	41.5
3000 Hz	30	42
3150 Hz	31	42.5
4000 Hz	35.5	43.5
5000 Hz	40	51
6000 Hz	40	51
6300 Hz	40	50
8000 Hz	40	50
9000 Hz		
10000 Hz		
11200 Hz		
12000 Hz		
12500 Hz		
14000 Hz		
16000 Hz		
18000 Hz		
20000 Hz		

- B81 uses ANSI S3.13 or IEC60318-6 2007 mechanical coupler and RETFL come from ANSI S3.6 2018 and ISO 389-3 1994. Force 5.4N ±0.5N

BONE VIBRATOR - ANSI SPEECH RETSPL

Transducer	B81	B81
Impedance	10 Ω	10 Ω
Coupler	Mastoid	Forehead
	RETFL	RETFL
Speech (Non-linear)	55	63.5
Speech noise Non-linear	55	63.5
White noise in speech	57.5	66

- ANSI Speech level 12.5 dB + 1 kHz RETSPL ANSI S3.6 2018 (acoustical linear weighting)
- ANSI Speech Equivalent free field level 12.5 dB + 1 kHz RETSPL - (GF-GC) from ANSI S3.6 2018(acoustical equivalent sensitivity weighting)
- ANSI Speech Not linear level 1 kHz RETSPL ANSI S3.6 2018 (DD45-TDH50-HDA200) and IP30- B81 12.5 dB + 1 kHz RETSPL ANSI S3.6 2018 (no weighting)

BONE VIBRATOR - IEC SPEECH RETSPL

Transducer	B81	B81
Impedance	10 Ω	10 Ω
Coupler	Mastoid	Forehead
	RETFL	RETFL
Speech (Non-linear)	55	63.5
Speech noise Non-linear	55	63.5
White noise in speech	57.5	66

- IEC Speech level IEC60645-2 2017 (acoustical linear weighting)
- IEC Speech Equivalent free field level (G_F-G_C) from IEC60645-2 2017 (acoustical equivalent sensitivity weighting)
- IEC Speech Not linear level 1 kHz RETSPL (DD45-TDH50-HDA200) and IP30 - B81 IEC60645-2 2017 (no weighting)

BONE VIBRATOR - PURE TONE MAX HL

Transducer	B81	B81
Impedance	10 Ω	10 Ω
Coupler	Mastoid	Forehead
	Max HL	Max HL
Tone 125 Hz		
Tone 160 Hz		
Tone 200 Hz		
Tone 250 Hz	50	35
Tone 315 Hz	60	45
Tone 400 Hz	70	55
Tone 500 Hz	70	55
Tone 630 Hz	75	60
Tone 750 Hz	75	60
Tone 800 Hz	75	60
Tone 1000 Hz	85	75
Tone 1250 Hz	90	80
Tone 1500 Hz	90	80
Tone 1600 Hz	90	75
Tone 2000 Hz	90	75
Tone 2500 Hz	85	70
Tone 3000 Hz	85	70
Tone 3150 Hz	85	70
Tone 4000 Hz	85	70
Tone 5000 Hz	70	55
Tone 6000 Hz	60	50
Tone 6300 Hz	55	45
Tone 8000 Hz	50	40
Tone 9000 Hz		
Tone 10000 Hz		
Tone 11200 Hz		
Tone 12000 Hz		
Tone 12500 Hz		
Tone 14000 Hz		
Tone 16000 Hz		
Tone 18000 Hz		
Tone 20000 Hz		

BONE VIBRATOR - NB NOISE EFFECTIVE MASKING LEVEL

Transducer	B81	B81
Impedance	10 Ω	10 Ω
Coupler	Mastoid	Forehead
	EM	EM
NB 125 Hz		
NB 160 Hz		
NB 200 Hz		
NB 250 Hz	71	83
NB 315 Hz	68	80.5
NB 400 Hz	65	78.5
NB 500 Hz	62	76
NB 630 Hz	57.5	71
NB 750 Hz	53.5	66.5
NB 800 Hz	52	64
NB 1000 Hz	48.5	57
NB 1250 Hz	45	55
NB 1500 Hz	42.5	53.5
NB 1600 Hz	41.5	52.5
NB 2000 Hz	37	48.5
NB 2500 Hz	35.5	47.5
NB 3000 Hz	36	48
NB 3150 Hz	37	48.5
NB 4000 Hz	40.5	48.5
NB 5000 Hz	45	56
NB 6000 Hz	45	56
NB 6300 Hz	45	55
NB 8000 Hz	45	55
NB 9000 Hz		
NB 10000 Hz		
NB 11200 Hz		
NB 12000 Hz		
NB 12500 Hz		
NB 14000 Hz		
NB 16000 Hz		
NB 18000 Hz		
NB 20000 Hz		
White noise	42.5	51
TEN noise		

BONE VIBRATOR - NB NOISE MAX HL

Transducer	B81	B81
Impedance	10 Ω	10 Ω
Coupler	Mastoid	Forehead
	Max HL	Max HL
NB 125 Hz		
NB 160 Hz		
NB 200 Hz		
NB 250 Hz	40	25
NB 315 Hz	50	35
NB 400 Hz	60	45
NB 500 Hz	60	45
NB 630 Hz	65	50
NB 750 Hz	65	50
NB 800 Hz	65	50
NB 1000 Hz	70	60
NB 1250 Hz	75	60
NB 1500 Hz	75	60
NB 1600 Hz	75	60
NB 2000 Hz	70	55
NB 2500 Hz	65	50
NB 3000 Hz	65	50
NB 3150 Hz	65	50
NB 4000 Hz	60	50
NB 5000 Hz	55	45
NB 6000 Hz	50	40
NB 6300 Hz	45	35
NB 8000 Hz	40	30
NB 9000 Hz		
NB 10000 Hz		
NB 11200 Hz		
NB 12000 Hz		
NB 12500 Hz		
NB 14000 Hz		
NB 16000 Hz		
NB 18000 Hz		
NB 20000 Hz		
White noise	70	60
TEN noise		

BONE VIBRATOR - ANSI SPEECH MAX HL

Transducer	B81	B81
Impedance	10 Ω	10 Ω
Coupler	Mastoid	Forehead
	Max HL	Max HL
Speech		
Speech Equ.FF.		
Speech Non-linear	60	50
Speech noise		
Speech noise Equ.FF.		
Speech noise Non-linear	50	40
White noise in speech	60	50

BONE VIBRATOR - IEC SPEECH MAX HL

Transducer	B81	B81
Impedance	10 Ω	10 Ω
Coupler	Mastoid	Forehead
	Max HL	Max HL
Speech		
Speech Equ.FF.		
Speech Non-linear	60	50
Speech noise		
Speech noise Equ.FF.		
Speech noise Non-linear	50	40
White noise in speech	60	50

FREE FIELD SPEAKERS – ANSI RETSPL AND MAX HL

Hz	Binaural RETSPL dB					Max HL dB ¹		
						Basic Speakers	GSI High Performance Speakers	GSI High Performance Speakers External Amplifier
	0° ⁴	45° ⁴	90° ⁴	135° ³	180° ³	45°	45°	45°
125	22	21.5	21	21.5	22	70	70	75
160	18	17	16.5	17	18	75	75	75
200	14.5	13.5	13	13.5	14.5	80	80	75
250	11.5	10.5	9.5	10	12	80	80	90
315	8.5	7	6	7	9	80	80	90
400	6	3.5	2.5	3.5	6.5	90	90	95
500	4.5	1.5	0	1.5	5	95	100	110
630	3	-0.5	-2	-0.5	3.5	95	100	110
750	2.5	-1	-2.5	-1.5	2.5	95	105	110
800	2	-1.5	-3	-2	1.5	95	105	110
1000	2.5	-1.5	-3	-2.5	0.5	95	105	110
1250	3.5	-0.5	-2.5	-2.5	1	95	105	110
1500	2.5	1	-2.5	-2	1	100	105	110
1600	1.5	-2	-3	-2.5	0	100	105	110
2000	-1.5	-4.5	-3.5	-2	0	100	105	115
2500	-4	-7.5	-6	-2	-1.5	100	105	115
3000	-6	-11	-8.5	-4	-3.5	105	105	115
3150	-6	-11	-8	-4	-3.5	105	105	115
4000	-5.5	-9.5	-5	-0.5	-0.5	105	105	115
5000	-1.5	-7.5	-5.5	2.5	3	100	100	105
6000	4.5	-3	-5	3	9.5	95	95	105
6300	6	-1.5	-4	2.5	11	95	95	100
8000	12.5	7	4.1	8.5	16.5	80	80	75
9000	14	9	7	10.5	18.5	85	85	85
10000	14	9.5	8	10	18	80	80	80
11200	13	9	6	8	16	80	80	80
12500	12.5	11	4.5	6.5	15	75	75	75
Speech	14.5	12.5	11	10	13	90 ²	96 ²	102 ²

- 1) When using the Calibration App, the Custom Max procedure must be utilized.
- 2) To achieve the maximum output level the crest factor per ANSI S3.6 (2018) is reduced.
- 3) The angle 135 and 180 degrees are from E.A.G Shaw and M.M Vaillancourt publication 1985
- 4) Based on ANSI S3.6-2018 / ISO 389-7

APPENDIX D - WARRANTY AND REPAIR

We, Grason-Stadler, warrant that this product is free from defects in material and workmanship and, when properly installed and used, will perform in accordance with applicable specifications. This warranty is extended to the original purchaser of the instrument by GSI through the distributor from whom it was purchased and covers defects in material and workmanship for a period of one year from date of delivery of the instrument to the original purchaser. If within one year after original shipment, it is found not to meet this standard; it will be repaired, or at our option, replaced at no charge except for transportation costs, when returned to an authorized Grason-Stadler facility. If field service is requested, there will be no charge for labor or material; however, there will be a charge for travel expense at the service center's current rate.

NOTE: Opening the instrument case or changes to the product not approved in writing by Grason-Stadler shall void this warranty. Grason-Stadler shall not be responsible for any indirect, special, or consequential damages, even if notice has been given in advance of the possibility of such damage. The transducers may go out of calibration due to rough handling or impact (dropping).

THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

We advise you against attempting to rectify any faults by yourself or commissioning non-experts to do so. Equipment is not user repairable. Repairs must be performed by an authorized service representative only.

No modifications of the equipment are allowed by anyone other than a qualified GSI representative. Modification of the equipment could be hazardous. If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.

To ensure that your instrument works properly, the GSI Pello should be checked and calibrated at least once per year. This check must be carried out by your dealer or authorized GSI service facility.

When returning the instrument for repairs or calibration it is essential to send the acoustic transducers with the device. Send the device to an authorized service center only. Please include a detailed description of faults. To prevent damage in transit, please use the original packing material, if possible, when returning the instrument.

APPENDIX E – RECYCLING AND DISPOSAL

Many local laws and regulations require special procedures to recycle or dispose of electrical equipment and related waste including batteries, printed circuit boards, electronic components, wiring and other elements of electronic devices. Follow all local laws and regulations for the proper disposal of batteries and any other parts of this system.

Below is the contact address for proper return or disposal of electronic wastes relating to Grason-Stadler products in Europe and other localities. The contact information for the WEEE in Europe:

Grason-Stadler
c/o DGS Diagnostics A/S
Audiometer Alle 1
5500 Middelfart
Denmark

Within the European Union, it is illegal to dispose of electric and electronic items in unsorted municipal waste. Electric and electronic waste may contain hazardous substances and therefore must be collected separately. Such products will be marked with the crossed-out wheeled bin symbol. The cooperation of the user is important to ensure a high level of reuse and recycling of electric and electronic waste. Failing to recycle such waste products in an appropriate way may endanger the environment and consequently the health of human beings.

APPENDIX F – ELECTROMAGNETIC COMPATIBILITY (EMC)

Portable and mobile RF communications equipment can affect the Pello. Install and operate the system according to the EMC information presented in this chapter.

The Pello has been tested for EMC emissions and immunity as a stand-alone device. Do not use the system adjacent to or stacked with other electronic equipment. If adjacent or stacked use is necessary, the user should verify normal operation in the configuration.

The use of accessories, transducers, and cables other than those specified, except for servicing parts sold by GSI as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the device.

Anyone connecting additional equipment is responsible for making sure the system complies with the IEC 60601-1-2 standard.

CAUTIONS REGARDING EMC

This instrument is suitable in hospital environments except for near active HF surgical equipment and RF shielded rooms of systems for magnetic resonance imaging, where the intensity of electromagnetic disturbance is high.

Use of this instrument adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this instrument and the other equipment should be observed to verify that they are operating normally.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Pello, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

The Pello is intended for use in the electromagnetic environment specified below. The customer or the user of the Instrument should assure that it is used in such an environment.

GUIDANCE AND MANUFACTURER'S DECLARATION
ELECTROMAGNETIC EMISSIONS

Emissions Test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B Limits	The Pello is suitable for use in all commercial, industrial, business, and residential environments.
Harmonic emissions IEC 61000-3-2	Class A Category	
Voltage fluctuations /flicker emissions IEC 61000-3-3	Complies	

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT

The Pello is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the system can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the system as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum output power of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.23\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.70	3.70	7.37
100	11.70	11.70	23.30

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1 At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2 These guidelines may not apply to all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

GUIDANCE AND MANUFACTURER'S DECLARATION

ELECTROMAGNETIC IMMUNITY

Immunity Test	IEC 60601 Test level	Compliance	Electromagnetic Environment-Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be greater than 30%.
Electrical fast transient/burst IEC61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or residential environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Not applicable	Mains power quality should be that of a typical commercial or residential environment.
Voltage dips, short interruptions, and voltage variations on power supply lines IEC 61000-4-11	$< 5\%$ <i>UT</i> ($>95\%$ dip in <i>UT</i>) for 0.5 cycle 40% <i>UT</i> (60% dip in <i>UT</i>) for 5 cycles 70% <i>UT</i> (30% dip in <i>UT</i>) for 25 cycles $<5\%$ <i>UT</i> ($>95\%$ dip in <i>UT</i>) for 5 sec	Not applicable	Mains power quality should be that of a typical commercial or residential environment. If the user of the system requires continued operation during power mains interruptions, it is recommended that the system be powered from an uninterruptable power supply or its battery.
Power frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or residential environment.

Note: *UT* is the A.C. mains voltage prior to application of the test level.

Immunity test	IEC / EN 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC / EN 61000-4-6</p> <p>Radiated RF IEC / EN 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2,5 GHz</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any parts of the system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>$d = 1,2\sqrt{P}$</p> <p>$d = 1,2\sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d = 2,3\sqrt{P}$ 800 MHz to 2,5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, (a) should be less than the compliance level in each frequency range (b)</p>

			<p>Interference may occur near equipment marked with the following symbol:</p> 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.</p>			

^(a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the system is used exceeds the applicable RF compliance level above, the system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the system.

^(b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

NOTICE: There are no deviations from the collateral standard and allowances uses

NOTICE: All necessary instruction for maintaining compliance regarding EMC can be found in the general maintenance section in this instruction. No further steps required.

Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. To ensure compliance with the EMC requirements as specified in IEC 60601-1-2, it is essential to use only the following accessories.

Item	Manufacturer	Item #
IP30 Insert Earphone	RadioEar	8101884
DD45 Headphones	RadioEar	8503851
DD450 Headphones	RadioEar	8520321
DD65 v2 Headphones	RadioEar	8510115
B81 Bone Vibrator	RadioEar	8501119

Conformance to the EMC requirements as specified in IEC 60601-1-2 is ensured if the cable types and cable lengths are as specified below:

Description	Length	Screened/Unscreened
IP30 Insert Earphone	2.5 m	Unscreened
DD45 Headphones	2.5 m	Unscreened
DD450 Headphones	2.5 m	Unscreened
DD65 v2 Headphones	2.5 m	Unscreened
B81 Bone Vibrator	2 m	Unscreened
Power Supply Cable	2.5 m	Unscreened
USB Cable	2 m	Screened

NOTICE: The use of the accessories, transducers, and cables with medical equipment/system other than this equipment may result in increased emissions or decreased immunity of the medical equipment/system.

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The screenshot displays the GSI Pello software interface for a speech test. At the top, it shows 'Patient ...' and 'Examiners ...'. The main area is divided into three sections: Channel 1 (Right), Score - SRT, and Channel 2 (Left). Channel 1 is set to 40 dB HL and Channel 2 to 0 dB HL. The score for SRT is 0%. The word list 'Playground' is selected. Below the channels are two 'Speech Test Results' tables and a 'Basic Auditory Tests - Adult: Spondee A' table. The bottom status bar shows 'Test Type', 'Word Lists', 'Word Nav', 'Aided', '5 dB Step', 'Ch2 Mic', and the date/time '12:02 PM 6/29/2016'.

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BKB-SIN (VOLITELNÁ LICENCE SPEECH PLUS)

The screenshot displays the BKB-SIN software interface. At the top, it shows 'Patient ...', 'BKB-SIN', and 'Examiner 1'. The main area is divided into two channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Below the channels are sections for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' with 'R', 'B', and 'L' indicators. There are also fields for 'INT A' and 'INT B' with 'Right' and 'Left' labels. Below these are fields for 'PTA AC: BC:' and 'SII:'. The 'Reliability' field is empty. The 'Test Results Group 1' section contains three tables for 'Right' and 'Left' ears, with columns for 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. The 'Score' section shows a list of word pairs (A1-A10) and their corresponding scores (S/N 21 to S/N 9). The bottom bar includes 'Word Lists', 'Word Nav', 'Age', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a timestamp '10:03 AM 02-26-2019'.

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QUICKSIN (VOLITELNÁ LICENCE SPEECH PLUS)

The screenshot displays the QuickSIN software interface. At the top, it shows 'Patient ...' and 'QuickSIN'. Below this are two main channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Each channel has a section for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages', with sub-sections for 'Basic', 'HFE', and 'HFE-LP' across 'R', 'B', and 'L' ears. Below the channels are 'Test Results Group 1' tables for each ear, with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. A 'QuickSin : List 1 (Track 3)' section contains several sentences for testing. To the right of these sentences is a 'Score' table with columns for 'S/N' and 'Score'. The bottom status bar includes 'Word Lists', 'Word Nav', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and the date/time '10:53 AM 03-01-2019'.

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The screenshot displays the GSI Pello software interface for a speech test. The interface is divided into several sections:

- Channel 1 (Right):** Shows a score of 40 dB HL. The word list is 'Playground'.
- Channel 2 (Left):** Shows a score of 0 dB HL. The word list is 'Speech Noise'.
- Score - SRT:** Shows a score of 0%.
- Speech Test Results:** Two empty tables for recording test results.
- Basic Auditory Tests - Adult: Spondee A:** A table of word lists for testing.

Ear	Test Type	Int Ext Mic	Word List	Aid	%	dB HL	dB EM

Ear	Test Type	Int Ext Mic	Word List	Aid	%	dB HL	dB EM

Basic Auditory Tests - Adult: Spondee A						
Playground	Daybreak	Northwest	Mushroom	Doormat	Eardrum	Iceberg
Padlock	Sunset	Duck pond	Cowboy	Inkwell	Baseball	Whitewash
Oatmeal	Greyhound	Hot dog	Mousetrap	Airplane	Headlight	Hothouse
Stairway	Woodwork	Drawbridge	Armchair	Schoolboy	Horseshoe	Railroad
Workshop	Pancake	Hardware	Toothbrush	Grandson	Birthday	Sidewalk
Farewell						

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The screenshot displays the BKB-SIN software interface. At the top, it shows 'Patient ...', 'BKB-SIN', and 'Examiner 1'. The main area is divided into two channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Below the channels are sections for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' with 'R', 'B', and 'L' input boxes. There are also fields for 'INT A', 'Phone', 'PTA AC: BC:', and 'SII:'. The 'Test Results Group 1' section contains three empty tables for recording results. The 'BKB-SIN : List Pair 1' section shows a list of word pairs with corresponding scores. The bottom status bar includes 'Word Lists', 'Word Nav', 'Age', 'Aided', 'dB Step' (5), 'Group' (1), 'Research', and a timestamp of 10:03 AM on 02-26-2019.

Score	S/N	Word List	Score
-	21	A1 THEY are LOOKING AT the CLOCK	-
-	18	A2 The CAR ENGINE is RUNNING	-
-	15	A3 CHILDREN LIKE STRAWBERRIES	-
-	12	A4 THEY are BUYING some BREAD	-
-	9	A5 The GREEN TOMATOES are SMALL	-
-		A6 HE PLAYED with his TRAIN	-
-		A7 The BAG FELL to the GROUND	-
-		A8 The BOY DID a HANDSTAND	-
-		A9 The WATER BOILED QUICKLY	-
-		A10 The MAN is PAINTING a SIGN	-
-		Sum	-

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The screenshot displays the QuickSIN software interface. At the top, it shows 'Patient ...' and 'QuickSIN'. Below this, there are two main channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. The interface includes sections for SNR Loss Averages (Group 1 and Group 2) with tables for Basic, HFE, and HFE-LP tests across Right (R), Back (B), and Left (L) ears. There are also fields for PTA AC: BC, SII, and Reliability. Below these are three tables for Test Results Group 1, each with columns for Ear, Word List, dB HL, SNR 50, and SNR Loss. At the bottom, there is a section for 'QuickSin : List 1 (Track 3)' with a list of sentences and a 'Score' table showing S/N ratios from 25 to 0 and a Sum. The bottom status bar includes 'Word Lists', 'Word Nav', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a timestamp of 10:53 AM on 03-01-2019.

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BKB-SIN (LICENCE SPEECH PLUS EN OPTION)

The screenshot displays the BKB-SIN software interface. At the top, it shows 'Patient ...', 'BKB-SIN', and 'Examiner 1'. Below this, there are two main channels: Channel 1 (red header) and Channel 2 (blue header). Channel 1 shows a presentation level of 70 dB HL, while Channel 2 shows 0 dB HL. Each channel has a frequency scale from -20 to 3 kHz. In the center, there are two boxes for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages', each with 'R', 'B', and 'L' sub-sections. Below these are fields for 'PTA AC: BC:', 'SII:', and 'Reliability'. The bottom section contains three 'Test Results Group 1' tables with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. At the very bottom, there is a 'BKB-SIN : List Pair 1' table with 10 pairs of word lists (A1-A10) and their corresponding scores (S/N 21 to -6). A control bar at the bottom includes 'Word Lists', 'Word Nav', 'Age', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a timestamp '10:03 AM 02-26-2019'.

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The screenshot displays the GSI Pello software interface for a speech test. At the top, it shows 'Patient ...' and 'Examiners ...'. The main area is divided into three sections: Channel 1 (Right), Score - SRT, and Channel 2 (Left). Channel 1 is set to 40 dB HL and Channel 2 to 0 dB HL. The score for SRT is 0%. The word list is 'Playground'. Below the channels are two 'Speech Test Results' tables and a 'Basic Auditory Tests - Adult: Spondee A' table. The 'Basic Auditory Tests' table has 'Playground' highlighted in yellow. The bottom status bar shows 'Test Type', 'Word Lists', 'Word Nav', 'Aided', '5 dB Step', 'Ch2 Mic', and the date/time '12:02 PM 6/29/2016'.

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BKB-SIN (OPCIONÁLIS SPEECH PLUS LICENC)

The screenshot displays the BKB-SIN software interface. At the top, it shows 'Patient ...', 'BKB-SIN', and 'Examiner 1'. The main area is divided into two channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Below the channels are sections for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' with 'R', 'B', and 'L' indicators. There are also fields for 'INT A' and 'INT B' with 'Right' and 'Left' labels, and 'Phone' options. A scale from -20 to 3 is visible for both channels. Below these are fields for 'PTA AC: BC:' and 'SII:'. The 'Reliability' field is empty. The 'Test Results Group 1' section contains three tables for recording results. The 'BKB-SIN : List Pair 1' section shows a list of word pairs with scores. The bottom bar includes 'Word Lists', 'Word Nav', 'Age', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a timestamp '10:03 AM 02-26-2019'.

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QUICKSIN (OPCIONÁLIS SPEECH PLUS LICENC)

The screenshot displays the QuickSIN software interface. At the top, it shows 'Patient ...' and 'QuickSIN'. Below this, there are two main channels: Channel 1 (red header) and Channel 2 (blue header). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Each channel has a 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' section with tables for 'Basic', 'HFE', and 'HFE-LP' tests across 'R', 'B', and 'L' ears. Below the channels are 'PTA AC: BC:' and 'SII:' fields for both Right and Left ears. The 'Test Results Group 1' section contains three empty tables with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. The 'QuickSin : List 1 (Track 3)' section shows a list of test sentences. To the right of this list is a 'Score' table with columns for 'S/N' and 'Score'. The bottom of the interface has a control bar with 'Word Lists', 'Word Nav', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a timestamp '10:53 AM 03-01-2019'.

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The screenshot displays the GSI Pello software interface for a speech test. At the top, it shows 'Patient ...' and 'Examiners ...'. The main area is divided into three sections: Channel 1 (Right), Score - SRT, and Channel 2 (Left). Channel 1 shows a score of 40 dB HL, and Channel 2 shows 0 dB HL. The word 'Playground' is displayed in the center. Below the channels are sections for 'Speech Test Results' (empty tables) and 'Basic Auditory Tests - Adult: Spondee A' (a grid of word lists). The bottom status bar shows 'Test Type', 'Word Lists', 'Word Nav', 'Aided', '5 dB Step', 'Ch2 Mic', and the date/time '12:02 PM 6/29/2016'.

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
The screenshot displays the BKB-SIN software interface. At the top, it shows 'Patient ...', 'BKB-SIN', and 'Examiner 1'. The interface is divided into two main channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Below the channels, there are sections for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' with input fields for Right (R), Better (B), and Left (L) ears. The interface also includes 'INT A' and 'INT B' sections for Right and Left ears, respectively, with 'Phone' and 'Reliability' options. Below these are 'Test Results Group 1' tables for each ear, with columns for Ear, Word List, dB HL, SNR 50, and SNR Loss. At the bottom, there is a 'List Pair 1' section with a table of word lists and scores, and a 'Score' section on the right. The bottom status bar shows 'Word Lists', 'Word Nav', 'Age', 'Aided', 'dB Step' (5), 'Group' (1), 'Research', and a timestamp of 10:03 AM on 02-26-2019.

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The screenshot displays the QuickSIN software interface. At the top, it shows 'Patient ...' and 'QuickSIN'. Below this, there are two main channels: Channel 1 (red header) and Channel 2 (blue header). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Each channel has a 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' section with tables for 'Basic', 'HFE', and 'HFE-LP' across 'R', 'B', and 'L' ears. Below these are 'Test Results Group 1' tables for each channel, with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. A 'QuickSin : List 1 (Track 3)' section contains several sentences for testing. To the right of this list is a 'Score' table with columns for 'S/N' and 'Score'. At the bottom, there is a control bar with buttons for 'Word Lists', 'Word Nav', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a clock showing '10:53 AM 03-01-2019'.

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The screenshot displays the 'Pure Tone' test interface. At the top, it shows 'Patient ...' and 'Examiners ...'. The main area is divided into two channels:

- Channel 1 (Right Ear):** Shows a hearing level of 40 dB HL. Below it, there are options for 'Tone Steady', 'Right', and 'Phone'.
- Channel 2 (Left Ear):** Shows a hearing level of 0 dB HL. Below it, there are options for 'NBN', 'Left', and 'Phone'.

In the center, the 'Frequency' is set to 1000 Hz. Below the channels, there are fields for 'PTA AC: BC:', 'SII:', and 'Reliability' for both Right and Left ears.

Two large graphs are shown for 'Right Ear' and 'Left Ear'. Both graphs plot 'Hearing Level (dB HL)' on the y-axis (ranging from -10 to 120) against 'Frequency (Hz)' on the x-axis (logarithmic scale: 125, 250, 500, 1K, 2K, 4K, 8K). Below each graph are checkboxes for 'AC' and 'BC' and a 'dB EM' label.

At the bottom, there is a 'Signal' menu with options: '5 dB Step', 'HTL Test Type', 'Aided', and 'Auto Hz'. There is also a 'Ch2 Tone' button and a notes icon (pencil) in the bottom right corner. The time and date '2:10 PM 6/28/2016' are displayed in the bottom right corner.

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The screenshot displays the GSI Pello software interface for a speech test. At the top, it shows 'Patient ...' and 'Examiners ...'. The main area is divided into two channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 40 dB HL and Channel 2 to 0 dB HL. The score for SRT is 0%. The word list 'Playground' is displayed. Below the channels are two 'Speech Test Results' tables and a 'Basic Auditory Tests - Adult: Spondee A' table. The 'Playground' word is highlighted in yellow in the table. The bottom status bar shows 'Test Type', 'Word Lists', 'Word Nav', 'Aided', '5 dB Step', 'Ch2 Mic', and the date/time '12:02 PM 6/29/2016'.

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BKB-SIN : List Pair 1 Page 1/2

Score	Word List	Score	Word List
S/N 21	A1 THEY are LOOKING AT the CLOCK	A6	HE PLAYED with his TRAIN
S/N 18	A2 The CAR ENGINE is RUNNING	A7	The BAG FELL to the GROUND
S/N 15	A3 CHILDREN LIKE STRAWBERRIES	A8	The BOY DID a HANDSTAND
S/N 12	A4 THEY are BUYING some BREAD	A9	The WATER BOILED QUICKLY
S/N 9	A5 The GREEN TOMATOES are SMALL	A10	The MAN is PAINTING a SIGN

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The screenshot displays the QuickSIN software interface. At the top, it shows 'Patient ...' and 'QuickSIN'. Below this, there are two main channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Each channel has a volume control slider and a 'Insert Phone' button. The interface also includes sections for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages', each with a grid for 'Basic', 'HFE', and 'HFE-LP' tests across 'R', 'B', and 'L' ears. Below these are 'Test Results Group 1' tables for each channel, with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. A 'QuickSin : List 1 (Track 3)' section contains several sentences for listening. To the right of this list is a 'Score' table with columns for 'S/N' and 'Score'. The bottom of the interface features a control bar with buttons for 'Word Lists', 'Word Nav', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a clock showing '10:53 AM 03-01-2019'.

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BKB-SIN (VOLITEL'NÁ LICENCIA SPEECH PLUS)

The screenshot displays the BKB-SIN software interface. At the top, it shows 'Patient ...', 'BKB-SIN', and 'Examiner 1'. The main area is divided into two channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Below the channels are sections for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' with input fields for Right (R), Back (B), and Left (L) ears. There are also fields for 'INT A' and 'INT B' with 'Right' and 'Left' options, and 'Phone' checkboxes. A scale from -20 to 3 is visible for both channels. Below these are fields for 'PTA AC: BC:' and 'SII:'. The 'Reliability' field is empty. The 'Test Results Group 1' section contains three empty tables with columns for Ear, Word List, dB HL, SNR 50, and SNR Loss. The 'BKB-SIN : List Pair 1' section shows a list of word pairs with scores. The 'Score' columns on either side show scores for S/N 21, 18, 15, 12, 9, 6, 3, 0, -3, -6, and a Sum. The bottom bar includes 'Word Lists', 'Word Nav', 'Age', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a timestamp '10:03 AM 02-26-2019'.

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QUICKSIN (VOLITEĽNÁ LICENCIA SPEECH PLUS)

The screenshot displays the QuickSIN software interface. At the top, it shows 'Patient ...' and 'QuickSIN'. Below this, there are two main channels: Channel 1 (red header) and Channel 2 (blue header). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Each channel has a frequency scale from -20 to 3 dB HL. In the center, there are two 'SNR Loss Averages' sections for Group 1 and Group 2, each with a grid for Basic, HFE, and HFE-LP tests across Right (R), Back (B), and Left (L) ears. Below these are 'PTA AC: BC:' and 'SII:' fields for both Right and Left ears. The bottom section contains three 'Test Results Group 1' tables with columns for Ear, Word List, dB HL, SNR 50, and SNR Loss. A 'QuickSin : List 1 (Track 3)' section shows a list of sentences for testing. To the right of this list is a 'Score' table with columns for S/N (25, 20, 15, 10, 5, 0) and Sum. At the very bottom, there is a control bar with buttons for 'Word Lists', 'Word Nav', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a clock showing 10:53 AM on 03-01-2019.

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BKB-SIN (OPCIJSKA LICENCA SPEECH PLUS)

The screenshot displays the BKB-SIN software interface. At the top, it shows 'Patient ...', 'BKB-SIN', and 'Examiner 1'. The main area is divided into two channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Below the channels are sections for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' with 'R', 'B', and 'L' indicators. There are also fields for 'INT A' and 'INT B' with 'Right' and 'Left' labels, and 'Phone' indicators. A scale from -20 to 3 is visible for both channels. Below these are fields for 'PTA AC: BC:' and 'SII:' for both Right and Left ears, along with a 'Reliability' field. The 'Test Results Group 1' section contains three empty tables with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. At the bottom, there is a 'BKB-SIN : List Pair 1' section with a table of word pairs and a 'Score' column. The interface also includes a 'Word Lists' button, 'Word Nav', 'Age', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a timestamp '10:03 AM 02-26-2019'.

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QUICKSIN (OPCIJSKA LICENCA SPEECH PLUS)

The screenshot shows the QuickSIN software interface. At the top, there's a header with 'Patient ...' and 'QuickSIN'. Below this, there are two main channels: Channel 1 (Right ear) and Channel 2 (Left ear). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Each channel has a section for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' with tables for 'Basic', 'HFE', and 'HFE-LP' across 'R', 'B', and 'L' ears. Below the channels are fields for 'PTA AC: BC:', 'SII:', and 'Reliability'. The middle section contains three 'Test Results Group 1' tables with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. Below these is a list of speech tracks under 'QuickSin : List 1 (Track 3)'. To the right of the list is a 'Score' table with columns for 'S/N' and 'Score'. The bottom section is a control panel with buttons for 'Word Lists', 'Word Nav', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a clock showing '10:53 AM 03-01-2019'.

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The screenshot displays the GSI Pello software interface for a speech test. At the top, it shows 'Patient ...' and 'Examiners ...'. The main area is divided into three sections: Channel 1 (Right), Score - SRT, and Channel 2 (Left). Channel 1 shows a score of 40 dB HL, and Channel 2 shows a score of 0 dB HL. The word 'Playground' is displayed in the center. Below the channels are two 'Speech Test Results' tables and a 'Basic Auditory Tests - Adult: Spondee A' table. The 'Basic Auditory Tests' table has 'Playground' highlighted in yellow. The bottom status bar shows 'Test Type', 'Word Lists', 'Word Nav', 'Aided', '5 dB Step', 'Ch2 Mic', and the date/time '12:02 PM 6/29/2016'.

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BKB-SIN : List Pair 1 Page 1/2

Score	S/N	Word List	Score
-	A1	THEY are LOOKING AT the CLOCK	A6 HE PLAYED with his TRAIN
-	A2	The CAR ENGINE is RUNNING	A7 The BAG FELL to the GROUND
-	A3	CHILDREN LIKE STRAWBERRIES	A8 The BOY DID a HANDSTAND
-	A4	THEY are BUYING some BREAD	A9 The WATER BOILED QUICKLY
-	A5	The GREEN TOMATOES are SMALL	A10 The MAN is PAINTING a SIGN
-	S/N 6		
-	S/N 3		
-	S/N 0		
-	S/N -3		
-	S/N -6		
-	Sum		

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The screenshot displays the QuickSIN software interface. At the top, it shows 'Patient ...' and 'QuickSIN'. Below this, there are two main channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Each channel has a volume control slider and a 'Insert Phone' button. The interface also includes sections for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages', each with a grid for 'Basic', 'HFE', and 'HFE-LP' across 'R', 'B', and 'L' ears. Below these are 'Test Results Group 1' tables for each ear, with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. A 'QuickSin : List 1 (Track 3)' section contains several sentences for the test. To the right of this list is a 'Score' table with columns for 'S/N' and 'Score'. The bottom status bar shows 'Word Lists', 'Word Nav', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and the time '10:53 AM 03-01-2019'.

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Copyright © 2021 Grason-Stadler. Všechna práva vyhrazena. Žádná část této publikace nesmí být bez předchozího písemného povolení společnosti Grason-Stadler reprodukována nebo přenášena v žádné podobě ani žádným způsobem. Informace uvedené v této publikaci jsou majetkem společnosti Grason-Stadler.
Copyright © 2021 Grason-Stadler. Με επιφύλαξη όλων των δικαιωμάτων. Κανένα μέρος της παρούσας δημοσίευσης δεν επιτρέπεται να αναπαραχθεί ή να διαβιβαστεί σε οποιαδήποτε μορφή ή με οποιοδήποτε μέσο χωρίς την προηγούμενη γραπτή άδεια της Grason-Stadler. Οι πληροφορίες σε αυτή τη δημοσίευση αποτελούν ιδιοκτησία της Grason-Stadler.
Copyright © 2021 Grason-Stadler. Tous droits réservés. Aucune partie de cette publication ne peut être reproduite ou transmise sous quelque forme ou par quelque moyen que ce soit sans l'autorisation écrite préalable de Grason-Stadler. Les informations contenues dans cette publication sont la propriété de Grason-Stadler.
Copyright © 2021 Grason-Stadler. Minden jog fenntartva. A kiadvány semmilyen formában és eszközzel nem sokszorosítható vagy továbbítható a Grason-Stadler előzetes írásbeli engedélye nélkül. A jelen kiadványban szereplő információk a Grason-Stadler tulajdonát képezik.
Copyright © 2021 Grason-Stadler. Tutti i diritti riservati. Nessuna parte di questa pubblicazione può essere riprodotta o trasmessa in alcuna forma o con alcun mezzo senza la previa autorizzazione scritta di Grason-Stadler. Le informazioni contenute nella presente pubblicazione sono di proprietà di Grason-Stadler.
Copyright © 2021 Grason-Stadler. Wszelkie prawa zastrzeżone. Żadna część niniejszej publikacji nie może być powielana ani przekazywana w jakiegokolwiek formie lub w jakikolwiek sposób bez uprzedniej pisemnej zgody firmy Grason-Stadler. Informacje zawarte w niniejszej publikacji są własnością firmy Grason-Stadler.