

OTICON | More

Give your patients the power to go

life-changing
technology

Introducing the new, easy-to-use Oticon SmartCharger.* This SmartCharger joins our current desktop charger to give patients freedom to enjoy life on the go. It features a built-in power bank, a lid for additional protection, and an automatic drying feature. The SmartCharger offers convenient, reliable charging without access to a power outlet.

Patients now have the option to purchase the sleek portable SmartCharger to give them the freedom to go wherever life takes them. The SmartCharger keeps Oticon More hearing aids powered up, dry and protected.



The new Oticon SmartCharger:

- Recharges hearing aids without access to a power outlet—minimum of three full charges
- Charges hearing aids in **three hours for a full day of power****
- **30-minute quick charge** provides an additional six hours of power
- **Designed for travel** with a protective lid and built-in power bank
- Drying function **automatically removes moisture** from hearing aids



*Oticon SmartCharger is compatible with Oticon More miniRITE R and Oticon CROS PX only.

**Lithium-ion battery performance varies depending upon hearing loss, lifestyle and streaming behavior.

Charging times and light status



Charging times

Hearing aids and SmartCharger power bank charge simultaneously.



Oticon Hearing Aids		
3 hours	1 hour	0.5 hours
Fully charged	50% charged	25% charged

Oticon SmartCharger	
4 hours	Up to 8 hours
Fully charged power bank battery	Fully charged power bank and hearing aids

Charging the SmartCharger

When the charger is connected to a power source, the LEDs on the back of the charger indicate the charging status of the power bank.

1	Charging (battery too low to recharge the hearing aid)
2	Charging
3	
4	
5	Fully charged

1 ORANGE blink 2 ORANGE 3 GREEN

Power bank function

When the charger is disconnected from power, the LEDs on the back of the charger indicate the charger's power bank battery level. Your patient can check the remaining battery level by opening the SmartCharger's lid.

5	Three full charging cycles of two hearing aids
4	Two full charging cycles of two hearing aids
3	One full charging cycles of two hearing aids
2	Battery too low for a full charging cycle

2 ORANGE blink 3 GREEN

Please Note: Oticon More instruments will need to be updated to FW 1.1.1 in order to work with the SmartCharger

21-176909 15500-0568/07.21