

# OTICON | Play PX

## Technical data sheet

### miniRITE R

60 85 100 105



	Play PX 1	Play PX 2	
Speech Understanding	MoreSound Intelligence™	Level 1	Level 3
	- Environment configuration	5 Options	3 Options
	- Virtual Outer Ear	3 Configurations	1 Configuration
	- Spatial Balancer	100%	60%
	- Neural Noise Suppression, Difficult / Easy	10 dB / 4 dB	6 dB / 0 dB
	- Sound Enhancer	3 Configurations	1 Configuration
	MoreSound Amplifier™	•	•
	Feedback Prevention	MoreSound Optimizer™ & Feedback shield	MoreSound Optimizer™ & Feedback shield
	Spatial Sound™	4 Estimators	2 Estimators
	Soft Speech Booster	•	•
Sound Quality	Frequency lowering	Speech Rescue™	Speech Rescue™
	Clear Dynamics	•	-
	Better-Ear Priority	•	-
	Fitting Bandwidth*	10 kHz	8 kHz
	Bass Boost (streaming)	•	•
Listening Comfort	Processing Channels	64	48
	Transient Noise Management	4 configurations	3 configurations
Optimizing Fitting	Wind Noise Management	•	•
	Fitting Bands	24	18
	REM Autofit	Verifit@LINK, IMC 2**	Verifit@LINK, IMC 2**
	Pediatric Fitting Mode	•	•
	DSL Fitting Range***	•	•
Designed for children	Fitting Formulas	DSL v5.0, NAL-NL 1/ NAL-NL 2, VAC+	DSL v5.0, NAL-NL 1/ NAL-NL 2, VAC+
	LED	•	•
	Biological safe	•	•
	Nano coating	•	•
	Color options	12	12
	Hands-free communication****	•	•
	Direct streaming*****	•	•
	Edumic	•	•
Oticon ON app	•	•	

\* Bandwidth accessible for gain adjustments during fitting

\*\* Inter Module Communication 2

\*\*\* Available in this Technical Data sheet and Oticon Play PX Product Guide

\*\*\*\* Available for Oticon Play PX from FW 1.1 with selected iPhone models

\*\*\*\*\* From iPhone®, iPad®, iPod touch®, and selected Android™ devices

#### Operating and charging conditions

Temperature: +5°C to +40°C (41°F to 104°F)  
Relative humidity: 5% to 93%, non-condensing  
Atmospheric pressure: 700 hPa to 1060 hPa

#### Storage and transportation conditions

Temperature and humidity should not exceed the below limits for extended periods during transportation and storage.

#### Transport

Temperature: -20°C to +60°C (-4°F to 140°F)  
Relative humidity: 5% to 93%, non-condensing  
Atmospheric pressure: 700 hPa to 1060 hPa

#### Storage

Temperature: -20°C to +30°C (-4°F to 86°F)  
Relative humidity: 5% to 93%, non-condensing  
Atmospheric pressure: 700 hPa to 1060 hPa

Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries.

Oticon Play PX miniRITE R offers a discreet design powered by a rechargeable lithium-ion battery. The style features telecoil and a double push-button. It is a Made for iPhone® hearing aid and compatible with the new Android™ protocol for Audio Streaming for Hearing Aids (ASHA) - making it possible to stream directly from iPhone, iPad®, iPod touch® and select Android devices.

MoreSound Intelligence™ creates a more precise and natural representation of individual sounds with clearer and more distinct contrasts providing access to all relevant sounds.








MoreSound Amplifier™ analyzes details in sound, and optimally amplifies them for the brain to have access to relevant information.

Oticon Play PX is built on the innovative Polaris™ platform, which uses a Deep Neural Network to rapidly and optimally manage incoming sounds based on individual needs.










For information on compatibility, please visit [www.oticon.com/support/compatibility](http://www.oticon.com/support/compatibility)



		Ear Simulator Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	2CC Coupler Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006
 <p>60</p>  <p>DSL Fitting Range</p> <p>■ Mold, Bass &amp; Power dome</p> <p>□ OpenBass dome</p> <p><b>Technical information</b> Omnidirectional mode is used unless otherwise stated.</p>		<p><b>OSPL90</b></p>  <p><b>Full-on gain</b></p>  <p><b>Frequency response</b></p>  <p>— Acoustic input: 60 dB SPL - - - Magnetic input: 31.6 mA/m</p>	<p><b>OSPL90</b></p>  <p><b>Full-on gain</b></p>  <p><b>Frequency response</b></p> 
		<p>Peak OSPL90 116 dB SPL 1600 Hz 110 dB SPL HFA-OSPL90 110 dB SPL</p> <p>Peak Full-on gain<sup>1</sup> 46 dB 1600 Hz 37 dB HFA-FOG 38 dB</p> <p>Reference test gain 31 dB</p> <p>Frequency range 100-9600 Hz</p> <p>1 mA/m field 68 dB SPL Telecoil output (1600 Hz) 10 mA/m field 88 dB SPL</p> <p>SPLITS L/R -</p> <p>500 Hz &lt;2 % Total harmonic distortion (Input 70 dB SPL) 800 Hz &lt;3 % 1600 Hz &lt;2 %</p> <p>Omni 18 dB SPL Equivalent input noise level Dir 26 dB SPL</p> <p>Battery Lithium-Ion</p> <p>Expected operating time, hours<sup>2</sup> 24</p>	<p>Peak OSPL90 106 dB SPL 1600 Hz 102 dB SPL HFA-OSPL90 103 dB SPL</p> <p>Peak Full-on gain<sup>1</sup> 36 dB 1600 Hz 29 dB HFA-FOG 30 dB</p> <p>Reference test gain 26 dB</p> <p>Frequency range 100-9400 Hz</p> <p>- Telecoil output (1600 Hz) -</p> <p>83/83 dB SPL</p> <p>&lt;2 % Total harmonic distortion (Input 70 dB SPL) &lt;2 % &lt;2 %</p> <p>17 dB SPL Equivalent input noise level 28 dB SPL</p> <p>Battery Lithium-Ion</p>

1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+A1:1994 but without influence of feedback.








2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.

		<b>Ear Simulator</b> Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	<b>2CC Coupler</b> Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006
 <p>60</p> <p>DSL Fitting Range</p> <p> <input checked="" type="checkbox"/> Mold, Bass &amp; Power dome  <input type="checkbox"/> OpenBass dome                 </p> <p><b>Technical information</b> Omnidirectional mode is used unless otherwise stated.</p>		<b>OSPL90</b> 	<b>OSPL90</b> 
		<b>Full-on gain</b> 	<b>Full-on gain</b> 
		<b>Frequency response</b>  <p>                     — Acoustic input: 60 dB SPL                      - - - Magnetic input: 31.6 mA/m                 </p>	<b>Frequency response</b> 
		<b>OSPL90</b> Peak 1600 Hz HFA-OSPL90	116 dB SPL 110 dB SPL 110 dB SPL
<b>Full-on gain<sup>1</sup></b> Peak 1600 Hz HFA-FOG	46 dB 37 dB 38 dB	36 dB 29 dB 30 dB	
<b>Reference test gain</b>	31 dB	26 dB	
<b>Frequency range</b>	100-7500 Hz	100-7500 Hz	
<b>Telecoil output (1600 Hz)</b> 1 mA/m field 10 mA/m field SPLITS L/R	68 dB SPL 88 dB SPL -	- - 83/83 dB SPL	
<b>Total harmonic distortion (Input 70 dB SPL)</b> 500 Hz 800 Hz 1600 Hz	< 2 % < 3 % < 2 %	< 2 % < 2 % < 2 %	
<b>Equivalent input noise level</b> Omni Dir	19 dB SPL 26 dB SPL	17 dB SPL 29 dB SPL	
<b>Battery</b>	Lithium-Ion	Lithium-Ion	
<b>Expected operating time, hours<sup>2</sup></b>	24		

1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+A1:1994 but without influence of feedback.  
 2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.








		<b>Ear Simulator</b> Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	<b>2CC Coupler</b> Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006
<p><b>85</b></p> <p>DSL Fitting Range</p> <p> <input checked="" type="checkbox"/> Mold, Bass &amp; Power dome  <input type="checkbox"/> OpenBass dome                 </p> <p><b>Technical information</b>                      Omnidirectional mode is used unless otherwise stated.</p>		<b>OSPL90</b> 	<b>OSPL90</b> 
		<b>Full-on gain</b> 	<b>Full-on gain</b> 
		<b>Frequency response</b> <p>                     — Acoustic input: 60 dB SPL                      - - - Magnetic input: 31.6 mA/m                 </p>	<b>Frequency response</b> 
OSPL90	Peak 1600 Hz HFA-OSPL90	127 dB SPL 121 dB SPL 122 dB SPL	117 dB SPL 113 dB SPL 114 dB SPL
Full-on gain <sup>1</sup>	Peak 1600 Hz HFA-FOG	66 dB 53 dB 56 dB	55 dB 45 dB 48 dB
Reference test gain		46 dB	37 dB
Frequency range		100-9500 Hz	100-8900 Hz
Telecoil output (1600 Hz)	1 mA/m field 10 mA/m field SPLITS L/R	84 dB SPL 104 dB SPL -	- - 94/94 dB SPL
Total harmonic distortion (Input 70 dB SPL)	500 Hz 800 Hz 1600 Hz	< 2 % < 4 % < 5 %	< 2 % < 2 % < 2 %
Equivalent input noise level	Omni Dir	21 dB SPL 29 dB SPL	18 dB SPL 28 dB SPL
Battery		Lithium-Ion	Lithium-Ion
Expected operating time, hours <sup>2</sup>		24	

1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+A1:1994 but without influence of feedback.  
 2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.

		Ear Simulator Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	2CC Coupler Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006
 <p>85</p>  <p>DSL Fitting Range</p> <p>■ Mold, Bass &amp; Power dome</p> <p>□ OpenBass dome</p> <p><b>Technical information</b> Omnidirectional mode is used unless otherwise stated.</p>		<p><b>OSPL90</b></p>  <p><b>Full-on gain</b></p>  <p><b>Frequency response</b></p>  <p>— Acoustic input: 60 dB SPL - - - Magnetic input: 31.6 mA/m</p>	<p><b>OSPL90</b></p>  <p><b>Full-on gain</b></p>  <p><b>Frequency response</b></p> 
	OSPL90	Peak 1600 Hz HFA-OSPL90	127 dB SPL 121 dB SPL 122 dB SPL
Full-on gain <sup>1</sup>	Peak 1600 Hz HFA-FOG	66 dB 53 dB 56 dB	55 dB 45 dB 48 dB
Reference test gain		46 dB	37 dB
Frequency range		100-7500 Hz	100-7500 Hz
Telecoil output (1600 Hz)	1 mA/m field 10 mA/m field SPLITS L/R	84 dB SPL 104 dB SPL -	- - 94/94 dB SPL
Total harmonic distortion (Input 70 dB SPL)	500 Hz 800 Hz 1600 Hz	< 2 % < 4 % < 5 %	< 2 % < 2 % < 2 %
Equivalent input noise level	Omni Dir	22 dB SPL 29 dB SPL	18 dB SPL 27 dB SPL
Battery		Lithium-Ion	Lithium-Ion
Expected operating time, hours <sup>2</sup>		24	








1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+A1:1994 but without influence of feedback.

2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.








		<b>Ear Simulator</b> Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	<b>2CC Coupler</b> Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006
 <p>DSL Fitting Range</p> <p>Power receiver mold, Bass &amp; Power dome</p>		<b>OSPL90</b> 	<b>OSPL90</b> 
		<b>Full-on gain</b> 	<b>Full-on gain</b> 
<b>Technical information</b> Omnidirectional mode is used unless otherwise stated.		<b>Frequency response</b> 	<b>Frequency response</b> 
<b>Warning to the hearing aid dispenser</b> The maximum output capability of the hearing aid may exceed 132 dB SPL (IEC 711). Special care should be exercised in selecting and fitting the hearing aid, as there may be risk of impairing the remaining hearing of the hearing aid user.		— Acoustic input: 60 dB SPL - - - Magnetic input: 31.6 mA/m	
OSPL90	Peak 1600 Hz HFA-OSPL90	132 dB SPL 130 dB SPL 127 dB SPL	124 dB SPL 122 dB SPL 120 dB SPL
Full-on gain <sup>1</sup>	Peak 1600 Hz HFA-FOG	66 dB 60 dB 61 dB	57 dB 52 dB 53 dB
Reference test gain		53 dB	42 dB
Frequency range		100-8900 Hz	100-7500 Hz
Telecoil output (1600 Hz)	1 mA/m field 10 mA/m field SPLITS L/R	91 dB SPL 111 dB SPL -	- - 100/100 dB SPL
Total harmonic distortion (Input 70 dB SPL)	500 Hz 800 Hz 1600 Hz	<9 % <6 % <3 %	<2 % <2 % <2 %
Equivalent input noise level	Omni Dir	17 dB SPL 26 dB SPL	16 dB SPL 28 dB SPL
Battery		Lithium-Ion	Lithium-Ion
Expected operating time, hours <sup>2</sup>		24	

1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+A1:1994 but without influence of feedback.

2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.








		<b>Ear Simulator</b> Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	<b>2CC Coupler</b> Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006
 <p>DSL Fitting Range</p> <p>Power receiver mold, Bass &amp; Power dome</p>		<b>OSPL90</b> 	<b>OSPL90</b> 
		<b>Full-on gain</b> 	<b>Full-on gain</b> 
<b>Technical information</b> Omnidirectional mode is used unless otherwise stated.		<b>Frequency response</b> 	<b>Frequency response</b> 
<b>Warning to the hearing aid dispenser</b> The maximum output capability of the hearing aid may exceed 132 dB SPL (IEC 711). Special care should be exercised in selecting and fitting the hearing aid, as there may be risk of impairing the remaining hearing of the hearing aid user.		— Acoustic input: 60 dB SPL - - - Magnetic input: 31.6 mA/m	
OSPL90	Peak 1600 Hz HFA-OSPL90	132 dB SPL 130 dB SPL 127 dB SPL	124 dB SPL 122 dB SPL 120 dB SPL
Full-on gain <sup>1</sup>	Peak 1600 Hz HFA-FOG	66 dB 60 dB 61 dB	57 dB 52 dB 53 dB
Reference test gain		53 dB	42 dB
Frequency range		100-7500 Hz	100-7500 Hz
Telecoil output (1600 Hz)	1 mA/m field	91 dB SPL	-
	10 mA/m field	111 dB SPL	-
	SPLITS L/R	-	100/100 dB SPL
Total harmonic distortion (Input 70 dB SPL)	500 Hz	< 9 %	< 2 %
	800 Hz	< 6 %	< 2 %
	1600 Hz	< 3 %	< 2 %
Equivalent input noise level	Omni	17 dB SPL	17 dB SPL
	Dir	26 dB SPL	29 dB SPL
Battery		Lithium-Ion	Lithium-Ion
Expected operating time, hours <sup>2</sup>		24	

1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+A1:1994 but without influence of feedback.  
 2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.

		Ear Simulator Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	2CC Coupler Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006
 <p>105</p>  <p>DSL Fitting Range</p> <p>Power receiver mold</p>		<p>OSPL90</p>  <p>Full-on gain</p>  <p>Frequency response</p> 	<p>OSPL90</p>  <p>Full-on gain</p>  <p>Frequency response</p> 
	<p><b>Technical information</b> Omnidirectional mode is used unless otherwise stated.</p> <p><b>Warning to the hearing aid dispenser</b> The maximum output capability of the hearing aid may exceed 132 dB SPL (IEC 711). Special care should be exercised in selecting and fitting the hearing aid, as there may be risk of impairing the remaining hearing of the hearing aid user.</p> <p>— Acoustic input: 60 dB SPL - - - Magnetic input: 31.6 mA/m</p>		
OSPL90	Peak 1600 Hz HFA-OSPL90	135 dB SPL 133 dB SPL 131 dB SPL	127 dB SPL 126 dB SPL 123 dB SPL
Full-on gain <sup>1</sup>	Peak 1600 Hz HFA-FOG	72 dB 66 dB 65 dB	64 dB 59 dB 58 dB
Reference test gain		58 dB	47 dB
Frequency range		100-9100 Hz	100-7900 Hz
Telecoil output (1600 Hz)	1 mA/m field 10 mA/m field SPLITS L/R	96 dB SPL 116 dB SPL -	- - 105/105 dB SPL
Total harmonic distortion (Input 70 dB SPL)	500 Hz 800 Hz 1600 Hz	< 2 % < 2 % < 4 %	< 2 % < 2 % < 2 %
Equivalent input noise level	Omni Dir	16 dB SPL 25 dB SPL	16 dB SPL 28 dB SPL
Battery		Lithium-Ion	Lithium-Ion
Expected operating time, hours <sup>2</sup>		24	

1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+A1:1994 but without influence of feedback.

2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.

		Ear Simulator Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010	2CC Coupler Measured according to ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006
 <p><b>105</b></p> <p>DSL Fitting Range</p> <p>Power receiver mold</p>		<p><b>OSPL90</b></p> 	<p><b>OSPL90</b></p> 
		<p><b>Full-on gain</b></p> 	<p><b>Full-on gain</b></p> 
		<p><b>Frequency response</b></p>  <p>— Acoustic input: 60 dB SPL - - - Magnetic input: 31.6 mA/m</p>	<p><b>Frequency response</b></p> 
OSPL90	Peak	135 dB SPL	127 dB SPL
	1600 Hz	133 dB SPL	126 dB SPL
	HFA-OSPL90	131 dB SPL	123 dB SPL
Full-on gain <sup>1</sup>	Peak	72 dB	64 dB
	1600 Hz	66 dB	59 dB
	HFA-FOG	65 dB	58 dB
Reference test gain		58 dB	47 dB
Frequency range		100-7500 Hz	100-7500 Hz
Telecoil output (1600 Hz)	1 mA/m field	96 dB SPL	-
	10 mA/m field	116 dB SPL	-
	SPLITS L/R	-	104/104 dB SPL
Total harmonic distortion (Input 70 dB SPL)	500 Hz	< 2 %	< 2 %
	800 Hz	< 2 %	< 2 %
	1600 Hz	< 4 %	< 2 %
Equivalent input noise level	Omn	16 dB SPL	16 dB SPL
	Dir	25 dB SPL	28 dB SPL
Battery		Lithium-Ion	Lithium-Ion
Expected operating time, hours <sup>2</sup>		24	

1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+A1:1994 but without influence of feedback.

2) Expected operating time for rechargeable battery depends on use pattern, active feature set, hearing loss, sound environment, battery age and use of wireless accessories.





**Headquarters**  
Oticon A/S  
Kongebakken 9  
DK-2765 Smørum  
Denmark



SBO Hearing A/S  
Kongebakken 9  
DK-2765 Smørum  
Denmark