

B-81

High Output Bone Transducer

#8506737

Features

- Reliable performance
- Meets industry standards
- High maximum output, low distortion
- Mechanically robust
- Sensitivity: 119 dB re.1 μ N @ 1 VRMS & 1 kHz
- Total harmonic distortion: 1.1% @ 1 VRMS & 1 kHz
- Impedance: 12.5 ohm @ 1 kHz
- Secured plug concept
- RoHS compliant



The Gold
Standard

of Audiometric Transducers

Products:

- B-81 Bone Transducer, item #8104108
- Custom Safety Cable, item #8104110
- Standard Transducer Headband, item #8011098
- Pre-assembled packaged (incl. #8104108, #8104110, #8011098), item #8506491

RadioEar B-81 Bone Transducer

The B-81 is the new audiometric Bone Conductor from RadioEar. This new and enhanced bone conductor achieves higher output levels at low frequencies with a superior distortion performance. With the B-81 it is now possible, for example, to reliably measure bone conduction thresholds up to 50 dBHL at 250 Hz.

The B-81 is based on the Balanced Electro-magnetic Separation Transducer (BEST principle), where static forces are counter-balanced so that non-linear distortion forces are reduced and maximum output levels can be increased. Furthermore, the robust mechanical construction results in a significantly improved shock resistance compared to conventional designs.

The B-81 is compatible with all standard headbands and high quality cables. It also has the capability of using a custom cable designed to securely attach the plug to the transducer body. This safety feature eliminates the possibility of accidentally unplugging the device while in use.

Audiometric Benefit

"The B-81 bone vibrator provides higher output and lower harmonic distortion compared to the B71. This may offer a significant clinical advantage. Conductive components of severe hearing losses that cannot be measured with the B71 vibrator may be measurable with the B-81. This is important for medical diagnosis and treatment, programming hearing aids, and determining cochlear implant candidacy."

Robert H. Margolis

Professor Emeritus, University of Minnesota



Technical Specifications

Technical Performance

- The maximum output of the B-81 can be increased by 5-20 dB over the whole frequency range and still not exceed 5.5% THD or 6 VRMS drive voltage in reference to the IEC-ANSI Type 1 standard values. At 250 Hz, the maximum output for B-81 (median 48 dB HL) meets the standard IEC 60645-1.
- Measurements performed on Brüel & Kjær 4930, Artificial Mastoid with static force 5.4N.
- Reference equivalent threshold force levels (RETFLs) for bone vibrators.


Compliance Standards

- IEC 606045-1:2001: Electroacoustics – Audiological equipment. Part 1: Pure-tone audiometers
- ANSI/ASA S3.6-2010 American National Standard Specification for Audiometers
- ISO 389-3:1994: Acoustics – Reference zero for the calibration of audiometric equipment. Part 3: Reference equivalent threshold force levels for pure tones and bone vibrators.


| Frequency Hz | Mastoid (dB re 1uN) | Forehead (dB re 1uN) | Forehead minus mastoid |
|--------------|---------------------|----------------------|------------------------|
| 250 | 67.0 | 79.0 | 12.0 |
| 315 | 64.0 | 76.5 | 12.5 |
| 400 | 61.0 | 74.5 | 13.5 |
| 500 | 58.0 | 72.0 | 14.0 |
| 630 | 52.5 | 66.0 | 13.5 |
| 750 | 48.5 | 61.5 | 13.0 |
| 800 | 47.0 | 59.0 | 12.0 |
| 1000 | 42.5 | 51.0 | 8.5 |
| 1250 | 39.0 | 49.0 | 10.0 |
| 1500 | 36.5 | 47.5 | 11.0 |
| 1600 | 35.5 | 46.5 | 11.0 |
| 2000 | 31.0 | 42.5 | 11.5 |
| 2500 | 29.5 | 41.5 | 12.0 |
| 3150 | 31.0 | 42.5 | 11.5 |
| 4000 | 35.5 | 43.5 | 8.0 |
| 5000 | 40.0 | 51.0 | 11.0 |
| 6000 | 40.0 | 51.0 | 11.0 |
| 6300 | 40.0 | 50.0 | 10.0 |
| 8000 | 40.0 | 50.0 | 10.0 |
| Speech | 55.0 | 63.5 | 8.5 |

Table taken from ANSI standard (ANSI/ASA S3.6-2010) levels (RETFLs) for bone vibrators.


Graph 1: Vibratory Force Level (VFL)



Graph 2: Total Harmonic Distortion (THD) at Maximum Output



Graph 3: Maximum Output HL



Dimensions

- Height – 16 mm
- Length – 31.7 mm
- Width – 18.2 mm
- Weight: 20 g

