RADIOEAR B81 HEADSET

Instructions for Use – UK/US

The B81 is a bone conductor for audiometry testing for patients in all ages. The device converts electrical signals from the audiometer into vibratory i.e., bone transmitted signals and couple these signals to the skull of the patient. The device is intended to be used by an audiologist, hearing healthcare professional, or trained technician.

Warning: Only connect to IEC 60601-1 and IEC 60601-1-2 approved devices. Do not use a defective device.

Caution: Before use, make sure the headset is calibrated to the audiometric equipment it is to be used with.

Cleaning instruction: To clean the headset, use a soft, damp cloth. The use of organic solvents and aromatic oils must be avoided. If necessary, clean the device with a lint free cloth lightly dampened in cleaning solution.

Specifications:

- Standards: IEC 60645-1:2012, ANSI S3.6-2010 and ISO 389-3:1994
- Sensitivity: 119 dB re.1 μ N @ 1 VRMS & 1 kHz
- Total harmonic distortion: 1.1% @ 1 VRMS & 1 kHz
- Impedance: 10 Ω or 12.5 Ω
- Cable length: 2 m to 2.9 m, depending on variant.

Compliant to RoHS 2011/65/EU incl. 2015/863

See www.radioear.us for range of spare parts and accessories

In case of serious incident in relation to the use of the device, the incident must immediately be reported to DGS Diagnostics A/S and the local national competent authority.

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