

# BETTER HEARING FOR ME AND YOU



NEVARA 1



**bernafon**<sup>®</sup>  
Your hearing · Our passion

# BLOSSOMING CONVERSATIONS FOR ME AND YOU



Nevara is a solid solution built on Bernafon's unique Audio Efficiency™ technology. The latest developments in powerful and effective features are no longer reserved for the upper performance categories only. That's why we call Nevara a hearing aid family for everyone.

Like a florist who has a variety of flowers, plants, and decorations to choose from, there is a variety of Nevara styles and options too. Add to that an easy fitting process and remote control capabilities, and you've got a winning combination for you and your clients.

**No matter how diverse your clients are, Nevara offers a sound solution for almost every hearing loss.**

# Hearing Aids for Me and You

At Bernafon, we design hearing aids with high performance, discreetness, and durability in mind. Nevara is no exception. With peak MPOs between 119 and 138 dB SPL, mild to severe hearing losses can all be covered.

The all-new Power BTE model is the ideal companion for people who need the power in challenging listening situations. A well-designed modern hearing aid with lots of acoustic options makes the P BTE easy to fit and comfortable to wear.

Nevara hearing instruments are compatible with the RC-N remote control for additional discreetness and comfort when changing the programs or the volume. Also, Nevara is compatible with the wireless FittingLINK programming device, which makes fitting sessions easier for you and more relaxed for your clients.



## HARDWARE

	CPx	P	N	NR	ITED	ITCD	ITC	CICP	CICx
Directionality	dual omni	dir	dir	dir	dir	dir	omni	omni	omni
Battery Size	13	13	312	312	13	312	312	10	10
Program Selector	●	●	●	●	○	○	○	○	○
Volume Control	●	●	*	*	*	*	○	*	*
IPxx Rating	IP58	IP58	IP57	IP57	–	–	–	–	–

## ACCESSORIES

RC-N Remote Control	○	○	○	○	○	○	–	○	○
FittingLINK	○	○	○	○	○	○	–	○	○
FM/DAI Adapter	○	○	–	–	–	–	–	–	–

## PERFORMANCE (EAR SIMULATOR)















OSPL 90, Peak dB SPL	137	138	128	132	129	128	128	125	119
OSPL 90, 1600 Hz dB SPL	135	134	127	131	121	119	119	117	108
Full-On Gain, Peak dB	77	72	55	73	60	60	60	57	52
Full-On Gain, 1600 Hz dB	73	67	54	69	54	48	49	49	40
Reference Test Gain	60	58	47	54	47	41	42	42	34

● standard ○ optional

\*Push button can be programmed for volume control use

# Attractive Color Combinations for Me and You

Together with your client, choose an attractive combination from the following color palette. There is something suitable for everyone with Nevara.

CPx, N, NR	P		ITEs
Base Shell Graphite	Base Shell Metallic Anthracite	Base Shell Metallic Silver	
 <p data-bbox="293 919 347 942">beige</p>	 <p data-bbox="560 919 655 942">sand beige</p>	 <p data-bbox="847 919 943 942">sand beige</p>	 <p data-bbox="1158 919 1209 942">beige</p>
 <p data-bbox="268 1164 368 1187">dark brown</p>	 <p data-bbox="547 1164 668 1187">metallic silver</p>		 <p data-bbox="1131 1164 1235 1187">light brown</p>
 <p data-bbox="268 1406 368 1430">gray brown</p>	 <p data-bbox="541 1406 671 1430">antique bronze</p>		 <p data-bbox="1115 1406 1251 1430">medium brown</p>
 <p data-bbox="300 1651 336 1674">gray</p>	 <p data-bbox="525 1651 691 1674">metallic anthracite</p>		 <p data-bbox="1131 1651 1235 1674">dark brown</p>
 <p data-bbox="236 1896 400 1919">metallic anthracite</p>			

# CUSTOM-MADE SOLUTIONS FOR ME AND YOU



# A Winning Feature Set for Me and You



Audio Efficiency™

In every Nevara hearing aid, you'll find Bernafon's reliable Audio Efficiency™ technology. Important features that enhance speech intelligibility for your clients without compromising on comfort make Nevara a good choice for those who simply want to hear better.

**Speech Cue Priority™** is one of the latest inventions engineered by Bernafon. While some of your clients may benefit from comprehensive amplification of all phonemes, others might find it difficult to access the detailed speech cues. To address these different needs, Speech Cue Priority™ provides two alternative amplification strategies. Phoneme Priority mode focuses on maximum audibility of all the fine details of the speech signal. This strategy adjusts the gain to apply the precise amount of amplification to each phoneme. Envelope Priority applies less alteration to the amplification of the speech signal, preserving the crucial envelope information.

**AFC Plus** and **ANR Plus** are highly beneficial when noise is present. AFC Plus detects and removes the acoustic feedback signals effectively. The key to this feature is speed. Your clients will notice the improvement in situations such as putting on a hat, giving someone a hug, answering the telephone, or adjusting the volume control. ANR Plus manages background noise by measuring and responding to the signal-to-noise ratio. In highly noisy situations, the system applies maximum attenuation. When speech begins to dominate noise, ANR Plus is smart enough to apply some or even no attenuation to the signal to keep the speech information clear.

**Automatic Directional** fades between the full omni and a fixed directional mode. Depending on the listening environment detected, the feature selects the best setting. This is done automatically, making it more convenient for your clients.

# Feature Overview

## NEVARA 1

### SIGNAL PROCESSING

ChannelFree™	●
Speech Cue Priority™	●
Frequency Bandwidth	8 kHz

### LISTENING COMFORT

Adaptive Noise Reduction Plus (ANR Plus)	2 ctr
Adaptive Feedback Canceller Plus (AFC Plus)	●
Soft Noise Management	2 ctr

### BINAURAL COORDINATION

VC, Program Change	●
--------------------	---

### DIRECTIONALITY CONTROLS

Fixed Omni	●
Fixed Directional	●
Automatic Directional	●

### CONVENIENCE FEATURES

VC Clicks	●
Mute Via Push Button	●
Configurable Start-Up Delay	●

### INDIVIDUALIZATION

Program Options/Memories	7/4
Data Logging	●
Language Specific Targets	●
REMfit™	●
Client Interactive	●


● standard

Since 1946, we have been passionate about developing quality hearing systems that enable people with hearing difficulties to enjoy authentic listening experiences. With Swiss Engineering, precision technology, and our commitment to individual service, we strive to exceed our customer's expectations. Our aim is to provide extra value to our partners every day. Bernafon representatives and employees in over 70 countries fulfill our vision to help people with impaired hearing communicate again without limitation.

**Manufacturer:**  
Bernafon AG  
Morgenstrasse 131  
3018 Bern  
Switzerland  
[www.bernafon.com](http://www.bernafon.com)

**Local Manufacturer  
& Distributor:**  
Bernafon Canada  
500 Trillium Drive, Unit 15  
Kitchener, ON, N2R 1A7  
[www.bernafon.ca](http://www.bernafon.ca)

**SWISS**   
Engineering

 follow us on facebook  
[www.facebook.com/bernafon](http://www.facebook.com/bernafon)

**Bernafon Companies**

Australia · Canada · China · Denmark · Finland · France · Germany · Italy · Japan · Korea · Netherlands · New Zealand · Poland · Spain · Sweden · Switzerland · Turkey · UK · USA

[www.bernafon.com](http://www.bernafon.com)

**bernafon**   
Your hearing · Our passion