

POWERFUL  
PERSPECTIVES



SUPREMIA



**bernafon**®  
Your hearing • Our passion

# Developing New Perspectives

Sounds are important to all of us, but for children they are indispensable for their auditory and cognitive development. As a parent, you would like your child to benefit from the many sights and sounds all around them. For children with a hearing loss, there are many promising solutions provided by today's technology.



Hearing is critical to speech and language development, communication, and learning. Hearing aids can help a child learn to speak clearly by helping them hear other people's voices better. This makes imitating voices easier, something that is necessary in the development of speech. They can also help a child hear their own voice better, addressing an additional barrier for natural speech development.



# POWERFUL PERSPECTIVES

There are different types of hearing aids for different degrees of hearing losses.

If your child has been diagnosed with a severe to profound hearing loss, then Bernafon's Supremia Super Power BTE (behind-the-ear) might be the ideal solution.

Supremia is designed to assist children in every facet and stage of their lives.





# Specially Designed Features for Children

Supremia Super Power BTE is Bernafon's most powerful device for a pediatric fitting.

Children's lives are full of big and small events. Supremia is highly robust and durable to keep up with their daily activities. It is also equipped with a range of components which makes it a true pediatric hearing aid.





### **Volume Control and Program Button**

The easy-to-operate volume control wheel and program button are designed in a way that small fingers can easily find them and make adjustments.

### **LED Status Light**

A small status light is located below the volume control wheel. It can be programmed to visually indicate the various settings of your child's hearing aid. This will help you to more easily and discreetly monitor the operation of the hearing aid (e.g., hearing aid ON/OFF, low battery status, program change, etc.).

### **Tamper-Resistant Battery Door**

Small children are curious. In order to prevent them from accidentally losing or even swallowing the small battery, Supremia comes with an optional tamper-resistant battery door.



# The Benefits of Modern Technology

Bernafon's unique Audio Efficiency™ technology found in Supremia focuses on both speech understanding and listening comfort at the same time.

Supremia not only amplifies the sounds to make them audible. With sophisticated technology, Supremia strives to make all the speech sounds understandable. This will support your child's speech and language development.





Supremia hearing aids are also capable of adapting to many different listening environments. At school, at a birthday party, at home, and elsewhere, your child will benefit from the powerful feature set. Listening will be easier and more enjoyable with Supremia.



# Wireless Connectivity Options

Sometimes even the best hearing aids need assistance. In the classroom, an FM adapter can be attached to Supremia to connect to the school's FM system. The teacher's voice will then be transmitted directly into your child's hearing aids. This can help your child by making it easier to hear and understand the teacher.

Bernafon also offers the SoundGate 3 streaming device. Via SoundGate 3, your child can stream sound from his or her mobile phone, computer, and other audio sources directly into their hearing aids. The SoundGate App makes it possible to easily control the hearing aids via the SoundGate 3 through their iPhone®.<sup>1)</sup>

## SoundGate 3



## SoundGate App



Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.



The possibilities are better than ever with Supremia. To learn more, speak with your hearing care professional or go to our website [www.bernafon.com](http://www.bernafon.com).

<sup>1)</sup> The SoundGate App is for use with an iPhone (SoundGate App supports iPhone 6 Plus, iPhone 6, iPhone 5s, iPhone 5c, iPhone 5, iPhone 4s, and iPhone 4). The SoundGate App requires Bernafon's SoundGate 3 (or newer versions) and wireless hearing aids.

Hearing and understanding are basic human needs. Hearing problems and the resulting social consequences represent a challenge to society. We believe in a world where people with restricted hearing can communicate again thanks to advanced technology.

Visit **[www.bernafon.com](http://www.bernafon.com)** for more information about hearing and hearing loss.

500 Trillium Drive, Unit 15  
Kitchener, Ontario  
Canada N2R 1A7  
Phone +1 519 748 6669  
Toll Free +1 800 265 8250  
[contact@bernafon.ca](mailto:contact@bernafon.ca)

**SWISS**   
Engineering

**bernafon**   
*Your hearing · Our passion*

**Berna fon Companies:** Australia · Canada · China · Denmark · Finland · France · Germany · Italy · Japan · Korea · Netherlands · New Zealand · Poland · Spain · Sweden · Switzerland · Turkey · UK · USA