

MA 25-2

Air Conduction Audiometer

The MA 25 is engineered to be a portable, tough, and reliable air conduction audiometer with full range of frequencies.

Features At a Glance

- Portable air conduction screener (125-8000 Hz)
- AC adapter or battery operated (3 AA batteries)
- Easy fingertip operation with frequency and hearing level controls
- Continuous, pulse, warble tones
- Exceptionally lightweight, only 3.5 lbs
- Soft carrying case included

Upgrade Option: MA 25e

- Automated Hughson-Westlake
- Results recalled from internal memory following test
- PC compatible with SNAP software
- Includes patient response switch & USB cable

Accurate, Dependable Engineering

Take the MA 25 anywhere with confidence. User-friendly knobs make adjusting frequency and hearing levels easy. Simple push-button operation for right/left, pulse tone and warble.

Sleek, Small Design

The super-lightweight MA 25 weighs in at only 3.5 lbs. Small enough to hold in the palm of your hand! Tucked into its carrying case (included), the MA 25 is one of our smallest, lightest portables.



Technical Data

Standards:

Audiometer	EN 60645-1/ANSI S3.6, Type 4
Safety	EN 60601-1
EMC	EN 60601-1-2

Inputs/Outputs Patient response switch USB Power Supply

Rise/Fall Times ~35 msec

Distortion 0.5% typical, 2.5% maximum

Crosstalk -100 dB minimum

On/Off Ratio Greater than 80 dB

Frequency Modulation ±5% frequency modulating rate at 5Hz

Frequency Range 125 to 8000 Hz

Hearing Level Range -10 to 100 dB HL in 1 or 5 dB steps

Power Supply 3 AA batteries/External Power Supply
5V DC 100, 100 – 240V AC

Dimensions 8.9" W x 7.1" D x 2.2" H

Weight 3.5 lbs
(with accessories and carrying case)

Stimulus Continuous, Pulse, Warble

STANDARD ACCESSORIES

- DD 45 Audiometric Headset with HB7 Headband
- Audiogram Pad
- Operation Manual
- 3 AA Batteries
- Carrying Bag
- AC Power Adapter

OPTIONAL ACCESSORIES

- Patient Resonse Switch
- Audiocup Headset
- Patch Cord

CPT CODES

- 92551: Screening, pure tone
- 92552: Threshold, air

Specifications subject to change without notice.

Accessories



Standard Accessories

Carrying Case

Audiogram Pad

AC Adapter

