



The **Bio-logic Navigator Pro AEP** is flexible, portable and adaptable, providing results you can trust.

**Designed by audiologists to grow with you and your patients' needs**


The Bio-logic Navigator Pro offers a full range of AEP functionality.

- Portable and light weight
- Multiple waveform display configurations
- Active data analysis during collection
- Continuous EEG display
- Customizable protocols and reports
- Windows® 7 compatible




# Bio-logic Navigator Pro AEP – designed by audiologists for audiologists.

One small package meeting your needs...



Multiple protocols available or create your own

Dual, scrolling displays for unlimited waveforms



## Advanced data collection and analysis features

- Automatic electrode switching – eliminating the need to manually switch electrodes during testing
- Residual noise calculation and FsP – optimizing response detection
- Waveforms saved automatically – no lost data
- Flexible data analysis – mark and compare waveforms during collection
- Collection protocols – select from several built in protocols with commonly used parameters or easily create your own
- Quantifiable reproducibility – mathematically correlate waveforms to ensure reliability
- GraphMaster – create charts to compare latencies to published age referenced normative data

## Flexible report software



- Save reports as pdf file – allowing for simple transfer to the patient’s electronic medical record (EMR)
- Customizable report templates – creating professional reports to meet the needs of your practice
- Comprehensive patient database – allowing for creation of optional folders to organize patient records

## Full AEP functionality includes:

- Click evoked ABR
- Tone burst ABR
- Electocochleography (ECoHG)\*
- Middle and late latency response (MLR, ALR)
- P300
- Electrical ABR (EABR)\*
- Vestibular Evoked Myogenic Potential (VEMP)\*\*

\*required cables available at additional cost

\*\*not available in the U.S.



## The Bio-logic Navigator Pro family of products

### AEP

Standard auditory evoked potentials including ABR, ECoChG, MLR, ALR, P300 and EABR

### Scout®

Full range of OAE test capability including DPOAE and TEOAE

### ABaer®

Newborn hearing screening system including ABR, DPOAE and TEOAE

### MASTER II

Multiple Steady-State Evoked Response (ASSR) technology

### ENoG

External SEP stimulator option supports electroneuronography

### VEMP\*

EMG rectification for Vestibular Evoked Myogenic Potential test

\*not available in the U.S.



Scout



ABaer



MASTER II



ENoG



VEMP



**Special Olympics**

Natus Medical incorporated  
is an Official Healthy Hearing  
Sponsor of Special Olympics

# natus

**Natus Medical Incorporated**  
1501 Industrial Road  
San Carlos, CA 94070 USA  
1-800-303-0306  
+1-650-802-0400  
[www.natus.com](http://www.natus.com)