



## Why do professionals choose **Bio-logic Scout Sport** as their PC-based OAE system of choice?


Scout Sport's Windows®-based software features easy-to-use menus with a library of diagnostic OAE protocols:

- Diagnostic and screening DP and TE protocols
- Customizable DP protocols with a frequency range of 500 Hz - 10,000 Hz
- Diagnostic DP protocols meet CPT requirements
- Input/Output functionality
- Multiple choices for normative data comparison, included in the software
- Customizable report templates
- User interface common to the Bio-logic family of products

Scout Software can be added to a Bio-logic Navigator Pro device for added functionality of your AEP system.



Scout OAE software adds flexibility and complements existing Bio-logic brand hearing products such as the Bio-logic Navigator® Pro AEP diagnostic system.



DPOAE

TEOAE

I/O Function

### Is an OAE Test Reimbursable?

There are three CPT codes for OAE testing: 92587 (Screening otoacoustic emissions; limited), 92558 (Automated pass/refer screening) and 92588 (Evoked otoacoustic emissions; comprehensive or diagnostic evaluation).

Scout software is compatible with Windows® 7

**natus.**  
**Natus Medical Incorporated**  
 1501 Industrial Road  
 San Carlos, CA 94070 USA  
 1-800-303-0306  
 +1-650-802-0400  
[www.natus.com](http://www.natus.com)

© 2013 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.