

WIRELESS VIDEO VRA

(VDS-1500 INFRA-RED)



VDS-2500 SPEAKER/CONSOLE



CENTERING
MONITOR VDS-510

- Infra-Red Controlled with Multifunction Remote
- Control Image, Sound, LED Lights, Center Box
- Multiple Systems without Cross Interference
- Selectable "ON" Time Interval
- Dual Remote IR Sensors Included
- Optional Center Monitor or Animated Toy Ctr
- Optional VIDEO VRA Console/Speakers

www.pehratek.com



PEHRATEK PRODUCTS 9285 Pierson Lake Drive, Chaska, MN


pehratek

THE WIRELESS VIDEO VRA-RF



This combination remote control permits the user to activate the right/left/and center box (optional) functions from one remote control.

USE ENCLOSED LCD *THIN* REMOTE TO MUTE SOUND OR ADJUST VOLUME ONLY



*LED blinks when infra-red signals are detected





DUAL ARM MOUNTING BRACKETS



POWER SUPPLY

14 1/4"



11"

2 1/4"



www.pehratek.com

SYSTEM CONTENTS

- 2- 15" VIDEO LCD Monitors
- 1-Wireless IR Remote Control
(Uses 2 AA Battery)
- 2-Power Packs 12Volts DC
- 2-Dual-Arm Mounting Brackets
- 1-Instruction Manual

OPTIONAL ACCESSORIES

- Centering Monitor/Toy Box
- Rolling Stand
- LCD Centering Monitor (VDS-510)