XM-5000 Animated VRA

WIRELESS INFRA-RED or WIRED





- 10 Selectable Sounds
- 6 Selectable LED Patterns
- 4 Selectable Movement Patterns
- Infra-Red Wireless Control
- Corner/Dual Arm Mounting





FEATURES

- Two Animated 3D Toys with Movements, Sounds, Lights
- Bright 40 LED Selectable Light Patterns /10 Selectable Sounds
- Infra-Red Remote Controlled –No interference
- Wired or Wireless Operation-Field Selectable
- QUAD Operation 4 Boxes with Wired Kit
- Compact Sized Enclosure Fits Flat to Wall (11"-11"-6")
- Hi Intensity Backlights Brightly Illuminate Toys
- No Interconnecting Cable-Independent Operation
- Selectable ON-TIME Interval
- Optional Corner Mounting Bracket-Installation Kit
- Optional Center Boxes (controlled from same remote)

www.pehratek.com

pehratek

THE WIRELESS XM-5000 VRA

SOUND

Activates Sound

CYCLE TIME

Provides a 2 second on time

LEFT

Activates Left



LIGHTS

Turns on Back

ACTION

Turns on toy

RIGHT

Activated right

CENTER

Activate Center Box

CONNECTOR BLOCK

BATTERY COMPARTMENT

(UNDERSIDE) 2-AA BATTERIES



LEGEND

Volume UP Button Volume DOWN Button

POWER-Coax Silver

IR SENSOR-BLUE

CENTER BOX -RED

WIRED OUTPUT-BLACK

WIRED INPUT-BLACK





SYSTEM CONTENTS


- 2- Reinforcer Plexi Boxes
- 1-Wireless Remote Control (Uses two AA Batteries)
- 2- Power Packs 9Volts DC
- 2-Remote IR Sensors/Cables
- 1-Instruction Manual


OPTIONAL ACESSORIES

Mounting Brackets Corner Mounting Bracket Centering Boxes Wired Conversion Kit

11" WIDE







11" Pre-Drilled Mounting Holes on Loudspeaker Back of Box

TYPICAL SOUND ROOM INSTALLATIONS

www.pehratek.co