

Science **made** smarter

# StaticVR

Rehabilitation  
**made** objective

Static  
posturography  
force plates for  
functional balance  
assessments and  
training



  
**Interacoustics**

Audiometry

Tympanometry

ABR

OAE

Hearing Aid Fitting

**Balance**  




# Rehabilitate with objective balance control data

With the StaticVR posturography force plates, you collect real-time data on your patient's balance control during assessments and training, adding another layer of performance data to monitor your patient's progress and tailor the intensity to your patient.

## Expanding your assessment capabilities

The StaticVR offers a variety of functional

balance assessments such as CTSIB and Limits of Stability (LoS). This provides a solid foundation for functional balance assessments, with crucial balance control data delivered by the posturography force plates.

## Personalized training experiences




With the two independent force plates of the StaticVR and training modules such as LoS Rehab and BirdVR, you have

the opportunity to deliver personalized training while collecting real-time objective performance data.

## Tailored to your needs

As your balance rehabilitation program grows, StaticVR can grow alongside you. You can expand your solution with the PhysioVR or BalanceVR software for more comprehensive training and assessment when the need arises.

What do you get?



# Tailor your balance rehabilitation solutions to your patients' needs

BalanceVR

PhysioVR



## The power of rehabilitation with virtual reality

- Fully immersive, no visual reference in peripheral view.
- Simulated real-life training modules promoting adaptation, substitution and habituation.



Module examples:  
SVV, Target Tracking,  
Optokinetics, Supermarket


StaticVR

BalanceVR Smart

PhysioVR Smart

## ...Or with a StaticVR force plate

- Static posturography force plates for basic functional balance assessment.
- Force plates with objective data on balance control for tailored balance training.



Module examples:  
CTSIB, LOS, LOS Rehab,  
BirdVR



MotionVR (+) (OH)

BalanceVR Premium (+) (OH)

PhysioVR Premium (+) (OH)

## ...Combined with a MotionVR dynamic force plate

- CDP (Computerized Dynamic Posturography) for full functional balance assessment.
- Dynamic 360 degree force plate to simulate real-life surface and stimulate the full balance system, including the otoliths.



Module examples:  
SOT, ADT, MCT, LOS,  
Motion Program, SkiVR

SOT: Sensory Organization Test, ADT: Adaptation Test, MCT: Motor Control Test, LOS: Limits of Stability, CTSIB: Clinical Test for Sensory Interaction on Balance, SVV: Subjective Visual Vertical.

# Science made smarter

## Interacoustics is more than state-of-the-art solutions

Our mission is clear. We want to lead the way in audiology and balance by translating complexity into clarity:

- Challenges made into clear solutions
- Knowledge made practical
- Invisible medical conditions made tangible and treatable

Our advanced technology and sophisticated solutions ease the lives of healthcare professionals.

We will continue to set the standard for an entire industry. Not for the sake of science. But for the sake of enabling professionals to provide excellent treatment for their millions of patients across the globe.

[Interacoustics.com](http://Interacoustics.com)

**Interacoustics USA**  
10393 West 70th Street  
Eden Prairie, MN 55344

T +1 800 947 6334  
F +1 952 903 4200

[info@interacoustics-us.com](mailto:info@interacoustics-us.com)  
[interacoustics.com](http://interacoustics.com)



Rehabilitate balance disorders with the Virtualis solutions, which are a part of the Interacoustics balance portfolio.

Go online to  
explore our  
full product  
range

## Related products



**EyeSeeCam vHIT**  
Video Head Implants Test



**VisualEyes™ 525**  
Complete VNG solution for  
balance assessment



**TRV Chair**  
Diagnosing and treating  
Benign Paroxysmal  
Positional Vertigo (BPPV)

### Product specifications

All technical and hardware specifications concerning all products can be downloaded from our website.

