



SEAMLESS AUDIOMETER OPTIONS

AUDIOSTAR PRO VS. PELLO COMPARISON CHART

EXAM OPTIONS	AUDIOSTAR PRO	PELLO
Full Routing Flexibility	✓	
QuickSIN/BKB-SIN (Basic)	✓	Speech Plus License
QuickSIN/BKB-SIN (Basic + Research)	✓	
AZ BIO Adult Sentences and Pediatric	✓	Speech Plus License
ABLB/SISI/TEN/Tone Decay	✓	Special Test License
Pediatric Noise	✓	Special Test License
Remote Keyboard Shortcuts	✓	Special Test License
High Frequency Range (8-20 kHz)	HF Headphones Required	High Hz License
Fine Hz Functionality	✓	High Hz License
GSI AMTAS	GSI AMTAS License	GSI AMTAS License
Ototoxic Monitoring	✓	High Hz License
CAPD	✓	✓
Cochlear Implant Evaluations	✓	✓
Tinnitus Evaluations	✓	High Hz License

ROUTING ABILITY**AUDIOSTAR PRO****PELLO**

Channel 1 Narrow Band Noise	✓	
Same Ear Speech Stimuli and Speech Noise	✓	
Same Ear Tone and Narrowband Noise	✓	
Channel 1 Tone/Channel 2 Speech	✓	
Channel 1 Speech/Channel 2 Tone	✓	
Masking Stimuli in Channel 1	✓	
Binaural Masking	✓	

HARDWARE DIFFERENCES**AUDIOSTAR PRO****PELLO**

Asst Headset Monitor/Aux Intercom	✓	
Internal Amplifier	✓	
1 Calibrated Headphones + Inserts	Three Connections	Two Connections
2 Calibrated Headphones + Inserts	Three Connections	
Embedded VRA Controls	✓	
HDMI Port	✓	
Monitor Size	6 ¼" x 8 ¼"	5" x 6 ¾"
Presentation Bar	2	1
Store Button	2	1
Monitor Speaker (Internal)	✓	
Gooseneck Microphone	✓	
Dedicated Examiner Button	✓	

SOUND FIELD REQUIREMENTS**AUDIOSTAR PRO****PELLO**

SOUND FIELD REQUIREMENTS	AUDIOSTAR PRO	PELLO
90 dB HL	Basic Speakers	Amplified Speakers
96 dB HL	High Performance Speakers Plus Internal Amp	High Performance Speakers Plus External Amp
102 dB HL	High Performance Speakers Plus External Amp	High Performance Speakers Plus External Amp
2 Speakers	✓	✓
4 Speakers	✓	

	AudioStar Pro	Pello Basic	Speech Plus	Special Tests	High Hz
Air + Bone + Speech	✓	✓	✓	✓	✓
Speech Stenger	✓	✓	✓	✓	✓
Pure Tone Stenger	✓	✓	✓	✓	✓
Standard Word Lists	✓	✓	✓	✓	✓
QuickSIN Basic	✓		✓		
BKB-SIN Basic	✓		✓		
QuickSIN Advanced/Research	✓				
BKB-SIN Advanced/Research	✓				
AZ Bio Adult/Pediatric	✓		✓		
TEN Test	✓			✓	
Pediatric Noise	✓			✓	
Keyboard Shortcuts	✓			✓	
ABLB	✓			✓	
SISI	✓			✓	
Tone Decay	✓			✓	
High/Full Hz Audiometry	✓				✓
1 Hz Steps	✓				✓
Fine Hz Functionality	✓				✓
Gooseneck Microphone	✓				
Monitor Speaker (Internal)	✓				
HDMI Port	✓				
Aux Intercom	✓				
Internal Amplifier	✓				
2 Calibrated Headphones + Inserts	✓				
VRA Controls	✓				
Full Routing Capabilities	✓				