



GSI OAE DEVICE COMPARISON CHART

FEATURE	AUDERA PRO	CORTI	NOVUS
Dimensions	12 in x 15 in x 3 in	2.8 in x 1.3 in x 7 in	3.3 in x 6.2 in x 0.8 in
Weight	4.4 lbs.	0.4 lbs.	0.6 lbs.
Handheld		✓	✓
Device Type	Clinical	Diagnostic/Screening	Screening
DPOAE	✓	Per License	Per License
TEOAE	✓	Per License	Per License
TEOAE Stimulus Type	Click	Click	Click
DPOAE Levels	0 - 80 dB SPL	40 - 70 dB SPL	65/55 dB SPL
TEOAE Stimulus Range	250 - 5000 Hz	700 - 4000 Hz	1000 - 4000 Hz
DPOAE Stimulus Range	500 - 12000 Hz	1500 - 12000 Hz	1500 - 6000 Hz
Configurable Protocols	✓	✓	
Memory	Unlimited, PC-Based	500 Tests on Unit	500 Tests on Unit
Power	AC-Power	Battery, 20 hours	Battery, 150 OAE Screens
Evoked Potentials Capabilities	Comprehensive		Screening Only

FEATURE	AUDERA PRO	CORTI	NOVUS
GSI Suite Compatible		✓	
Other Data Management	Audera Pro Database	Corti Data Manager, HearSIM	HearSIM

IDEAL ENVIRONMENT	AUDERA PRO	CORTI	NOVUS
Private Practice	✓	✓	
ENT	✓	✓	
Hospital	✓	✓	✓
University	✓	✓	
Newborn Screening Program		✓	✓
Schools		✓	
VA	✓	✓	



DESIGNED SMART. BUILT STRONG.

800-700-2282 | www.grason-stadler.com