

COMPREHENSIVE PATIENT DATA MANAGEMENT

TECHNICAL SPECIFICATIONS

MINIMUM SYSTEM REQUIREMENTS

- CPU: 1.6 GHz
- RAM: 2 GB
- Hard Disk: 20 GB free space
- USB Ports: 1 for each connected GSI instrument
- Minimum Screen Resolution: 1024 x 768
- Network LAN Card (if connecting to database via network)

COMPATIBLE PRODUCTS

- AudioStar Pro
- Pello
- GSI 39
- TympStar Pro
- Corti
- Select GSI Legacy Products

SUPPORTED OPERATING SYSTEMS EDITIONS

- Windows® 10 Pro

LANGUAGES

- English, Simplified Chinese, French, German, Italian, Russian, Spanish, Thai

OTHER REQUIRED APPLICATIONS

- .NET 4.6 automatically installed

OPTIONAL APPLICATIONS

- Otoaccess™
- Noah 4.8 or higher