



# RELIABLE CLINICAL TYMPANOMETER

## TYMPSTAR PRO **NEW WIDEBAND**

The standard TympStar Pro includes all the test types a busy clinician expects. WideBand Tympanometry, along with Multi-Hz and ARLT, may be added with the advanced license.

TEST TYPE	STANDARD	ADVANCED
Screening Tymps	✓	✓
Screening Reflexes	✓	✓
Diagnostic Tymps	✓	✓
226, 678, 1000 Hz Probe Tones	✓	✓
Acoustic Reflex Thresholds	✓	✓
Acoustic Reflex Decay	✓	✓
ETF - Intact	✓	✓
ETF - Perforated	✓	✓
Patulous ET Test	✓	✓
ARLT		✓
Multi-Hz Test		✓
WideBand Tympanometry		✓



DESIGNED SMART. BUILT STRONG.

800-700-2282 | [www.grason-stadler.com](http://www.grason-stadler.com)