

Audiological Diagnostic & Video Imaging Equipment



AVANT™ A2D+ Audiometer

Powerful PC-based Dual Channel Diagnostic Audiometer

- AC, BC, Speech, Masking & Free Field Testing
- Integrated Speech Audiometry Wordlists
- Optional Fully Automated Audiometry
- Optional QuickSIN™ Speech In Noise Testing
- Powerful 3rd Party Counseling Tools
- Small Footprint, Approx 5" x 5" x 1.25"

AVANT REM Speech+

Complete Binaural Live Speech Mapping, Real Ear Measurement System

- Dual Sure Probe™ Microphone System
- USB Powered
- Results Can Be Viewed In Either LTASS Or Peak Hold Mode
- Features DSL (I/O), NAL-NL1 And NAL-NL2 Targets
- Hearing Loss Simulator For 3rd Party Demonstrations
- Small Footprint, Approx 5" x 5" x 1.25"



AVANT REMsp

Monaural Live Speech Mapping, Real Ear Measurement System

- USB Powered
- Results Can Be Viewed In Either LTASS Or Peak Hold Mode
- Features DSL (I/O), NAL-NL1 And NAL-NL2 Targets
- Hearing Loss Simulator For 3rd Party Demonstrations
- Lightweight, Very Small Size, Approx 1" x 3" x 0.5"

AVANT "Stealth" Audiometer

PC-Based Clinical High-Frequency Audiometer

- AC, BC, Speech, Masking And Free Field Testing
- Designated Transducer Outputs
- Dual Internal Amplifiers
- Powerful 3rd Party Counseling Tools
- Optional Fully Automated Audiometry
- Optional QuickSIN™ Speech In Noise Testing
- Optional High-Frequency Audiometry



Powerful PC-based Instrumentation for Audiometry, Real Ear Measurement, Live Speech Mapping, Hearing Instrument Testing & Evaluation and Video Otoscopy



Otowave Tympanometer

Handheld Portable Tympanometer

- Compatible With All MedRx Software Modules Via An Infrared Link
- Programmable Reflex Testing
- Fast And Accurate Middle Ear Measurement

AVANT Polar HIT+ Chamber

Unique USB Powered Directional Hearing Instrument Test Box

- Conducts 11 selectable Automated ANSI & IEC Hearing Instrument Tests
- Acoustic Isolation For Accurate Testing
- Provides A Motor Driven Rotating Test Pad To Automatically Generate Polar Plots
- Omni And Directional Microphone Verification



Video Companion

Bringing Video Otoscopy To Your Computer

- Superior Image Quality
- Battery Operated LED Or Halogen Light Source
- Compact, Color Video Camera
- Focus Always Feature
- Optional Trackball Mouse Including Operator Foot-Pedal

Ultra Vac

Hearing Instrument Restoration System, Includes:

- Pulse / Constant Power Feature
- Durable Vacuum Wand
- Vacuum Drying Chamber
- Selectable Wand (Vacuum Or Pressure)



(MedRx Products Are NOAH™ and TIMS® Compatible)

MedRx

www.medrx-usa.com • Toll Free 888-392-1234

MedRx

www.medrx-usa.com • Toll Free 888-392-1234