

MedRx software Windows 11 compatibility overview

February 2022 status

Product	Current Version Software Available	Compatible with WIN 11 from beginning	Compatible updating WIN 10 to 11
A2D+	Avant AUD 3.4	Yes	Yes
	Studio Version 1.1	Yes	Yes
Stealth	Avant AUD 3.4	Yes	Yes
	Studio Version 1.1	Yes	Yes
ARC	Avant AUD & REM 3.4	Yes	Yes
	Studio Version 1.1	Yes	Yes
AIR	Avant AUD 3.4	Yes	Yes
	Studio Version 1.1	Yes	Yes
Tinnometer	Tinnometer 1.3	Yes	Yes
REM+	Avant REM 3.4	Yes	Yes
	Studio Version 1.1	Yes	Yes
HIT+	Avant HIT 3.4	Yes	Yes
	Studio Version 1.1	Yes	Yes
Video Otoscope	Avant 3.4	Yes	Yes
	Studio Version 1.1	Yes	Yes
KIOSK	Studio Version 1.1	Yes	Yes
VRA	Studio Version 1.1	Yes	Yes
TYMP	Avant 3.4	Yes	Yes
	Studio Version 1.1	Yes	Yes