

**MedRx**<sup>®</sup>

www.medrx-usa.com



Ultra Vac+

## Vacuum & Drying Chamber

Repair, Restore and Maintain  
Hearing Instruments

**MedRx**<sup>®</sup>

**888-392-1234**

To Schedule An Online Demo



# UltraVac+

## Technical Specs

**Power Requirements:** 100 – 240 V, 50/60 Hz

**Power Consumption:** <100 Watts

**Fuse:** 2.5A @ 100-240V

**Operating Temperature:** 32° to 120° F • 0° to 50° C

**Weight:** Approx 11 lbs • Approx 5 kg

**Dimensions:** Approx 11.25" x 12" x 6.25" (L x W x H) Approx 29 cm x 31 cm x 16 cm (L x W x H)

**Debris Filter:** I-63S

**Vacuum Tips:** #14, #15, 2-#18 and 2-#20 gauge with Safety Lock attachment

**Drying Chamber Cycle Time:** 4 minutes with auto shut-off

**Vacuum Wand:** Continuous Use

**Pressure Wand:** Continuous Use

**Standard Accessories:** Drying Chamber Cup & Filter, Vacuum Wand Filter with O-Ring, 6 Vacuum Wand Tips, Cleaning Tool Kit, Instruction Manual

## Functions of your UltraVac+

**The Drying Chamber** – Displaces moisture trapped within the hearing instruments. You can place two ITEs, or one BTE with an earmold, into the chamber at one time.

**The Vacuum Wand** – Used to extract ear wax and debris from the receiver tubes, microphone tubes, vents and battery compartments of the hearing instruments.

**The Pressure Function** – Used to clear obstructed Vacuum Wand Tips and to blow debris from the receiver tubes, microphone tubes, vents and battery compartments of the hearing instruments.

*The Vacuum & Pressure Wands work simultaneously.*

# MedRx<sup>®</sup>

Good Things Come In Small Packages

### MedRx, Inc

1200 Starkey Rd. #105

Largo, FL 33771 U.S.A.

888-392-1234 • 727-584-9600

FAX: 727-584-9602

Email: [sales@medrx-usa.com](mailto:sales@medrx-usa.com)

[www.medrx-usa.com](http://www.medrx-usa.com)

### MedRx Minimum Computer Specs:

Windows<sup>®</sup> PC compatible computer, Intel<sup>™</sup> i5 Dual Core, 2.0 GHz or better. 4 GB RAM. 20 GB free hard drive space. Available 2.0 USB port. Windows 7, 8 or 10 Professional (32 or 64-bit).



MX-UV-DS-1

