

**MedRx**<sup>®</sup>

www.medrx-usa.com



**AVANT AIR+**  
Screening Audiometer

## Ideal Audiometer for Traveling

Use for Community Screenings where  
Only Air Thresholds Are Needed

**MedRx**<sup>®</sup>

**888-392-1234**

To Schedule An Online Demo



# AVANT AIR+

## Technical Specs

**Standards:** Screening Audiometer as per ANSI S3.6-2010, IEC 60645-1:2012, Type 4; Tone Audiometry

**Outputs:** DD65 Headphones

**Frequency Range:** Air: 125 Hz – 8 kHz

**Max Output:** Air Conduction: 100 dB HL For Mid-Range Frequencies

**Attenuation:** 1 dB Step Or 5 dB Step, User Selectable

**Optional Accessories:** DD450 Headphones, DD45 Headphones & Insert Earphones

**Compatible with:** NOAH™, TIMS® and SycleNet™

**Power Requirements:** USB-powered

**Dimensions:** Approx. 5" x 5" x 1.25" (L x W x H)

**Weight:** < 1 lb.

**Standard Accessories:** DD65 Headphones, Patient Response Switch, Talkback Microphone, Operator Mic, Monitor Headset, Software & Manuals, Carrying Case

# MedRx®

Good Things Come In Small Packages

### MedRx, Inc

1200 Starkey Rd. #105

Largo, FL 33771 U.S.A.

888-392-1234 • 727-584-9600

FAX: 727-584-9602

Email: [sales@medrx-usa.com](mailto:sales@medrx-usa.com)

[www.medrx-usa.com](http://www.medrx-usa.com)

### MedRx Minimum Computer Specs:

Windows® PC compatible computer, Intel™ i5 Dual Core, 2.0 GHz or better. 4 GB RAM. 20 GB free hard drive space. Available 2.0 USB port. Windows 7, 8 or 10 Professional (32 or 64-bit).



AVANT AIR+ comes standard with the DD65 Headset

MX-AIR-DS-1

