

# MedRx<sup>®</sup>

www.medrx-diagnostics.com



**AVANT Stealth**  
Clinical Audiometer


## The Full Clinical Solution

Allows Testing Up To 20 kHz  
with High Frequency Option.

# MedRx<sup>®</sup>

## 888-392-1234

To Schedule An Online Demo



# AVANT Stealth

## Technical Specs

**Standards:** 2-Channel Clinical Audiometer As Per ANSI S3.6-2018, IEC 60645-1:2012, IEC 60645-2:1993, IEC 60645-4:1994, Type 1 HFAE; Tone Audiometry, Speech Audiometry, Stenger Test, QuickSIN™, ABLB, SISI, Tone Decay, Hughson Westlake Automated Audiometry

**Outputs:** Insert Earphones, Headphones, Bone Conductor, Free Field Via High Power Internal Amplifiers, 2 x 20 Watts Into 4 Ohms

**Frequency Range:** Air: 125 Hz – 8000 Hz, Bone: 250 Hz – 8000 Hz, Optional: High Frequency Range With High Frequency Headphones: 8000 Hz – 20,000 Hz

**Max Output:** Air Conduction: 120 dB HL For Mid-Range Frequencies, Bone Conduction: 70 dB HL, Sound Field: 95 dB HL (Depends On Speakers)

**Attenuation:** 1 dB Step Or 5 dB Step, User Selectable

**Speech Input:** Live Microphone, MP3/Wave Files, CD

**Masking Signals Tone Audiometry:** Narrow Band Noise (Default), Speech Weighted Noise, White Noise

**Speech Audiometry:** Speech Weighted Noise (Default), White Noise, CD/File, Opposite Channel

**Hearing Loss Simulator And Hearing Instrument Simulator:** Frequency Range: 125 Hz – 8000 Hz, 13 Band Equalizer

**Communication Port:** USB 2.0 (Backward Compatible With 1.1) Compatible With: NOAH™ And TIMS®

**Power Requirements:**

USB-Powered Or External Power DC 15 V/2A

**Power Supply:** 100 V – 240 V, 50/60 Hz

**Operating Temperature:** 10°C To 35°C

**Operating Humidity:** 30% To 90%

**Storage Temperature:** -20°C To 50°C

**Storage Humidity:** 10% To 90%

**Dimensions:**

Approx. 20 cm x 12 cm x 3 cm (L x W x H)

Approx. 8" x 5" x 1.25" (L x W x H)

**Weight:** < 1 kg < 2 lbs

**Standard Accessories:** Transducers: DD65v2, DD45 Or IP30, Bone Conductor, Operator Mic / Monitor Headset, Patient Response Switch, Talkback Microphone, External Power Supply, Auditec Sound File License, QuickSIN License, USB Cable, Software & Manuals, Carrying Case

**Optional Accessories:** DD450 Headphones

**MedRx Minimum Computer Specs:**

Windows® PC Compatible Computer, Intel™ i5 Dual Core, 2.0 GHz Or Better. 4 GB RAM. 20 GB Free Hard Drive Space. Available 2.0 USB Port. Windows 7, 8 Or 10 Professional (32 Or 64-Bit).



Good Things Come In Small Packages

**MedRx, Inc**

1200 Starkey Rd. #105

Largo, FL 33771 U.S.A.

888-392-1234 • 727-584-9600

Fax: 727-584-9602

Email: sales@medrx-diagnostics.com

www.medrx-diagnostics.com



Stealth Audiometer the Full Clinical Solution.

