

MedRx[®]

www.medrx-diagnostics.com



AVANT Stealth
Clinical Audiometer


The Full Clinical Solution

Allows Testing Up To 20 kHz
with High Frequency Option.

MedRx[®]

888-392-1234

To Schedule An Online Demo



AVANT Stealth

Technical Specs

Standards: 2-Channel Clinical Audiometer As Per ANSI S3.6, IEC 60645-1, IEC 60645-2, IEC 60645-4, Type 1 HFAE; Tone Audiometry, Speech Audiometry, Stenger Test, QuickSIN™, ABLB, SISI, Tone Decay, Hughson Westlake Automated Audiometry

Outputs: Insert Earphones, Headphones, Bone Conductor, Free Field Via High Power Internal Amplifiers, 2 x 20 Watts Into 4 Ohms

Frequency Range: Air: 125 Hz – 8000 Hz, Bone: 250 Hz – 8000 Hz, Optional: High Frequency Range With High Frequency Headphones: 8000 Hz – 20,000 Hz

Max Output: Air Conduction: 120 dB HL For Mid-Range Frequencies, Bone Conduction: 70 dB HL, Sound Field: 95 dB HL (Depends On Speakers)

Attenuation: 1 dB Step Or 5 dB Step, User Selectable

Speech Input: Live Microphone, MP3/Wave Files, CD

Masking Signals Tone Audiometry: Narrow Band Noise (Default), Speech Weighted Noise, White Noise

Speech Audiometry: Speech Weighted Noise (Default), White Noise, CD/File, Opposite Channel

Hearing Loss Simulator And Hearing Instrument Simulator: Frequency Range: 125 Hz – 8000 Hz, 13 Band Equalizer

Communication Port: USB 2.0 (Backward Compatible With 1.1) Compatible With: Noah And TIMS

Power Requirements:

USB-Powered Or External Power DC 15 V/2A

Power Supply: 100 V – 240 V, 50/60 Hz

Operating Temperature: 10°C To 35°C

Operating Humidity: 30% To 90%

Storage Temperature: -20°C To 50°C

Storage Humidity: 10% To 90%

Dimensions:

Approx. 20 cm x 12 cm x 3 cm (L x W x H)

Approx. 8" x 5" x 1.25" (L x W x H)

Weight: < 1 kg < 2 lbs

Standard Accessories: Transducers: DD65v2, DD45 Or IP30, Bone Conductor, Operator Mic / Monitor Headset, Patient Response Switch, Talkback Microphone, External Power Supply, Auditec Sound File License, QuickSIN License, USB Cable, Software & Manuals, Carrying Case

Optional Accessories: DD450 Headphones

MedRx Minimum Computer Specs:

Windows® PC compatible computer, Intel™ i5, 2.0 GHz or better. 4 GB RAM. 20 GB free hard drive space. Available 2.0 USB Port. Windows 10 or 11 Professional (32 Or 64-Bit).



Good Things Come In Small Packages

MedRx, Inc

1200 Starkey Rd. #105

Largo, FL 33771 U.S.A.

888-392-1234 • 727-584-9600

Fax: 727-584-9602

Email: sales@medrx-diagnostics.com

www.medrx-diagnostics.com



Stealth Audiometer the Full Clinical Solution.

