

Oticon CROS

A revolutionary sound experience for single-sided deafness



oticon
life-changing technology



Oticon CROS/BiCROS

A new and innovative solution for single-sided deafness

We understand the struggle of living with single-sided deafness (SSD). When you can hear from only one ear, sounds coming from other directions are difficult to understand. Missing sounds and conversations on your “bad side” is a constant frustration. Relying on just one ear to hear, especially in difficult listening environments, is tiring and challenging.




Oticon CROS is the solution you've been hoping for, providing valuable and life-changing benefits:


- Hear more sounds from all directions - even on your poorer hearing side
- Stay focused and understand conversations, even in noisy environments
- Ability to be aware of sound on poorer ear side when watching TV or streaming music and audio
- No surgical procedures necessary

You don't have to compromise on your hearing anymore. With Oticon CROS, you can actively participate in conversations all around you.

The Oticon CROS hearing device consists of two parts:



A transmitter with a microphone that wirelessly picks up and sends the sound from your poorer ear to the better ear.



A hearing aid that picks up the sound from the transmitter and plays it in your better ear.

Connect to the world around you

Innovative TwinLink™ dual-streaming technology in Oticon CROS makes it possible to connect to external audio streams like music or TV, while enjoying sound transmission from your poorer ear to your better ear at the same time.



Oticon CROS/BiCROS

A perfect companion with our latest technology*

Oticon Xceed / Xceed Play



BTE SP



BTE UP

Oticon Opn S™ / Opn Play™



miniRITE



miniRITE T



miniRITE R



BTE PP

Oticon Ruby



miniRITE



miniRITE T



miniRITE R



BTE PP



BTE

Available in seven colors:



*Please note: Oticon CROS is not compatible with Oticon Opn S 3 and Oticon Opn Play 2

How will Oticon CROS **change your life?**

Oticon CROS gives you access to speech and sounds from all around you, opening up your world to a more natural, comfortable listening experience. But how does this help you in real life? Imagine more hearing ease in everyday situations like...



Getting to know your Oticon CROS

Ready to get started with your Oticon CROS? Adjusting to any new hearing device can be challenging. Keep these tricks in mind and you'll be on the road to better hearing health in no time!


- Work with your hearing care professional to **fine tune** your Oticon CROS.
- **Try new sounds and environments** with your Oticon CROS - join in the conversation!
- **Give your brain a break** - and give yourself time to adjust to your new devices.
- Note the **difference in sounds and experiences** you have with your Oticon CROS - and make sure you share with your hearing care professional for any necessary adjustments at your next visit.

See what Oticon CROS can do for you

We know you're not willing to compromise on your hearing health, and with Oticon CROS, you don't have to.

To find out more about what Oticon CROS can do for you, talk to your hearing care professional today.

Contact your local hearing care professional:



20-18300 15555-9943/02.20



www.oticon.com

oticon
life-changing technology