

# The Oticon CROS family

A revolutionary sound experience for single-sided deafness



**oticon**  
life-changing technology



## **Oticon CROS/BiCROS**

Innovative solutions for  
single-sided deafness

We understand the struggle of living with single-sided deafness (SSD). When you can hear from only one ear, sounds coming from other directions are difficult to understand. Missing sounds and conversations on your “bad side” is a constant frustration. Relying on just one ear to hear, especially in difficult listening environments, is tiring and challenging.




Oticon CROS PX miniRITE R and Oticon CROS miniRITE T are the solutions you've been hoping for, providing valuable and life-changing benefits:

- Hear more sounds from all directions – even from your poorer hearing side
- Stay focused and understand conversations, even in noisy environments
- Ability to be aware of sound from poorer ear side when watching TV or streaming music and audio
- No surgical procedures necessary
- Rechargeable option with Oticon CROS PX


You don't have to compromise on your hearing anymore. With Oticon CROS/BiCROS, you can actively participate in conversations all around you.

# The Oticon CROS PX/CROS

hearing devices consist of two parts:



A transmitter with a microphone that wirelessly picks up and sends the sound from your poorer ear to the better ear.



A hearing aid that picks up the sound from the transmitter and plays it in your better ear.

## Connect to the world around you

Innovative TwinLink™ dual-streaming technology in Oticon CROS PX/CROS makes it possible to connect to external audio streams like music or TV, while enjoying sound transmission from your poorer ear to your better ear at the same time.



# Oticon CROS PX/CROS

Perfect companions with our latest hearing aid technology

**Oticon Real™, Oticon More™, Oticon Zircon, and Oticon Play PX**



miniRITE T



miniRITE R



miniBTE T



miniBTE R

**Oticon Xceed and Oticon Xceed Play**



BTE SP



BTE UP

**Oticon Opn S™ and Oticon Opn Play™**



miniRITE



miniRITE T



miniRITE R



BTE PP

**Oticon Ruby**



miniRITE



miniRITE T



miniRITE R



BTE PP



BTE

**Our CROS family is available in a variety of colors:**



*Oticon CROS PX only*

\*Please note: Oticon CROS PX/CROS are compatible with all hearing devices listed except Oticon Opn S 3 and Oticon Opn Play 2

# How will Oticon CROS PX/CROS **change your life?**

The Oticon CROS family gives you access to speech and sounds from all around you, opening up your world to a more natural, comfortable listening experience. But how does this help you in real life? Imagine more hearing ease in everyday situations like...



# Getting to know your Oticon CROS/BiCROS

Ready to get started? Adjusting to any new hearing device can be challenging. Keep these tricks in mind and you'll be on the road to better hearing health in no time!

- Work with your hearing care professional to **fine tune** your Oticon device.
- **Try new sounds and environments** with your Oticon device – join in the conversation!
- **Give your brain a break** – and give yourself time to adjust to your new devices.
- Note the **difference in sounds and experiences** you have with your Oticon device – and make sure you share with your hearing care professional for any necessary adjustments at your next visit.

## See what the Oticon CROS family can do for you

We know you're not willing to compromise on your hearing health, and with Oticon CROS PX/CROS, you don't have to.

To find out more about what the Oticon CROS family can do for you, talk to your hearing care professional today.

Contact your local hearing care professional:



20-18300 15555-9943/03.23



[oticon.com](https://www.oticon.com)

**oticon**  
life-changing technology