

OTICON | CROS

## The Oticon CROS family

Introducing Oticon CROS PX miniRITE R - a rechargeable solution for single-sided deafness

life-changing  
technology


**The Oticon CROS family is expanding with the addition of the new rechargeable Oticon CROS PX.** The new Oticon CROS PX is powered by the Polaris™ platform and offers patients the benefits of a rechargeable device. With advanced TwinLink™ dual-streaming technology, the Oticon CROS family offers patients access to all relevant sounds in challenging listening environments while also providing the ability to stream sound simultaneously.



### The Oticon CROS Family

- Offers patients the choice of a rechargeable (Oticon CROS PX) or non-rechargeable (Oticon CROS) transmitter to match their personal preference and hearing aid style.
- Both Oticon CROS and Oticon CROS PX offer forwards and backwards compatibility with Oticon More™, Opn S™\*, Opn Play™\*, Ruby, Xceed and Xceed Play.
- Includes TwinLink an innovative technology featuring simultaneous Near-Field Magnetic Induction (NFMI) and 2.4 GHz Bluetooth® Low Energy streaming, making it possible to carry on a conversation while streaming sound.
- Can be used as either a CROS or BICROS fitting.
- Compatible with miniRITE, miniRITE T, miniRITE R, BTE, BTE PP, BTE SP and BTE UP styles.
- Both styles available in the 7 standard colors with Hear Pink also available for Oticon CROS PX.

\*Except Opn S 3 and Opn Play 2





## Oticon CROS family compatibility



Oticon CROS  
miniRITE T



Oticon CROS PX  
miniRITE R

**Oticon More 1\***: miniRITE R, miniRITE T

**Oticon More 2\***: miniRITE R, miniRITE T

**Oticon More 3\***: miniRITE R, miniRITE T

**Oticon Opn S 1\*\***: miniRITE, miniRITE T, miniRITE R, BTE PP

**Oticon Opn S 2\*\***: miniRITE, miniRITE T, miniRITE R, BTE PP

**Oticon Opn S 3\*\***: miniRITE, miniRITE T, miniRITE R, BTE PP

**Oticon Xceed 1\*\***: BTE SP, BTE UP

**Oticon Xceed 2\*\***: BTE SP, BTE UP

**Oticon Xceed 3\*\***: BTE SP, BTE UP

**Oticon Ruby 1\*\*\***: miniRITE, miniRITE T, miniRITE R, BTE, BTE PP

**Oticon Ruby 2\*\*\***: miniRITE, miniRITE T, miniRITE R, BTE, BTE PP

**Oticon Opn Play 1\*\***: miniRITE, miniRITE T, miniRITE R, BTE PP

**Oticon Opn Play 2\*\***: miniRITE, miniRITE T, miniRITE R, BTE PP

**Oticon Xceed Play 1\*\***: BTE SP, BTE UP

**Oticon Xceed Play 2\*\***: BTE SP, BTE UP

• Indicates feature included

\* Firmware 1.1.0 and onwards

\*\* Firmware 8.0 and onwards

\*\*\* Firmware 1.0 and onwards