

OTICON | CROS

The Oticon CROS family

Help your patients stay in the conversation

life-changing
technology


Oticon CROS and CROS PX transmitters are the perfect transmitters to give your patients with single-sided deafness confidence and peace of mind throughout their day. With advanced TwinLink™ dual-streaming technology, the Oticon CROS family offers patients access to all relevant sounds in challenging listening environments while also providing the ability to stream sound simultaneously.



The Oticon CROS Family

- Offers patients the choice of a rechargeable (Oticon CROS PX) or non-rechargeable (Oticon CROS) transmitter to match their personal preference and hearing aid style.
- Both Oticon CROS and Oticon CROS PX offer forwards and backwards compatibility with Oticon More™, Oticon Zircon, Oticon Play PX, Opn S™*, Ruby, Opn Play™*, Xceed and Xceed Play.
- Includes TwinLink, an innovative technology featuring simultaneous Near-Field Magnetic Induction (NFM) and 2.4 GHz Bluetooth® Low Energy streaming, making it possible to carry on a conversation while streaming sound.
- Can be used as either a CROS or BICROS fitting.
- Compatible with miniRITE, miniRITE T, miniRITE R, miniBTE T, miniBTE R, BTE, BTE PP, BTE SP and BTE UP.
- Both styles available in the 7 standard colors with Hear Pink also available for Oticon CROS PX.

**Except Opn S 3 and Opn Play 2*





Oticon CROS family compatibility



Oticon CROS
miniRITE T



Oticon CROS PX
miniRITE R

Oticon More 1*: miniRITE R, miniRITE T, miniBTE R, miniBTE T

•

Oticon More 2*: miniRITE R, miniRITE T, miniBTE R, miniBTE T

•

Oticon More 3*: miniRITE R, miniRITE T, miniBTE R, miniBTE T

•

Oticon Opn S 1**: miniRITE, miniRITE T, miniRITE R, BTE PP

•

Oticon Opn S 2**: miniRITE, miniRITE T, miniRITE R, BTE PP

•

Oticon Opn S 3**: miniRITE, miniRITE T, miniRITE R, BTE PP

Oticon Xceed 1**: BTE SP, BTE UP

•

Oticon Xceed 2**: BTE SP, BTE UP

•

Oticon Xceed 3**: BTE SP, BTE UP

•

NEW **Oticon Zircon 1***: miniRITE R, miniRITE T, miniBTE R, miniBTE T

•

NEW **Oticon Zircon 2***: miniRITE R, miniRITE T, miniBTE R, miniBTE T

•

Oticon Ruby 1***: miniRITE, miniRITE T, miniRITE R, BTE, BTE PP

•

Oticon Ruby 2***: miniRITE, miniRITE T, miniRITE R, BTE, BTE PP

•

NEW **Oticon Play PX 1***: miniRITE T, miniRITE R, miniBTE R, miniBTE T

•

NEW **Oticon Play PX 2***: miniRITE T, miniRITE R, miniBTE R, miniBTE T

•

Oticon Opn Play 1**: miniRITE, miniRITE T, miniRITE R, BTE PP

•

Oticon Opn Play 2**: miniRITE, miniRITE T, miniRITE R, BTE PP

•

Oticon Xceed Play 1**: BTE SP, BTE UP

•

Oticon Xceed Play 2**: BTE SP, BTE UP

•

• Indicates feature included

* Firmware 1.1.0 and onwards

** Firmware 8.0 and onwards

*** Firmware 1.0 and onwards