

life-changing
technology

Oticon Government Services

Selection Guide Book



Phone compatibility check



oticon
life-changing technology

RITE Styles			
Oticon Intent 1 miniRITE R (lithium-ion)	Oticon Real 1 miniRITE R (lithium-ion)	Oticon Play PX 1 miniRITE R (lithium-ion) Oticon Opn S 1 miniRITE R (lithium-ion) Meets SCIF requirements!	Oticon Real 1 miniRITE T (312 zinc-air)

60 85 100 105 60 85 100 105 60 85 100 105 60 85 100 105



BTE Styles				
Oticon Real 1 miniBTE R (lithium-ion)	Oticon Real 1 miniBTE T (312 zinc-air)	Opn S 1 BTE Plus Power 13	Xceed 1 BTE Super Power 13	Xceed 1 BTE Ultra Power 675

85 85 105 110 120



Custom Styles				
Oticon Own 1 IIC	Oticon Own 1 CIC	Oticon Own 1 ITC	Oticon Own 1 ITE HS	Oticon Own 1 ITE FS

75 90 75 90 75 90 100 75 90 100 75 90 100



Oticon CROS Family Compatibility			
Oticon CROS (312 zinc-air) Oticon CROS PX (rechargeable)			
Oticon Intent 1	Oticon Real 1	Oticon More 1 Oticon Play PX 1	Opn S 1 Xceed 1




Oticon CROS



Oticon CROS PX


Intent Fitting Ranges



60 85

□ OpenBass dome
■ Mold, Bass & Power dome


RITE Fitting Ranges




60 85

□ Open Bass/Open dome
■ Mold, Bass & Power dome

BTE Fitting Ranges




85 105




110 120

□ Corda miniFit standard thin tube □ Corda miniFit Power thin tube ■ Hook



100 105

■ MicroShell Detect 100, Bass & Power dome ■ MicroShell Detect 105





100 105

■ Bass, Power dome & Power Receiver Mold ■ Power Receiver Mold

Custom Fitting Ranges

IIC CIC ITC ITE HS ITE FS



Oticon RITE and BTE Colors

C044 Silver	C092 Steel Grey	C091 Silver Grey	C090 Chroma Beige	C116** Honey Beige	C094 Terracotta	C093 Chestnut Brown	C063 Diamond Black
		C079 Hear Pink*	C115** Sky Blue	C114† Olive Green			

Oticon Play PX Colors Only

C048 Emerald Green	C058 Aquamarine	C047 Cool Blue	C045 Purple	C046 Red	C057 Power Pink	C079 Baby Pink
------------------------------	---------------------------	--------------------------	-----------------------	--------------------	---------------------------	--------------------------

*Only available for Oticon Real and CROS PX.
**Only available for Oticon Intent.
†Only available for Oticon Real miniRITE R.

Oticon Custom Colors

C001 Beige	C002 Light Brown	C003 Medium Brown	C004 Dark Brown
C005 Black	Transparent*	Red*	Blue*

*Not available as a faceplate color.

¹Oticon Intent SmartCharger Only.
²Oticon Intent Charger only compatible with Oticon Intent.
³Oticon miniRITE SmartCharger is compatible with Oticon More devices on Firmware 1.1.1 or newer. Oticon miniRITE SmartCharger is not compatible with Opn S miniRITE R.
⁴Features for both Oticon SmartCharger and Desktop charger.
⁵Not compatible with Opn S.

Chargers

- Oticon Intent™ Chargers**
- Quick charge in 15 minutes
 - Advanced contact charging
 - Fitting during charging available
 - Power Bank¹



Oticon Intent SmartCharger²

NEW



Oticon Intent Charger²

Oticon SmartChargers for Oticon Real™, Oticon More™ and Oticon Play PX³

- Protective Lid
- Power Bank
- Drying Capabilities



Oticon miniRITE SmartCharger



Oticon miniBTE SmartCharger

Oticon Desktop Chargers for Oticon Real, More, Play PX and Opn S™

- Easy insertion and removal
- Contact Free⁴
- Quick charge⁴




Oticon miniRITE Charger



Oticon miniBTE Charger⁵

miniFit Earpiece Compatibility Overview




*Opn S only.

**Not compatible with Opn S except for 5mm OpenBass dome.

RITE molds

BTE Molds




*Compatible with miniFit receivers 60 and 85.

**Oticon Intent only. Compatible with Oticon Intent miniFit Detect receivers 60, 85, 100 and 105 Receivers.

***Not available for Oticon Intent.

Standard Earmold Styles*



*Reference pages 38-39 of pocket guide for more information on these earmold styles.

Special Earmold Styles*



*Reference pages 38-39 of pocket guide for more information on these earmold styles.

Wireless Accessories*



EduMic

- Automatic pairing
- 3.5mm input jack
- Telecoil
- Easy-to-use remote microphone



TV Adapter 3.0

- Seamlessly stream audio from a TV
- Wireless streaming connection up to 50 ft
- Digital/RCA connection



ConnectClip

- Make and receive hands-free calls to both Android™ and iPhone® devices
- Remote Microphone
- Control the volume and change programs



Remote Control 3.0

- Adjust volume
- Switch programs
- Mute



Phone Adapter 2.0

- Turn a landline phone into a wireless Bluetooth phone
- Wireless streaming connection up to 65 ft
- ConnectClip required



BTDB800 Bluetooth Dongle

- Small USB device to improve wireless signal on computers
- Stream video calls, movies and music from your computer
- ConnectClip required

*Only available in styles with 2.4 GHz.