

miniFit universe



| Domes and earpieces | | Polaris™ miniFit receivers and Corda miniFit compatibility: | | | | | | Velox S™ & Velox™ miniFit receivers and Corda miniFit compatibility: | | | | | |
|-------------------------------|--|----------------------------------------------------------------|----|-----|-----|-----------|-----------|-------------------------------------------------------------------------|----|-----|-----|-----------|-----------|
| | | 60 | 85 | 100 | 105 | Corda 0.9 | Corda 1.3 | 60 | 85 | 100 | 105 | Corda 0.9 | Corda 1.3 |
| OpenBass dome 5 mm | | ● | | | | ● | ● | | | | | | |
| OpenBass dome 6-12 mm | | ● | ● | | | ● | ● | | | | | | |
| Open dome 5 mm | | | | | | | | ● | | | | ● | |
| Open dome 6-10 mm | | | | | | | | ● | ● | | | ● | |
| Bass dome Double vent 6-12 mm | | ● | ● | ● | | ● | ● | ● | ● | ● | | ● | ● |
| Bass dome Single vent 6-12 mm | | | | | | | | ● | ● | ● | | ● | ● |
| Power dome 6-12 mm | | ● | ● | ● | | ● | ● | ● | ● | ● | | ● | ● |
| Grip Tip | | ● | ● | | | ● | ● | ● | ● | | | ● | ● |
| VarioTherm® LiteTip | | ● | ● | | | | | ● | ● | | | | |
| VarioTherm® MicroMold | | ● | ● | | | ● | ● | ● | ● | | | ● | ● |
| LiteTip | | ● | ● | | | | | ● | ● | | | | |
| MicroMold | | ● | ● | | | ● | ● | ● | ● | | | ● | ● |
| MicroShell | | ● | ● | | | | | | | | | | |
| Power receiver mold | | | | ● | ● | | | | | ● | ● | | |

● = Available Pictures not to scale

22-388363 15500-1124/09.22



oticon.com