

OTICON | **More**

More than



meets the ear.

More sound. More freedom. More life.



oticon
life-changing technology

It's time to get more out of life



DID YOU KNOW?

You hear with your brain, not your ears. Your ears capture sounds so your brain can identify them, locate them, and decide which ones to listen to.

At Oticon, we focus on understanding the brain and develop technology that supports it. New research tells us that the brain can handle access to more information from sounds.* Having a hearing aid that supports your brain, by delivering sounds in a precise and clear way, is all the more important.

* O'Sullivan, et al (2019). Puvvada, K. C., & Simon, J. Z. (2017).



Your brain needs the full perspective

The sounds around you are constantly changing in unpredictable ways - everywhere you go. And those sounds contain far more richness than just speech. Your brain recognizes all of this complexity and creates meaning from it.



OTICON | More Hear the sounds that shouldn't be missed



Designed to support your brain

The human brain is the gold standard for processing sound. So we developed Oticon More™ to work more like how the brain works because it learned through experience. This means Oticon More supports your brain by giving it more of what it needs to make better sense of sound.

Powered by a Deep Neural Network


Oticon More is the world's first hearing aid with a Deep Neural Network (DNN) embedded on the chip. The DNN is a technology we have trained with 12 million real-life sounds. This built-in intelligence ensures each individual sound is delivered with detail and clarity.



More than meets the ear

Oticon More gives your brain the clearest sound information, allowing you to hear sounds that shouldn't be missed. You benefit with improved speech understanding with less listening effort and the ability to remember more.





**Proven to give
you superior sound**

8
out of
10
people prefer
Oticon More

Eight out of 10 people preferred the sound quality of Oticon More over two leading competitors in situations including a conversation in a busy restaurant, café and speech through a facemask.*

*Man et al 2021



Oticon More delivers

30%

more sound to
the brain**

15%

increase in speech
understanding**

compared to industry leading Oticon Opn S™



**Santurette et al 2020

More options
More styles
More freedom



Why settle when you can get what you want?

Oticon More now gives you more choices along with superior sound quality. Choose the style that suits you—a rechargeable miniRITE R or miniRITE T with disposable batteries—in one of eight colors.



C044
Silver



C091
Silver Grey



C092
Steel Grey



C063
Diamond Black



C093
Chestnut Brown



C094
Terracotta



C090
Chroma Beige



C079
Hear Pink



A full day's power. Every day.

The rechargeable Oticon More miniRITE R comes with a stylish desktop charger that gives you a full day of battery life in just three hours* – even if you stream sound from other devices.



- No disposable batteries = less hassle & waste
- 30-minute recharge provides an additional 6 hours of power
- Convenient, easy-to-use, reliable

More power on the go

You now have the option to purchase our new easy-to-use SmartCharger to give you freedom to enjoy life on the go. Keep your Oticon More miniRITE R devices powered up, dry and projected when recharging, wherever life takes you.




*Lithium-ion battery behavior varies depending on hearing loss, lifestyle, and streaming behavior.

Connect to your favorite devices

With Oticon More you can enjoy high-quality sound streaming from your favorite devices.



Direct streaming from iPhone® and Android™ devices*



The Oticon ON app - Control at your fingertips

Oticon ON app lets you personalize your listening experience when connected to your device. You can adjust the volume, change the program, check the battery level, and control other Oticon devices. It even comes with a "Find my hearing aids" feature.



*Android devices need to support Audio Streaming for Hearing Aids (ASHA) to allow direct streaming to Oticon More. Please visit oticon.com/support/compatibility for more information.

Apple, the Apple logo, iPhone®, iPad®, and iPod touch® are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android™, Google Play, and the Google Play logo are trademarks of Google LLC.

A wide range of connectivity options



ConnectClip

Use ConnectClip as a discreet microphone accessory that wirelessly connects to your Bluetooth devices.



TV Adapter

Use TV Adapter to stream sound from the TV directly to Oticon More hearing aids.



Remote Control

Adjust volume, switch programs, or mute the hearing aids with the touch of a button.

Oticon More is music to your ears

Now, you can add richness to your overall listening experience with Oticon MyMusic. Whether you're listening to live music or streaming from a device, you can have a dedicated setting to ensure the best possible sound reproduction.



**To find out what Oticon More can do for you,
talk to your hearing care professional today.**

life-changing
technology

Contact your local hearing care professional:

21-205850 15500-0645/08.21



oticon.com

oticon
life-changing **technology**