

Hands-free communication is here

Now your patients can have seamless hands-free calls using their Oticon hearing aids as a headset

Your patients with an iPhone or iPad can now use their Oticon hearing aids* to make or receive hands-free phone or video calls.



Hands-free
communication
Now available



Convenient. High-quality. Hands-free

Your patients can now enjoy the simple convenience of hands-free calls whenever it suits them. With their hands free, they can talk to their families, friends and colleagues while they're on the go. Whether they're driving, walking, working or cooking, they can keep doing what they're doing and let their hearing aids take care of the rest.

- Video calls and phone calls with high-quality, direct audio
- Better, more present remote meetings
- Free hands to do other things during calls

Made for
iPhone | iPad

* Oticon More, Oticon Zircon and Oticon Play PX hearing aids on the Polaris platform



Hands-free communication works for:

- Mobile phone calls, for instance, FaceTime
- Video calls
- Online meetings

How it works:

1. The wearer makes or receives a call on their iPhone or iPad
2. Call audio is sent wirelessly to the wearer's hearing aids
3. The hearing aids microphones capture the wearer's voice
4. The wearer's voice is sent to the iPhone or iPad



Built on the latest Bluetooth Low Energy technology

Oticon's hands-free communication is at the forefront of technology, so it uses less of the hearing aids' batteries than the classic Bluetooth found in some other hearing aid brands.



Ready for everyone on the Polaris platform

Hands-free communication is now available for these hearing aids on the Polaris platform: Oticon More™, Oticon Zircon, and Oticon Play PX. For patients who already have Oticon More, you simply need to update their software to 1.3.0 to enable the new hands-free communication functionality.



Hands-free communication is available with iPhone 11 or later running iOS 15.2 or later, and iPad running iPadOS 15.2 or later.**

** See which hearing aids and devices are compatible here:
oticon.com/support/compatibility

22-266444 15500-9022/01.22



oticon.com

Oticon is part of the Demant Group.

oticon
life-changing technology