


OTICON | **Play PX**

Unleash the magic of childhood



oticon
life-changing technology



Introducing the **new Oticon Play PX**

Experiences are important to children's development

Childhood is a time of adventures and of learning and children develop from these experiences. Without assistance, children with hearing loss miss out on important details that are essential to development.

Supports development with access to the full sound environment

New Oticon Play PX is designed to offer your child full access to the hearing input needed to experience every magical moment of childhood. Advanced technology provides communication access to vital input from all meaningful sounds, which is important to development.



Powered by a Deep Neural Network

To best support the brain, Oticon Play PX is powered by Polaris™, our most intelligent chip ever. It is also the world's first pediatric hearing aid with an on-board Deep Neural Network that has learned to recognize sounds through experience, like children's brains do naturally. Trained with 12 million real-life sound scenes, the Deep Neural Network knows how virtually all sounds should sound.

Give your child a winning start and **unleash the magic of childhood.**



Backed by science. **Proven by children**



For children to develop, they must be able to understand and take part in the world around them. We tested the effect Oticon Play PX has on children's ability to hear and follow what is being said and happening in typical busy everyday environments like classrooms and playgrounds.

The results showed that Oticon Play PX helps children to recognize and recall words, even when speech comes from multiple directions.* This allows your child to participate freely in the many different activities and situations that children find themselves in – which is crucial to development, growth, and learning.

* Gordey, D. & Ng, E. 2021, Oticon Play PX: Supporting Communication, Learning and Inclusion for Children and Teens, Oticon Whitepaper



Better hearing, better learning

Did you know that we hear with our brains and not our ears?

That is why it is important to give the brain the best quality sound. In fact, the brain needs access to the full sound scene to work properly.**

But if your child has a hearing loss, their brain needs to work harder to figure out what others are saying. They also miss out on important sound details that help their brain learn and develop.

Oticon Play PX helps provide the brain with the best conditions to hear so it can make sense of sound more easily.

At Oticon, we call this unique approach BrainHearing™. This gives your child's brain an advantage that is especially important during the critical period of development.

** Gordey, D. & Ng, E. 2021, Pediatric BrainHearing, Oticon Whitepaper

Let children **connect** **with the world**



Oticon Play PX lets your child connect seamlessly to the world. They can use their hearing aids to connect to classroom electronics, smartphones, TVs, computers, and more.



EduMic

Oticon's EduMic is a small wireless microphone that teachers wear to ensure your child hears what is going on in educational settings like classrooms. It promotes participation and provides the best conditions to learn and develop.

Connect to their favorite devices

With Oticon Play PX your child can enjoy high-quality sound streaming from their favorite devices.



Stream directly from iPhone® and select Android™ devices

Enjoy high-quality streaming directly from iPhone, iPad® and select Android devices.*



Stream sound from your TV

With Oticon's TV Adapter your child can stream sound directly from the TV.



Connect a remote microphone

Wirelessly stream someone's voice to your child's hearing aids with EduMic or ConnectClip.

* Android devices need to support Audio Streaming for Hearing Aids (ASHA) to allow direct streaming to Oticon Play PX. Please visit oticon.com/support/compatibility for more information.

Gain control with the Oticon ON app

With the Oticon ON app you and your child can monitor and control Oticon Play PX from your smartphones - adjust the volume, change the program or check the battery status. The app also includes a "Find my hearing aid" feature if they are misplaced.



Made for

 iPhone | iPad | iPod

Works with

android 



Download on the
App Store



GET IT ON
Google Play

Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

A background image showing a child's hands holding a white fabric bag, possibly containing hearing aids, against a blurred outdoor setting.

A child-friendly and **rechargeable** **hearing aid**

Oticon Play PX offers two recharging options, the powerful desktop charger and the portable SmartCharger.

Powerful desktop charger

Oticon Play PX comes with a stylish desktop charger, ideal for charging at home.

Portable SmartCharger

With the option to purchase the portable SmartCharger, your child can always keep their hearing aids powered, dry, and protected.

It gives children the freedom and confidence to participate in camps, sleepovers, and excursions knowing the hearing aids can easily be recharged.

miniBTE R charging times

- Fully charges your child's hearing aids in three and a half hours*
- Quick charge of 30 minutes gives an additional five hours of power**



Tamper resistance

The miniBTE R style is tamper resistant because the battery is completely sealed in the hearing aid, and the miniBTE T style is available with a tamper resistant battery drawer.



* Or three hours for the miniRITE R style.

** Or six hours for the miniRITE R style.

Made for adventure



The robust Oticon Play PX has been designed to stand up to active young lives and with safety in mind.

IP68-certified

Oticon's IP68-certified design includes water, dust, and debris barriers and is nanocoated to repel water. It also protects casings and electrical enclosures from water.

LED light

The LED indicator on Oticon Play PX gives you visual confirmation that the hearing aid is functioning properly.

Free of allergens

Oticon Play PX meets international standards on safety and biocompatibility and contains no phthalates. The hearing aids and accompanying stickers are biologically safe for permanent skin contact.

Choose the hearing aid that suits your child

Oticon Play PX covers many different types of hearing loss and is available in four different styles. Select from 12 attractive colors and a choice of rechargeable or non-rechargeable batteries.



miniBTE R



miniBTE T




miniRITE R



miniRITE T




Terracotta



Chestnut Brown



Chroma Beige



Diamond Black



Silver



Baby Pink



Emerald Green



Aquamarine



Cool Blue



Purple



Cool Red



Power Pink

To find out what Oticon Play PX can do for you, talk to your hearing care professional today or go to [oticon.com/play-px](https://www.oticon.com/play-px)

life-changing
technology

Contact your local hearing care professional:

21-254012 15555-10535/08.22



[oticon.com](https://www.oticon.com)

oticon
life-changing **technology**