


Hands-free communication is here

Now your pediatric patients can have seamless hands-free calls using their Oticon hearing aids as a headset

Your pediatric patients with an iPhone® or iPad® can now use their Oticon hearing aids* to make or receive hands-free phone or video calls. They can also use their hearing aid buttons to answer, reject, and end calls.



**Hands-free
communication
Now available**



Convenient. High-quality. Hands-free

Your pediatric patients can now enjoy the simple convenience of hands-free calls whenever it suits them. With their hands free, they can talk to their families, friends and classmates while they're on the go. Whether they're walking, bicycling or running, they can keep doing what they're doing and let their hearing aids take care of the rest.

- Video calls and phone calls with high-quality, direct audio
- Easier to attend for remote learning
- Free hands to do other things during calls

Made for
iPhone | iPad

* Oticon More™, Oticon Own™ with 2.4 GHz, Oticon Zircon and Oticon Play PX hearing aids on the Polaris™ platform.

Hands-free communication works for:

- Mobile phone calls, for instance, FaceTime
- Video calls
- Online meetings

How it works:

1. The wearer makes or receives a call on their iPhone or iPad
2. Call audio is sent wirelessly to the wearer's hearing aids
3. The hearing aid's microphones capture the wearer's voice
4. The wearer's voice is sent to the iPhone or iPad



Built on the latest Bluetooth® Low Energy technology

Oticon's hands-free communication is at the forefront of technology, so it uses less of the hearing aids' batteries than the classic Bluetooth found in some other hearing aid brands.

Ready for everyone on the Polaris platform

Hands-free communication is available for a range of hearing aid styles on the Polaris platform - including Oticon More, Oticon Own*, Oticon Zircon, and Oticon Play PX. Oticon More requires 1.4.0 firmware update for push button control; Oticon Play PX and Oticon Zircon require 1.2.0 firmware update for push button control.

Hands-free communication is available with iPhone 11 or later running iOS 15.2 or later, and iPad running iPadOS® 15.2 or later.**

* Only available in styles with 2.4 GHz connectivity

** See which hearing aids and devices are compatible here:
oticon.com/support/compatibility

