

# Oticon Pediatrics Selection Guide



Feature	Play PX miniBTE R	Play PX miniBTE T	Play PX miniRITE R	Play Px miniRITE T	Opn Play BTE PP	Xceed Play SP*	Xceed Play UP*
MoreSound Intelligence™	•	•	•	•			
OpenSound Navigator™					•	•	•
Tamper resistant battery door	•		•		•	•	•
Telecoil	•	•	•	•	•	•	•
Double push-button			•	•	•	•	•
Single push-button	•	•				•	•
Corda compatible: 0.9mm	•	•			•		
Corda compatible: 1.3mm					•		
CROS/CROS PX	•	•	•	•	•††	•	•
Rechargeable	•		•				
Charging times	3.5 hours for full day** 30 minutes for 5 hours		3 hours for full day** 30 minutes for 6 hours				
Zinc air battery size		312		312	13	13	675
On-ear DAI/FM Compatibility with FM10					•	•	•
EduMic (compatible with Roger X receiver)	•	•	•	•	•	•	•
Made for iPhone®	•	•	•	•	•	•	•
Direct streaming to Android™***	•	•	•	•			
Hands-free communication†	•	•	•	•			
Compatible with ConnectClip	•	•	•	•	•	•	•

\*Xceed Play has single and double push-buttons. Single push-button changes program, double push-button changes volume and mutes.

\*\*Lithium-ion battery performance varies depending on hearing loss, lifestyle, and streaming behavior.

\*\*\*Android devices need to support Audio Streaming for Hearing Aids (ASHA) to allow direct streaming to Oticon Play PX.


†Hands-free communication is available with iPhone 11 or later running iOS 15.2 or later, and iPad® running iPadOS 15.2 or later. Please visit [oticon.com/support/compatibility](https://oticon.com/support/compatibility) for more information.

††CROS and CROS PX are not compatible with Opn Play 2 BTE PP

# DSL fitting ranges

Oticon Play PX miniBTE


85



Hook  
Corda miniFit (thin tube)

Opn Play BTE PP


105



Hook  
Corda miniFit

Oticon Xceed Play SP


110



DSL Fitting Range  
Hook, undamped

Oticon Xceed Play UP


120



DSL Fitting Range  
Hook, undamped


Oticon Play PX miniRITE

60




Mold, Bass & Power dome  
OpenBass dome

85




Mold, Bass & Power dome  
OpenBass dome

100



Power receiver mold,  
Bass & Power dome

105




Power receiver mold

## 12 different color options




C094  
Terracotta



C093  
Chestnut Brown



C090  
Chroma Beige



C063  
Diamond Black




C044  
Silver



C079  
Baby Pink




C048  
Emerald Green




C058  
Aquamarine




C047  
Cool Blue



C045  
Purple



C046  
Cool Red



C057  
Power Pink