

miniFit universe



Domes and earpieces		Polaris R™ & Polaris™ miniFit receivers and Corda miniFit compatibility:						Velox S™ & Velox™ miniFit receivers and Corda miniFit compatibility:					
		60	85	100	105	Corda 0.9	Corda 1.3	60	85	100	105	Corda 0.9	Corda 1.3
OpenBass dome 5 mm		●				●	●	●				●	●
OpenBass dome 6-12 mm		●	●			●	●						
Open dome 6-10 mm								●	●			●	
Bass dome Double vent 6-12 mm		●	●	●		●	●	●	●	●		●	●
Bass dome Single vent 6-12 mm								●	●	●		●	●
Power dome 6-12 mm		●	●	●		●	●	●	●	●		●	●
Grip Tip		●	●			●	●	●	●			●	●
VarioTherm® LiteTip		●	●			●	●	●	●			●	●
VarioTherm® MicroMold		●	●					●	●				
LiteTip		●	●			●	●	●	●			●	●
MicroMold		●	●					●	●				
MicroShell		●	●										
Power Receiver Mold				●	●					●	●		

● = Available Pictures not to scale

22-388363 15500-1124/01.23

