

OTICON | Opn Play

Product comparison

	Oticon Opn Play 1	Oticon Opn Play 2
Speech Understanding		
OpenSound Navigator™	Level 1	Level 3
- Balancing power effect	100%	50%
- Max. noise removal	9 dB	3 dB
OpenSound Optimizer™	•	•
Speech Guard™ LX	Level 1	Level 3
Spatial Sound™ LX	4 estimators	2 estimators
Speech Rescue™ LX	•	•
Sound Quality		
Clear Dynamics	•	-
Spatial Noise Management	•	-
Fitting Bandwidth*	10 kHz	8 kHz
Processing Channels	64	48
Bass Boost (streaming)	•	•
Listening Comfort		
Transient Noise Management	4 configurations	On/Off
Feedback shield LX	•	•
Wind Noise Management	•	•
Optimizing Fitting		
YouMatic™ LX	3 configurations	1 configurations
Fitting Bands	16	12
REM AutoFit	Verifit®LINK, IMC 2**	Verifit®LINK, IMC 2**
Pediatric Fitting Mode	•	•
DSL Fitting Range ***	•	•
Fitting Formulas	DSL v5.0, NAL-NL1 + 2, VAC+	DSL v5.0, NAL-NL1 + 2, VAC+
Designed for Children		
LED	•	•
Tamper Resistant Battery Drawer	•	•
Hypo Allergenic	•	•
IP Rating	IP68	IP68
Nano Coating	•	•
Color Options	12	12
Integrated 2.4 GHz Receiver	•	•
DAI/FM Compatibility	•	•



IP68

Made for iPhone | iPad | iPod

For information on compatibility, please visit www.oticon.global/compatibility

oticon
PEDIATRICS

* Bandwidth accessible for gain adjustments during fitting

** Inter Module Communication 2

*** Available in the Technical Data sheets and Opn Play Product Guide

Oticon Opn Play™ style options

BTE PP, miniRITE R, miniRITE T, miniRITE



Oticon Opn Play is a Made for iPhone hearing aid and all styles offer 2.4 GHz wireless technology.*

Oticon Opn Play BTE PP - an all-round pediatric instrument with double push button and LED light for monitoring hearing aid status.

Oticon Opn Play miniRITE R - discreet rechargeable hearing aid powered by a lithium-ion battery. Featuring a telecoil, a double push button for easy volume and program control and LED light for monitoring hearing aid status.

Oticon Opn Play miniRITE T - discreet style featuring a telecoil, a double push button for easy volume and program control and LED light for monitoring hearing aid status.

Oticon Opn Play miniRITE - small, discreet, with single push button for volume and program changes.

* Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.


© VarioTherm is a registered trademark of Dreve

miniFit receivers and Corda miniFit (thin tube) options


Select between three receivers with different output performance, with lengths 0-5 for non power speakers and 1-5 for power.



Power receiver molds



Corda miniFit



miniRITE and Corda miniFit earpieces

miniFit earpieces include standard domes, Micro mold, and LiteTip. Custom earpieces are available in acrylic, soft silicone or VarioTherm® material.



DSL fitting ranges

