

Oticon Amigo and Oticon hearing aids

DAI/FM overview

Hearing aids	DAI adapter		FM adapter	FM receiver	Modes of operation
<p>Automatically adds 2 programs to the hearing aid: - DAI/FM + HA (default) - DAI/FM only</p>	<p>Low Level Input Bodyworn FM solutions e.g. Amigo R5</p>	<p>High Level Input TV, stereo, MP3 player, CD player, computer</p>	<p>Adapter for Amigo receivers and other ear-level receivers</p>	<p>Universal and dedicated Amigo receivers</p>	<p>Default modes of operation when using DAI adapters and FM receivers</p>
 <p>Oticon Alta2/Nera2/Ria2 Oticon Alta/Nera/Ria Oticon Dynamo Oticon Sensei Oticon Sensei SP Oticon Safari Oticon Safari SP Oticon Chili Oticon Agil Oticon Acto Oticon Ino Oticon Geno Oticon Get** (BTE & RITE styles except designRITE, miniRITE and miniBTE)</p>	<p>AP 900 399-50-520-08 Grey</p> 	<p>FM 9 399-50-590-00 Transparent</p> 	<p>Amigo R2*</p> 	<p>Amigo R2 switch positions: ○ or ○○ = HA + FM receiver ● = HA only/FM off</p>	
			<p>Amigo R12 17 colours</p> 		<p>Amigo Arc, short 100637 Amigo Arc, long 140-07-410-00 5 colours</p> 
			<p>Streamer Pro white 120178 Streamer Pro black 120179</p> 	<p>Amigo R2</p> 	
			 <p>Oticon Opn S™ BTE PP Oticon Opn Play™ BTE PP Oticon Opn™ BTE13 PP (FW 6.0) Oticon Xceed Oticon Xceed Play Oticon Siya BTE PP & BTE Oticon Alta2/Nera2/Ria2 BTE13 PP Oticon Geno BTE PP</p>	<p>AP 1000 142207 Grey</p> 	<p>FM 10 142328 Transparent</p> 
<p>Amigo R12G2 14 colors</p> 	<p>Amigo Arc 5 colours</p> 	<p>Via telecoil</p>			
<p>Amigo R2</p> 					
<p>Amigo R7</p> 	<p>Amigo R7 switch positions: ○ = FM ON ● = FM OFF</p>				
<p>Amigo Arc 5 colours</p> 		<p>Via telecoil</p>			
<p>Oticon Sumo family</p>	<p>AP 700 399-50-900-08 Grey</p>  <p>Input levels are set manually</p>	<p>FM 7 399-50-920-02 Grey</p> 	<p>Amigo R2</p> 	<p>Amigo R2 switch positions: ○ or ○○ = HA + FM receiver ● = HA only/FM off</p>	
<p>Oticon Sumo family</p>	<p>Amigo R7</p> 	<p>Amigo R7 switch positions: ○ = FM ON ● = FM OFF</p>			
			<p>Amigo Arc 5 colours</p> 		<p>Via telecoil</p>

* On 312 battery style, the LED on Amigo R2 must be set to blink.

** Streamer Pro is not compatible with Oticon Get family.


*** Other configurations are possible using Genie.

See Amigo
compatibility guide
www.oticon.global/amigo

oticon
PEOPLE FIRST

Order information

Amigo Star

Colours	169-187 MHz	216-218 MHz	
● C090 Chroma beige	118016	124404	
● C093 Chestnut brown	128687	128692	
● C047 Cool Blue	128690	128695	
● C046 Cool Red	128689	128694	
○ C044 Silver	128688	128693	
● C063 Diamond black	128686	128691	

Amigo R12

Amigo R12G2

Colours	169-187 MHz	216-218 MHz	169-187 MHz	216-218 MHz
● Champagne	140-07-000-00	140-07-100-00		
● C090 Chroma beige	140-07-010-00	140-07-110-00	142137	162482
● C092 Steel grey	140-07-011-00	140-07-111-00	188314	188336
● C091 Silver grey	140-07-012-00	140-07-112-00	188313	188335
● C093 Chestnut brown	140-07-013-00	140-07-113-00	188310	188332
● C043 Baby pink	140-07-014-00	140-07-114-00	189051	189055
● C042 Baby blue	140-07-015-00	140-07-115-00		
● C047 Cool Blue	140-07-016-00	140-07-116-00	188317	188339
● C046 Cool Red	140-07-017-00	140-07-117-00	188324	188344
● C034 Purple	140-07-018-00	140-07-118-00	188323	188343
○ C044 Silver	140-07-019-00	140-07-119-00	188312	188334
● C094 Terracotta	102404	102408	188311	188333
● C063 Diamond black	101821	101827	188315	188337
● C048 Emerald green	131228	131234	188318	188340
● C058 Aquamarine	162058	162070	188316	188338
● C057 Power pink	162057	162069	188322	188342
○ C098 Pure white	162056	162068		

Amigo R2

Amigo R7

Colours	169-187 MHz	216-218 MHz	169-187 MHz	216-218 MHz
● Champagne	140-06-200-00	140-06-300-00	140-06-600-00	140-06-700-00

DAI Cords

Cord, Monaural	CD or MP3 Players, computer, etc	383-58-940-09
Cord, Binaural	CD or MP3 Players, computer, etc.	383-58-941-00
DAI Y cord	60 cm (24") white	384-01-066-00
DAI Y Cord	60 cm (24") grey	384-01-063-00
DAI Y Cord	75 cm (30") white	384-01-070-00
DAI Y Cord	75 cm (30") grey	384-01-068-00
DAI Cord	60 cm (24") white	384-01-065-00
DAI Cord	60 cm (24") grey	384-01-062-00
DAI Cord	75 cm (30") white	384-01-069-00
DAI Cord	75 cm (30") grey	384-01-067-00

See Amigo
compatibility guide
www.oticon.global/amigo

oticon
PEOPLE FIRST