

Pediatric earmold colors

Silicone (MicroMold and BTE Mold)

Shore 60



Transparent

Pink

Shore 40



Transparent

Rose

White

Black

Yellow


Orange

Red

Green

Blue

Acrylic (MicroShell Detect, MicroShell and Power Receiver Mold)*



Acrylic/OtoTherm™ (LiteTip, MicroMold and BTE Mold)

Only available in transparent



*MicroShell Detect is for Sirius™ platform instruments and later. Power Receiver Mold is available for Polaris R™ platform instruments and prior platforms.