



**People First**

People First is our promise to empower people to communicate freely, interact naturally and participate actively

**Power your senses**

**Let your child get the most out of speech**

**Your child deserves the best**

Oticon Sensei Super Power has been developed with one goal in mind: **to increase your child's understanding of speech.**

Supported by BrainHearing™ technology, Sensei Super Power provides your child with access to the full spectrum of sounds and processes speech in the clearest way possible. With Oticon Sensei Super Power, your child will get the audibility to hear more speech sounds, especially the important high frequency sounds, such as “s” and “th”. Giving your child more speech details than ever before will help your child perceive speech more clearly.

All this care will allow for the best spoken language development. It will also support your child's brain when it comes to making sense of speech sounds in order to improve speech understanding - which is so important for children with severe to profound hearing losses.

**child friendly hearing care**

Our pediatric audiological mission is to ensure a better future for every child with hearing loss. We will deliver solutions, tools and techniques that optimize auditory and cognitive habilitation, embrace the complexities of growing up with hearing loss and empower you to adapt solutions to each child's developmental stage on their journey to adulthood.

 [youtube.com/OticonPeopleFirst](https://youtube.com/OticonPeopleFirst)

 [facebook.com/OticonPeopleFirst](https://facebook.com/OticonPeopleFirst)

163245US / 07.15 Printed on 100% recycled paper



PEDIATRIC HEARING SOLUTIONS



**More speech details**  
- more power to achieve your goals





# Ready for the journey of life - with Sensei SP

## Get and give hugs without annoying whistling sounds

- Inium Sense feedback shield<sup>SP</sup> effectively controls whistling sounds

## Listen comfortably all day long

- Continuously high sound quality in a small-size hearing aid

## Hearing "s", "th" and "f" makes speech understanding easier

- Speech Rescue™ uniquely rescues high frequency sounds to make them audible

## Get more speech details than ever

- Speech Guard E protects the important details of speech

## Hear the teacher's voice in a noisy classroom

- FM compatible with dedicated FM options for school

## Keeping speech in focus

- Automatic TriState Noise Management dampens loud noise in the surrounding environment

## Stimulated to learn

- BrainHearing technology stimulates the brain to learn and support the child's development in making sense of sound

## Be ready for an active day outside

- IP58 classified for water- and dust-resistance

## Free of allergens

- Completely hypoallergenic, phthalates-free and biocompatible case colors

## Ideal for active kids

- Shock and drop resistant to withstand children's rough treatment

## Easy talking on the phone

- ConnectLine wireless connectivity makes it easy to use smartphones and listen to music

## Extra safety for parents

- Unique LED status indicator for battery function certainty  
- for parents and caregivers. Tamper-resistant battery door to keep curious fingers away from small batteries

## A better future for every child with hearing loss

### Giving your child's developing brain what it needs

At Oticon, we have helped people with hearing loss for more than a hundred years and would love to help your child too. Sound allows us to learn to talk. It allows us to have conversations, to express needs and wants and to connect with others. To maximize speech and language development, your child deserves the best hearing solution possible.

Although hearing loss happens in the ears, the real effect is in the brain. It is the brain that makes sense of sound and so, the brain must be stimulated to develop. Oticon's unique BrainHearing™ technology ensures a consistent sound processing approach that will give your child the best opportunity to get the most out of their hearing.



# Power your child's senses

**More speech details**  
**More power to achieve your goals**

**Discover how Oticon Sensei Super Power can give your child's developing brain what it needs.**

Oticon Sensei Super Power delivers the best POWER performance in the smallest hearing aid size, with the clearest sound signal. Your child's developing brain will get what it needs to maximize speech, language and social development. Improved communication skills can positively affect all areas of your child's life, including school performance, social relationships and even self-esteem.

Contact your hearing care professional to learn more about Sensei Super Power.

## Give your child's brain the most complete and clearest signal it needs

- Speech Rescue™ and Speech Guard E deliver more speech details than ever before.

Unique to Oticon, Speech Rescue™ and Speech Guard E are Sensei Super Power's innovative "Rescue" and "Guard" sound processing technologies that work in harmony to increase your child's understanding of speech.

Being able to hear the important high frequency sounds, such as "s" and "th", makes speech and language development easier. Moreover, the ability to experience variations among speakers, speech intensity and distance, as well as capturing as many speech details as possible will help improve your child's overall communication skills.

Sensei Super Power ensures access to the full spectrum of sounds and secures a strong speech input that is so vital for your child's cognitive growth. The benefits that make it easier for your child's brain to recognize sounds and take meaning from sound information are found in Oticon's BrainHearing™ concept. It ensures the brain gets the most complete and clearest signal it needs to let your child get the most out of speech.



## Get and give hugs without annoying whistling sounds

- Inium Sense feedback shield<sup>SP</sup> effectively controls whistling sounds

Your child wants to interact with friends, family and peers like everybody else, without being embarrassed by whistling sounds and annoying squealing sounds. Inium Sense feedback shield<sup>SP</sup> in Oticon Sensei Super Power keeps whistling to an absolute minimum, ensuring that your child has maximum wearing comfort.



## Complete peace of mind

- Unique LED status indicator

Sensei Super Power's unique LED status indicator gives you, your child's caregivers and teachers immediate, visual confirmation of the status of the hearing aids and lets you know for sure that the battery is functioning properly.



## Be ready for an active day outside

- IP58 classification for dust- and water-resistance

Jumping in puddles, playing sports, sweating and summer humidity are just some of the childhood situations hearing aids must endure. Oticon Sensei Super Power's mechanical design features multiple seals to prevent water and dust from entering the case, while all key components are nano coated to repel water and moisture. All ensured by Sensei Super Power's IP58 classification for dust- and water-resistance.



## Free of allergens

- meeting international standards on biocompatibility

When anything comes into contact with skin for long periods of time there is a chance of an allergic reaction occurring. Oticon has minimized this risk by excluding more than 200 potentially harmful substances and materials, such as phthalates. All materials in Oticon Sensei Super Power are tested for biocompatibility meeting strict international child safety standards and are completely hypoallergenic.



## Don't miss out - at school or at home

- ConnectLine wireless connectivity

ConnectLine is Oticon's complete, integrated system that gives your child easy wireless connection to hearing aids, smart phones, tablets, gaming systems, mp3 players, FM systems and more. With Oticon Streamer Pro, children and students can get their favorite music, phone calls and audio from TV - or the teacher's voice - streamed instantly and directly into their hearing aids.



## Hear the teacher's voice in a noisy classroom

- FM compatibility with dedicated FM options for school

For the classroom, Oticon Amigo FM is an optimal choice. This family of FM devices ensures your child gets the teacher's voice directly into the hearing aids. This direct FM signal gives your child an improved opportunity to listen and learn in class. Sensei Super Power is also compatible with Oticon ConnectLine wireless Streamer Pro and wireless systems from other manufacturers.

## Matching the colors of childhood

- A fun collection of decorative stickers

We have made sure the Oticon Sensei family of products is available in a wide range of modern colors, guaranteed to please children of all ages. Your child can also personalize Sensei Super Power with a fun collection of complementary stickers.

