

Your Information	Step 1: Account Information					Step 2: Patient Information									
	Ship to Account:					First Name:									
	Account Name:					Last Name:									
	Address:					Age:									
	City: State: Zip:					Last 4 digits of Social Security #:									
	Phone #:					Audiometric Information (required)									
	Contact Name:					Hz: 250 500 1k 2k 3k 4k 6k 8k									
	E-mail (required):					Right:									
	Bill to Account:					Left:									
	Date: Fitting Date:					Step 3: Fitter's Information									
Purchase Order #:					Fitter's Name:										
					Fitter's E-mail:										

4: Select instrument quantity, color, and speaker options

Model	Colors										Speaker						Domes						
	Diamond Black 63	Steel Grey 92	Silver Grey 91	Silver 44	Chroma Beige 90	Terracotta 94	Chestnut Brown 93	Honey Beige 116	Sky Blue 115	Ear	Power Level**	Length					OpenBass Dome		Bass dome, Double Vent		Power dome		
Oticon Intent™									R	60 85	0*	1	2	3	4	5	R	L	R	L	R	L	
									L	60 85							5 mm	<input type="checkbox"/>	<input type="checkbox"/>	---N/A---	---N/A---	<input type="checkbox"/>	<input type="checkbox"/>
										100							6 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
										100							8 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
																	10 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
																	12 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

*100 speaker level not available in size 0 length. **105 only available with MicroShell Detect

5. Style and product options

(Removal string standard for Solid (MicroMold), Hollow (LiteTip) and Embedded (MicroShell Detect))

SOLID (MicroMold)		HOLLOW (LiteTip)		EMBEDDED ² (Microshell Detect)	
Material	R L	Material	R L	Retention Styles	R L
Hard Acrylic (Transparent)	<input type="checkbox"/> <input type="checkbox"/>	Hard Acrylic (Transparent)	<input type="checkbox"/> <input type="checkbox"/>	Micro (IIC) (9mm ²)	<input type="checkbox"/> <input type="checkbox"/>
Soft Silicone (See color options)	<input type="checkbox"/> <input type="checkbox"/>	OtoTherm (Transparent)	<input type="checkbox"/> <input type="checkbox"/>	Canal (CIC) (12mm ²)	<input type="checkbox"/> <input type="checkbox"/>
OtoTherm™ (Transparent)	<input type="checkbox"/> <input type="checkbox"/>	Retention Styles		Semi-Skeleton	<input type="checkbox"/> <input type="checkbox"/>
Retention Styles		Micro (IIC) (9mm ²)	<input type="checkbox"/> <input type="checkbox"/>	Canal Lock	<input type="checkbox"/> <input type="checkbox"/>
Micro (IIC) (9mm ²)	<input type="checkbox"/> <input type="checkbox"/>	Canal (CIC) (12mm ²)	<input type="checkbox"/> <input type="checkbox"/>	Half Skeleton	<input type="checkbox"/> <input type="checkbox"/>
Canal (CIC) (12mm ²)	<input type="checkbox"/> <input type="checkbox"/>	Semi-Skeleton	<input type="checkbox"/> <input type="checkbox"/>	Skeleton Wing Lock (Same retention points as skeleton, less material)	<input type="checkbox"/> <input type="checkbox"/>
Semi-Skeleton	<input type="checkbox"/> <input type="checkbox"/>	Canal Lock	<input type="checkbox"/> <input type="checkbox"/>	Skeleton Lock	<input type="checkbox"/> <input type="checkbox"/>
Canal Lock	<input type="checkbox"/> <input type="checkbox"/>	Half Skeleton	<input type="checkbox"/> <input type="checkbox"/>	MicroShell Detect Fitting Level	
Half Skeleton	<input type="checkbox"/> <input type="checkbox"/>	Skeleton Wing Lock (Same retention points as skeleton, less material)	<input type="checkbox"/> <input type="checkbox"/>	60	<input type="checkbox"/> <input type="checkbox"/>
Skeleton Wing Lock (N/A Soft) (Same retention points as skeleton, less material)	<input type="checkbox"/> <input type="checkbox"/>	Skeleton Lock	<input type="checkbox"/> <input type="checkbox"/>	85	<input type="checkbox"/> <input type="checkbox"/>
Skeleton Lock	<input type="checkbox"/> <input type="checkbox"/>	miniFit Detect Fitting Level		100	<input type="checkbox"/> <input type="checkbox"/>
miniFit Detect Fitting Level		60	<input type="checkbox"/> <input type="checkbox"/>	105	<input type="checkbox"/> <input type="checkbox"/>
60	<input type="checkbox"/> <input type="checkbox"/>	85	<input type="checkbox"/> <input type="checkbox"/>	Vent Size	
85	<input type="checkbox"/> <input type="checkbox"/>	Vent Size		Max Vent	<input type="checkbox"/> <input type="checkbox"/>
Vent Size		Max Vent	<input type="checkbox"/> <input type="checkbox"/>	Extra Large > 2.4	<input type="checkbox"/> <input type="checkbox"/>
Max Vent	<input type="checkbox"/> <input type="checkbox"/>	Extra Large 1.5 (Default)	<input type="checkbox"/> <input type="checkbox"/>	Large 2.4 (Default 60, 85)	<input type="checkbox"/> <input type="checkbox"/>
Extra Large > 2.4	<input type="checkbox"/> <input type="checkbox"/>	Large 1.2	<input type="checkbox"/> <input type="checkbox"/>	Medium Plus 1.8	<input type="checkbox"/> <input type="checkbox"/>
Large 2.4 (Default)	<input type="checkbox"/> <input type="checkbox"/>	Medium 0.7 (short)	<input type="checkbox"/> <input type="checkbox"/>	Medium 1.4 (Default 100, 105)	<input type="checkbox"/> <input type="checkbox"/>
Medium 1.4	<input type="checkbox"/> <input type="checkbox"/>	Small 0.7 (long)	<input type="checkbox"/> <input type="checkbox"/>	Small 0.8	<input type="checkbox"/> <input type="checkbox"/>
Small 0.8	<input type="checkbox"/> <input type="checkbox"/>	No vent	<input type="checkbox"/> <input type="checkbox"/>	No vent	<input type="checkbox"/> <input type="checkbox"/>
No vent	<input type="checkbox"/> <input type="checkbox"/>				

¹Approximate length dependent on ear anatomy and impression.

²Only compatible with Sirius™ instruments.

³Impression will be destroyed unless otherwise requested.

⁴N/A with Matte Finish.

Color options

Soft Material Colors (Shore 40)

<input type="checkbox"/> Transparent	<input type="checkbox"/> Yellow
<input type="checkbox"/> White	<input type="checkbox"/> Green
<input type="checkbox"/> Black	<input type="checkbox"/> Red
<input type="checkbox"/> Rose	<input type="checkbox"/> Blue
<input type="checkbox"/> Orange	

MicroShell Detect only

<input type="checkbox"/> Transparent	<input type="checkbox"/> Dark Brown
<input type="checkbox"/> Beige	<input type="checkbox"/> Black
<input type="checkbox"/> Light Brown	<input type="checkbox"/> Red/Blue
<input type="checkbox"/> Medium Brown	

6. Mold Special Options

<input type="checkbox"/> Canal Tips Red/Blue ⁴ (Hard Acrylic Only)	<input type="checkbox"/> L&R on mold
<input type="checkbox"/> Heavy Removal Cord	<input type="checkbox"/> Large Ball Removal Cord

Optional Finish (Hard Coat STD or, choose 1 alternative below)

<input type="checkbox"/> Soft Coat (Hard Acrylic Only)	<input type="checkbox"/> Matte Finish (N/A OtoTherm) (N/A MicroShell)
--	---

Order contains **contaminated impression**. Handlers use gloves.

Tip Consider this style of Mold for the deepest/smallest mold available.

Tip Consider Semi-Skeleton as a great first option for patients in need of retention (short/narrow canals, difficulty inserting domes). This style provides better retention/stability than canal locks alone and offers better ease of insertion.

Please Send:

<input type="checkbox"/> Impression Boxes	<input type="checkbox"/> Service Ord. Forms
<input type="checkbox"/> Return Labels	<input type="checkbox"/> Return Impressions ³

Special Instructions:

24-570749 15500-0306/04.24