

Your Information	Step 1: Account Information				Step 2: Patient Information											
	Ship to Account:				First Name:											
	Account Name:				Last Name:											
	Address:				Age:											
					Last 4 digits of Social Security #:											
	City:		State:		Zip:		Audiometric Information (required)									
	Phone #:				Hz:		250	500	1k	2k	3k	4k	6k	8k		
	Contact Name:				Right:											
	E-mail (required):				Left:											
	Bill to Account:				Step 3: Fitter's Information											
Date:				Fitting Date:		Fitter's Name:										
Purchase Order #:				Fitter's E-mail:												

Step 4: Select instrument quantity, color, and speaker options																		
Model	Battery	Colors									Speaker							
		Diamond Black 63	Steel Grey 92	Silver Grey 91	Silver 44	Chroma Beige 90	Terracotta 94	Chestnut Brown 93	Hear Pink / Baby Pink 79	Olive Green 114	Ear	Power Level		Size				
		1	2	3	4	5	60	85	100	105		1	2	3	4	5		
Oticon Real 1 miniRITE R	Lithium-ion										L	60	85					
											R	100	105					
Oticon Real 1 miniRITE T	312 zinc-air										L	60	85					
										--N/A--	R	100	105					
Oticon CROS PX	Lithium-ion										L	60	85					
											R	60	85					
Oticon CROS	312 zinc-air										L	60	85					
										--N/A--	R	100	105					

Domes	OpenBass Dome		Double Vent Bass (1.4)		Power	
	L	R	L	R	L	R
	5 mm	<input type="checkbox"/>	<input type="checkbox"/>	---N/A---	---N/A---	---
6 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
10 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
12 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Grip Tips	L	R
Small, no vent	<input type="checkbox"/>	<input type="checkbox"/>
Small, 0.9mm vent	<input type="checkbox"/>	<input type="checkbox"/>
Large, no vent	<input type="checkbox"/>	<input type="checkbox"/>
Large, 0.9mm vent	<input type="checkbox"/>	<input type="checkbox"/>

Accessories	
ConnectClip	<input type="checkbox"/>
Remote Control 3.0	<input type="checkbox"/>
TV Adapter 3.0	<input type="checkbox"/>
Phone Adapter 2.0 (Requires ConnectClip)	<input type="checkbox"/>
EduMic	<input type="checkbox"/>

105 receiver is only available for use with a custom mold. Dome fitting is not an option.

Order contains **contaminated impression**. Handlers use gloves.

Impression(s) will be destroyed unless otherwise requested.

Oticon Chargers	
Oticon SmartCharger	<input type="checkbox"/>
Oticon Desktop Charger	<input type="checkbox"/>

Special Instructions: _____

Custom Mold Styles

Standard Mold		Hollow Mold		Embedded Mold	
Canal Mold	MicroMold (CIC style)	LiteTip		MicroShell	Power Receiver Mold

Available with various retention options from canal lock to skeleton.

Please make custom earmold selection on reverse side.

Standard (Speakers 60, 85)	L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	OtoTherm™**	Venting	Removal String		
	MicroMold (CIC style)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Please choose Venting selection options below.	STD	
	Canal	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
	Canal with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
	Canal with Helix Lock	<input type="checkbox"/>	<input type="checkbox"/>	STD	Not Available	<input type="checkbox"/>			
	Half Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
	Semi Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	STD	Not Available	<input type="checkbox"/>			
	Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Hollow Mold (LiteTip style)	L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	OtoTherm™**	Venting	Removal String		
Hollow	<input type="checkbox"/>	<input type="checkbox"/>	STD	Not Available	<input type="checkbox"/>	Please choose Venting selection options below.	<input type="checkbox"/>		
Hollow with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	STD	Not Available	<input type="checkbox"/>		<input type="checkbox"/>		
Hollow with Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	STD	Not Available	<input type="checkbox"/>		<input type="checkbox"/>		
Embedded									
MicroShell for 60 & 85 Speakers (embedded speaker)	L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	OtoTherm™**	Venting	Options		
Canal	<input type="checkbox"/>	<input type="checkbox"/>	STD	Not Available	Not Available	Please choose Venting selection options below.	MicroShell (acrylic only)	<input type="checkbox"/> L <input type="checkbox"/> 60 <input type="checkbox"/> 85 <input type="checkbox"/> R <input type="checkbox"/> 60 <input type="checkbox"/> 85	
Canal with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	STD	Not Available	Not Available		Receiver Length	<input type="checkbox"/> L <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5 <input type="checkbox"/> R <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5	
Canal with Helix Lock	<input type="checkbox"/>	<input type="checkbox"/>	STD	Not Available	Not Available		Shell and Faceplate color	<input type="checkbox"/> Clear <input type="checkbox"/> Beige <input type="checkbox"/> Black <input type="checkbox"/> Dark Brown <input type="checkbox"/> Med. Brown <input type="checkbox"/> Lt. Brown <input type="checkbox"/> Blue <input type="checkbox"/> Red	
Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	STD	Not Available	Not Available				
Semi Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	STD	Not Available	Not Available				
Half Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	STD	Not Available	Not Available				
Power Receiver Mold (100 and 105 speakers)	L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	OtoTherm™**	Venting	Options		
Power Receiver Mold	<input type="checkbox"/>	<input type="checkbox"/>	STD	Not Available	Not Available	Please choose Venting selection options below.	Power Receiver Mold (acrylic only)	<input type="checkbox"/> L <input type="checkbox"/> 100 <input type="checkbox"/> 105 <input type="checkbox"/> R <input type="checkbox"/> 100 <input type="checkbox"/> 105	
Power Receiver Mold with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	STD	Not Available	Not Available		Power Receiver Wire Length	<input type="checkbox"/> L <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5 <input type="checkbox"/> R <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5	
							Shell and Cover Plate Color	<input type="checkbox"/> Clear <input type="checkbox"/> Beige <input type="checkbox"/> Lt. Brown <input type="checkbox"/> Med. Brown <input type="checkbox"/> Dark Brown <input type="checkbox"/> Black	
							Removal String	<input type="checkbox"/> Yes <input type="checkbox"/> No	

Venting Selection															
All RITE molds except Hollow (LiteTip) Mold								Hollow (LiteTip) Mold RITE							
Venting based on audio-canal size <input type="checkbox"/>								Venting based on audio-canal size <input type="checkbox"/>							
Vent Size	Small (0.8mm)	Medium (1.4mm)	Large (2.4mm)	Max**	Pressure	SAV	No vent	Vent Size	Small (0.7mm long)	Medium (0.7mm short)	Large (1.2mm)	Max**	Pressure	SAV	No vent
L	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	L	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
R	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	R	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Retention Locks

Available options vary by style and material selected.

**As large as possible.

Canal Lock	Helix Lock	Skeleton Mold	Half Skeleton	Semi Skeleton
	(Not available in soft material) 			(Not available in soft material) 

Wax Protection System

Custom molds (except OtoTherm™, MicroShell, and Soft Silicone MicroMold) use the ProWax system. All miniFit receivers use the ProWax miniFit wax system.



*OtoTherm™ is flexible material that is rigid (Shore A 90) at room temperature and flexible (Shore A 70) at body temperature. These characteristics allow for ease of insertion, increased comfort, and optimal sealing.

Please note: This material requires gentle warming/softening of the mold with a hair dryer prior to insertion or removal of the receiver.