

Your Information	Step 1: Account Information										Step 2: Patient Information																	
	Ship to Account:										First Name:																	
	Account Name:										Last Name:																	
	Address:										Age:																	
											Last 4 digits of Social Security #:																	
	City:				State:			Zip:			Audiometric Information (required)																	
	Phone #:										Hz:		250		500		1k		2k		3k		4k		6k		8k	
	Contact Name:										Right:																	
	E-mail (required):										Left:																	
	Bill to Account:										Step 3: Fitter's Information																	
Date:					Fitting Date:					Fitter's Name:																		
Purchase Order #:										Fitter's E-mail:																		

Step 4: Select instrument quantity, color, and speaker options																			
Model	Battery	Colors										Speaker							
		Diamond Black 63	Steel Grey 92	Silver Grey 91	Silver 44	Chroma Beige 90	Terracotta 94	Chestnut Brown 93	Hear Pink / Baby Pink 79	Olive Green 114	Ear	Power Level		Size					
														1	2	3	4	5	
Oticon Real 1 miniRITE R	Lithium-ion										L	60	85						
											R	100	105						
Oticon Real 1 miniRITE T	312 zinc-air									--N/A--	L	60	85						
											R	100	105						
Oticon More 1 miniRITE R	Lithium-ion									--N/A--	L	60	85						
											R	100	105						
Oticon More 1 miniRITE T	312 zinc-air									--N/A--	L	60	85						
											R	100	105						
Oticon CROS PX	Lithium-ion									--N/A--	L	60	85						
											R	60	85						
Oticon CROS	312 zinc-air								--N/A--	--N/A--	L	60	85						
											R	60	85						

Domes						Grip Tips			L	R	Accessories	
	OpenBass Dome		Double Vent Bass (1.4)		Power	Small, no vent					ConnectClip	
	L	R	L	R								
5 mm	☐	☐	---N/A---		---N/A---	Small, 0.9mm vent					Remote Control 3.0	
6 mm	☐	☐	☐	☐	☐							
8 mm	☐	☐	☐	☐	☐	Large, no vent					TV Adapter 3.0	
10 mm	☐	☐	☐	☐	☐							
12 mm	☐	☐	☐	☐	☐	Large, 0.9mm vent					Phone Adapter 2.0 (Requires ConnectClip)	
											EduMic	

105 receiver is only available for use with a custom mold. Dome fitting is not an option.

Venting Information (required)								Oticon Chargers	
Vent Size	Small	Medium	Large	Max	SAV	Pressure	No vent	Oticon SmartCharger	
L	☐	☐	☐	☐	☐	☐	☐		
R	☐	☐	☐	☐	☐	☐	☐	Oticon Desktop Charger	

Special Instructions:



☐ Order contains contaminated impression. Handlers use gloves.

Impression(s) will be destroyed unless otherwise requested.

Standard (Speakers 60, 85)		L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	OtoTherm™**	Removal String	
MicroMold (CIC style)	☐	☐	☐	☐	☐		STD	
Canal	☐	☐	☐	☐	☐		☐	
Canal with Canal Lock	☐	☐	☐	☐	☐		☐	
Canal with Helix Lock	☐	☐	STD	Not Available	☐		☐	
Half Skeleton	☐	☐	☐	☐	☐		☐	
Semi Skeleton	☐	☐	STD	Not Available	☐		☐	
Skeleton	☐	☐	☐	☐	☐		☐	
Hollow Mold (LiteTip style)		L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	OtoTherm™**	Removal String	
Hollow	☐	☐	STD	Not Available	☐		☐	
Hollow with Canal Lock	☐	☐	STD	Not Available	☐		☐	
Hollow with Skeleton	☐	☐	STD	Not Available	☐		☐	
Embedded								
MicroShell for 60 & 85 Speakers (embedded speaker)		L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	OtoTherm™**	Options	
Canal	☐	☐	STD	Not Available	Not Available	MicroShell (acrylic only)	☐ L ☐ 60 ☐ 85	
Canal with Canal Lock	☐	☐	STD	Not Available	Not Available		☐ R ☐ 60 ☐ 85	
Canal with Helix Lock	☐	☐	STD	Not Available	Not Available	Receiver Length	☐ L ☐ 1 ☐ 2 ☐ 3 ☐ 4 ☐ 5	
Skeleton	☐	☐	STD	Not Available	Not Available		☐ R ☐ 1 ☐ 2 ☐ 3 ☐ 4 ☐ 5	
Semi Skeleton	☐	☐	STD	Not Available	Not Available	Shell and Faceplate color	☐ Clear ☐ Beige ☐ Black ☐ Dark Brown Translucent ☐ Med. Brown Translucent ☐ Light Brown ☐ Blue ☐ Red	
Half Skeleton	☐	☐	STD	Not Available	Not Available			
Power Receiver Mold (100 and 105 speakers)		L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	OtoTherm™**	Options	
Power Receiver Mold	☐	☐	STD	Not Available	Not Available	Power Receiver Mold (acrylic only)	☐ L ☐ 100 ☐ 105	
							☐ R ☐ 100 ☐ 105	
						Power Receiver Wire Length	☐ L ☐ 1 ☐ 2 ☐ 3 ☐ 4 ☐ 5	
							☐ R ☐ 1 ☐ 2 ☐ 3 ☐ 4 ☐ 5	
Power Receiver Mold with Canal Lock	☐	☐	STD	Not Available	Not Available	Shell color	☐ Clear	
						Cover Plate color	☐ Clear ☐ Beige ☐ Light Brown ☐ Medium Brown ☐ Dark Brown ☐ Black	
						Removal String	☐ Yes ☐ No	

Custom Mold Styles

Standard Mold		Hollow Mold		Embedded Mold	
					
Canal Mold	MicroMold (CIC style)	LiteTip		MicroShell	Power Receiver Mold

Available with various retention options from canal lock to skeleton.

Retention Locks

Available options vary by style and material selected.


Canal Lock	Helix Lock (Available in hard material only)	Skeleton Mold	Half Skeleton	Semi Skeleton (Available in hard material only)
				

Wax Protection System

Custom molds (except LiteTip, MicroShell, and Soft Silicone MicroMold) use the ProWax system. All miniFit receivers use the ProWax miniFit wax system.

ProWax

ProWax miniFit



*OtoTherm™ is flexible material that is rigid (Shore A 90) at room temperature and flexible (Shore A 70) at body temperature. These characteristics allow for ease of insertion, increased comfort, and optimal sealing.

Please note: This material requires gentle warming/softening of the mold with a hair dryer prior to insertion or removal of the receiver.