

Venting Information (required)								Notes
Vent Size	Small	Medium	Large	Max	SAV	Pressure	No vent	Power Molds default to 1.4mm vents unless otherwise requested.
L	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	LiteTip and MicroMolds come with largest vent possible unless otherwise requested.
R	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

MicroShell for 60 & 85 Speakers (embedded speaker)							
MicroShell (acrylic only)	<input type="checkbox"/> L <input type="checkbox"/> 60 <input type="checkbox"/> 85			Receiver Length	<input type="checkbox"/> L <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5		
	<input type="checkbox"/> R <input type="checkbox"/> 60 <input type="checkbox"/> 85				<input type="checkbox"/> R <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5		
Shell and Faceplate color			<input type="checkbox"/> Clear <input type="checkbox"/> Beige <input type="checkbox"/> Black				
Removal String/STD							
Molds for 60, 85	L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	*VarioTherm® silicone (clear only)	Removal String	
LiteTip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-----Not Available-----	<input type="checkbox"/>	STD	
LiteTip with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-----Not Available-----	<input type="checkbox"/>	<input type="checkbox"/>	
LiteTip Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-----Not Available-----	<input type="checkbox"/>	<input type="checkbox"/>	
MicroMold	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	STD	
Canal	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Canal with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Canal with Helix Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-----Not Available-----	-----Not Available-----	<input type="checkbox"/>	
Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Semi Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-----Not Available-----	-----Not Available-----	<input type="checkbox"/>	
Half Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
miniFit Power Molds for 100/105 speaker	L	R					
Power Receiver Mold	<input type="checkbox"/>	<input type="checkbox"/>	STD	-----Not Available-----	-----Not Available-----	<input type="checkbox"/>	
Power Receiver Mold with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	STD	-----Not Available-----	-----Not Available-----	<input type="checkbox"/>	
Clear shell - Face plate color:	<input type="checkbox"/> Clear <input type="checkbox"/> Beige <input type="checkbox"/> Light Brown <input type="checkbox"/> Medium Brown <input type="checkbox"/> Dark Brown <input type="checkbox"/> Black						
Power Receiver Wire Length	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5				

Custom Mold Styles

LiteTip	MicroMold	Power Receiver Mold	MicroShell	Skeleton	Half Skeleton	Semi Skeleton
						
<i>(Hollow)</i>	<i>(Solid)</i>	<i>Opn, Opn S and Opn Play ONLY</i>	<i>Opn S and Opn Play ONLY</i>			<i>(Available in hard material only)</i>

Retention Locks

All mold styles except MicroShell, are offered with canal locks and skeleton locks for better retention. Helix locks, Half Skeleton, and Semi Skeleton styles are also available upon request. (Retention Lock options vary by style and material selected.)

Canal Lock	Helix Lock
	
	<i>(Available in hard material only)</i>

***VarioTherm Material**

VarioTherm is a body temperature-reactive plastic for soft flexible earmolds. The unique properties of this material allow the mold to remain hard at room temperature for ease of insertion and then soften (become more flexible) at body temperature for increased comfort and optimum sealing.

Please note: This material requires gentle warming/softening of the mold with a hair dryer before insertion or removal of the receiver.

*miniFit speakers only

Wax Protection System

Custom molds (except LiteTips) use the ProWax system. All miniFit receivers and MicroShells use the ProWax miniFit wax system.

