

Your Information	Step 1: Account Information								Step 2: Patient Information								
	Ship to Account:								First Name:								
	Account Name:								Last Name:								
	Address:								Age:								
									Last 4 digits of Social Security #:								
	City:				State:		Zip:		Audiometric Information (required)								
	Phone #:								Hz:	250	500	1k	2k	3k	4k	6k	8k
	Contact Name:								Right:								
	Email (required):								Left:								
	Bill to Account:								Step 3: Fitter's Information								
Date:				Fitting Date:				Fitter's Name:									
Purchase Order #:								Fitter's Email:									

New Instrument																						
Model	Colors														Speaker							
	Diamond Black 63	Steel Grey 92	Silver Grey 91	Silver 44	Chroma Beige 90	Terracotta 94	Chestnut Brown 93	Hear Pink / Baby Pink 79	Emerald Green 48	Blue 47	Purple 45	Red 46	Aquamarine 58	Power Pink 57	Ear	Power Level		Size				
	1	2	3	4	5																	
Oticon More™ miniRITE R (Lithium-ion)									-----Not Available-----					L	60	85						
									-----Not Available-----					R	100	105						
Oticon CROS									-----Not Available-----					L	60							
									-----Not Available-----					R	60							
Opn S 1 miniRITE									-----Not Available-----					L	60	85						
									-----Not Available-----					R	100	105						
Opn S 1 miniRITE R (Lithium-ion)									-----Not Available-----					L	60	85						
									-----Not Available-----					R	100	105						
Opn S 1 miniRITE T									-----Not Available-----					L	60	85						
									-----Not Available-----					R	100	105						
Opn Play 1 miniRITE T		N/A	N/A											L	60	85						
													R	100	105							

Grip Tips	L	R	Accessories	Domes	Notes
Small, no vent	<input type="checkbox"/>	<input type="checkbox"/>	ConnectClip <input type="checkbox"/>	OpenBass Dome** L R	105 receiver is only available for use with a custom mold. Dome fitting is not an option. ** For Oticon More only *** For Oticon Opn S only
Small, 0.9mm vent	<input type="checkbox"/>	<input type="checkbox"/>	Remote Control 3.0 <input type="checkbox"/>	Open*** L R	
Large, no vent	<input type="checkbox"/>	<input type="checkbox"/>	TV Adapter 3.0 <input type="checkbox"/>	Single Vent Bass (0.8)*** L R	
Large, 0.9mm vent	<input type="checkbox"/>	<input type="checkbox"/>	Phone Adapter 2.0 (Requires ConnectClip) <input type="checkbox"/>	Double Vent Bass (1.4) L R	
				Power L R	
				5mm	
				6mm	
				8mm	
				10mm	
				12mm	

Venting Information (required)								Notes
Vent Size	Small	Medium	Large	Max	SAV	Pressure	No vent	LiteTip and MicroMolds come with largest vent possible unless otherwise requested.
L	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
R	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Special Instructions: _____





MicroShell for 60 & 85 Speakers (embedded speaker)						
MicroShell (acrylic only)	<input type="checkbox"/> L <input type="checkbox"/> 60 <input type="checkbox"/> 85	Receiver Length		<input type="checkbox"/> L <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5		
	<input type="checkbox"/> R <input type="checkbox"/> 60 <input type="checkbox"/> 85			<input type="checkbox"/> R <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5		
Shell and Faceplate color	<input type="checkbox"/> Clear <input type="checkbox"/> Beige <input type="checkbox"/> Black		Removal String (Standard)			
Molds for 60, 85	L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	*VarioTherm® silicone (clear only)	Removal String
LiteTip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-----Not Available-----	<input type="checkbox"/>	Standard
LiteTip with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-----Not Available-----	<input type="checkbox"/>	<input type="checkbox"/>
LiteTip Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-----Not Available-----	<input type="checkbox"/>	<input type="checkbox"/>
MicroMold	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Standard
Canal	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Canal with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Canal with Helix Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-----Not Available-----	-----Not Available-----	<input type="checkbox"/>
Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Semi Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-----Not Available-----	-----Not Available-----	<input type="checkbox"/>
Half Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
miniFit Power Molds for 100/105 speaker	L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	*VarioTherm® silicone (clear only)	Removal String
Power Receiver Mold	<input type="checkbox"/>	<input type="checkbox"/>	Standard	-----Not Available-----	-----Not Available-----	<input type="checkbox"/>
Power Receiver Mold with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	Standard	-----Not Available-----	-----Not Available-----	<input type="checkbox"/>
Clear shell - Faceplate color:	<input type="checkbox"/> Clear <input type="checkbox"/> Beige <input type="checkbox"/> Light Brown <input type="checkbox"/> Medium Brown <input type="checkbox"/> Dark Brown <input type="checkbox"/> Black					
Power Receiver Wire Length	<input type="checkbox"/> L <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5					
	<input type="checkbox"/> R <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5					

Custom Mold Styles

LiteTip	MicroMold	Power Receiver Mold	MicroShell	Skeleton	Half Skeleton	Semi Skeleton
						
(Hollow)	(Solid)	More, Opn S, Opn and Opn Play ONLY	More, Opn S, Opn and Opn Play ONLY			(Available in hard material only)

Retention Locks

All mold styles except MicroShell, are offered with canal locks and skeleton locks for better retention. Helix locks, Half Skeleton, and Semi Skeleton styles are also available upon request. (Retention Lock options vary by style and material selected.)

Canal Lock	Helix Lock
	 (Available in hard material only)

*VarioTherm Material

VarioTherm is a body temperature-reactive plastic for soft flexible earmolds. The unique properties of this material allow the mold to remain hard at room temperature for ease of insertion and then soften (become more flexible) at body temperature for increased comfort and optimum sealing.

Please note: This material requires gentle warming/softening of the mold with a hair dryer before insertion or removal of the receiver.

*miniFit 60 and 85 speakers only

Wax Protection System

Custom molds (except LiteTips) use the ProWax system. All miniFit receivers and MicroShells use the ProWax miniFit wax system.



20-43150 15500-2094/03.21